

# FLUX Walker



**VERSATILE POSTERIOR  
WALKING AID**

Distributed by  
**medifab**

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

November 2023

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Size Options

Measurements	Size 1	Size 2	Size 3	Size 4
Standard Grab Rail Grip Height	48 - 61 cm	57 - 81 cm	69 - 87 cm	82 - 102 cm
Standard Grab Rail Grip Width	41 cm	46 cm	51 cm	51 cm
Universal Grip Grab Rail Grip Height	43 - 68 cm	51 - 88 cm	62 - 94 cm	75 - 102 cm
Universal Grip Grab Rail Grip Width	29 - 41 cm	34 - 46 cm	39 - 51 cm	39 - 51 cm
Total Frame Width	64 cm	69 cm	74 cm	78 cm
Total Frame Length	65 cm	80 cm	88 cm	91 cm
Total Load Capacity	60 kg	75 kg	75 kg	100 kg
Front / Rear Wheel Diameter	15 / 20 cm	15 / 20 cm	15 / 20 cm	20 / 30 cm

## Package & Size Options

Package Includes	Size 1	Size 2	Size 3	Size 4
<ul style="list-style-type: none"> <li>Standard Handlebar</li> <li>Front Swivel Wheels</li> <li>Polyurethan Tyres</li> <li>Bumper Wheels</li> <li>Friction Brakes</li> <li>Anti-Reverse Mechanism</li> </ul>	5641-0001-000	5642-0001-000	5643-0001-000	5644-0001-000
				
	Luminous Red	Sky Blue	Gentian Blue	Grey Aluminium

# Accessories

## Anti-Tippers without Castors

**Size 1-3**  
5640-8421-000



**Size 4**  
5640-8424-000



## Anti-Tippers with Castors

**Size 1-3**  
5640-8425-000




**Size 4**  
5640-8428-000



## Fold Up Seat

**Size 1-2**  
5640-0201-000

**Size 3-4**  
5640-0204-000



## Carrier

**Size 1**  
5640-8061-000

**Size 2**  
5640-8063-000



## Sling Seat

**Small**  
Hip Size: 57-77cm  
5640-1921-000

**Large**  
Hip Size: 73-103cm  
5640-1922-000



## Sling Seat Rail

**Size 1**  
5640-1931-000

**Size 2**  
5640-1932-000

**Size 3**  
5640-1933-000



## Flux Basket

5640-8051-000



## Hand Disc

5640-5550-000



## Driving Brake

**Size 1-3**  
Includes wheel lock  
5640-8451-000

**Size 4**  
Includes wheel lock  
5640-8454-000






## Swivel Lock



To hold castors straight when brisk walking.  
2880-8440-000



# Accessories Continued

Dynamic Hip Control (Hip Belt Required)		Hip Belt (Required)	
<b>Suspension Soft (Yellow)</b> 5640-1952-000		<b>Small, Hip 53-83cm</b> 5640-6201-000	
<b>Suspension Hard (Red)</b> 5640-1953-000		<b>Large, Hip 65-115cm</b> 5640-6202-000	
<b>Kit for Orthoses</b> 5640-1951-000			

Universal Handlebar			
<b>Size 1</b> 5641-5751-000	<b>Size 2</b> 5642-5752-000		
<b>Size 3</b> 5643-5753-000	<b>Size 4</b> 5644-5754-000		

Forearm Supports - 2 Axial Adjust		Hip Pad	
22 x 10.5cm 5640-5501-000		Multi-Adjustable 5640-1900-000	

## Notes

---



---



---



---



---



---



---



---



---



---

# medifab

Distributed by Medifab Ltd

22 Detroit Drive  
Rolleston 7675  
New Zealand

0800 543 343

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Flux\_20231110\_NZ