

Hero-XL




**SPECIAL
PURPOSE
CAR SEAT**

SCRIPTING AND
ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Standard Package

Hero-XL Car Seat Package

Hero-XL extended rear facing special purpose car seat in Grey/Black complete with accessory mounting frame, adjustable lateral supports and installation straps for lower tether points.

4182-0010-090



Specifications

Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Distance Between Laterals	Total Width	Capacity
27cm	24cm	25-34cm	28-45cm	13-20cm	43cm	36kg

See next page for accessory options.

Accessory Options

Seat Depth Adaptor

One Size
4182-0280-000



Medial Thigh Support

One Size
4182-1510-000



Combination Medial Thigh Support and Depth Adaptor

One Size
4182-0895-000



Tray Table

One Size
4182-5544-000



Additional Lower Tether Anchorage Straps

A set of these straps is provided as a part of the standard package, additional straps available for installation in a second vehicle.

4182-6953-000



Seat Protector Mat

One Size
4182-0190-000

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script Hero-XL_20240925