

Jenx Multistander



**PRONE OR SUPINE
STANDING AID**

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

October 2024


Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Notes

User Measurements

A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	



Basic Packages

Multistander Size 1, Lime (Prone)	Multistander Size 1, Lime (Supine)
-----------------------------------	------------------------------------

Order Code 5230-0011-054	Order Code 5230-0013-054
<p>Option one, for prone set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 1 • Manual frame • Size 1 sandals • Size 1 knee blocks • Small hip support pads • Small trunk support pads and tray • MUS01-P1-03 (add headrest of choice to convert to supine) 	<p>Option one, for supine set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 1 • Manual frame • Size 1 sandals • Size 1 knee blocks • Small hip support pads • Small trunk support pads • Tray and 1 link head support

Multistander Size 2, Lime (Prone)	Multistander Size 2, Lime (Supine)
-----------------------------------	------------------------------------

Order Code 5230-0022-054	Order Code 5230-0021-054
<p>Option two, for prone set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 2 • Manual frame • Size 2 sandals • Size 2 knee blocks • Medium hip support • Medium trunk support and tray 	<p>Option two, for supine set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 2 • Manual frame • Size 2 sandals • Size 2 knee blocks • Slim thoracic • Medium hips supports • Moulded head rest • Multigrip body support and tray • MUS02-P1-03

Product Details

Product Specifications	Size 1	Size 1 - With Abduction	Size 2
Average age	9 Months – 6 Years	9 Months – 6 Years	3 – 13 Years
Weight capacity (kg)	35	35	60
Hip width (cm)	14 – 27.5	14 – 27.5	20 – 32
Hip height (cm)	25 – 44	25 – 44	37.5 – 60
Chest width (cm)	14 – 27.5	14 – 27.5	20 – 32
Chest height (cm)	49 – 85	50 – 89	75 – 120
Height range (cm)	62 – 116	60 – 116	110 – 155
Tilt range (°) (90 = upright)	10 – 90	10 – 90	0 – 85 (Powered) / 10 – 85 (Manual)
Frame width (cm)	61.5	61.5	67.5
Frame length (cm)	76	76	107
Abduction range (°)	-	30 (Either Leg)	-
Height of surface when horizontal (cm)	56	56	77.5

*Max ranges are achieved with additional accessories


Key

Blue	Suitable for Prone and Supine
Green	Suitable for Supine only
Yellow	Suitable for Prone only





Items outlined in blue are recommended as part of standard package


Power & Size Options

Manual		Battery Powered	
Size 1 5230-0001-000		Size 2 5230-0008-000	
Size 2 5230-0002-000			


Colour Options

Blue		Lime	
-066		-054	



Abduction Accessory



Only Available on Size 1	
Supports 0-30° in each leg in 5° increments. Fitted prior to shipping. 5230-5000-000	

Support Boards

Size 1	Size 2	
PU chest & hip support boards. 450mm (W) 5230-2201-XXX	PU chest & hip support boards. 520mm (W) 5230-2202-XXX	



Head Supports

Moulded Head Support		Oval Head Support	
Moulded polyurethane headrest for children who require minimal to moderate head support. 280mm (W) x 210mm (H) 5230-4101-XXX		A simple support for children who require minimal head positioning. 190mm (W) x 100mm (H) 5356-4611-XXX	

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm (W) x 5cm (D) x 12cm (H) 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm (W) x 5cm (D) x 12cm (H) 5356-4603-000	
Wipeable Cover (optional) 5356-4671-000	Wipeable Cover (optional) 5356-4672-000	Wipeable Cover (optional) 5356-4673-000	

Pelvic Supports

Support Pad & Straps, Small	Support Pad & Straps, Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	

Flexi Thoracics (Pair)			
Size 1 A 180mm B 80mm 5356-3333-000	Size 1 Long  A 215mm B 80mm 5356-3334-000	Size 2 A 220mm B 110mm 5356-3335-000	
Covers (required) 5356-3363-XXX	Covers (required) 5356-3364-XXX	Covers (required) 5356-3365-XXX	

Hi-Rise Thoracic Laterals

Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX	
--	---

Trunk Supports

Support Pad & Straps, Small	Support Pad & Straps, Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	


Flexi Thoracics (Pair)			
Size 1 A 180mm B 80mm 5356-3333-000	Size 1 Long A 215mm B 80mm 5356-3334-000		Size 2 A 220mm B 110mm 5356-3335-000
Covers (required) 5356-3363-XXX	Covers (required) 5356-3364-XXX		Covers (required) 5356-3365-XXX

Hi-Rise Thoracic Laterals	Multigrip Body Support
Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX	 For Size 2 only. Provides extra support on shoulders or lumbar area. 420mm (W) 5230-2721-XXX

Sternal Prompt
(For Prone). 5230-2910-000
Cover (required) 5230-8801-XXX

Waistcoat Harness

Size 2	Size 3	Size 4	Mounting Kit
Provides additional trunk and shoulder stability.			
A 5cm B 32cm C 14cm D 15.5cm 5230-6712-000	A 6.5cm B 37cm C 17.5cm D 18cm 5230-6713-000	A 8.5cm B 39.5cm C 20.5cm D 18.5cm 5230-6714-000	Mounting bracket for Harness (Required for Size 1 Multistander with Abduction only). 5230-PART-056



Knee Supports


Size 1	Size 2	
Easy to attach and adjust for height, angle & depth.		
A 100mm / B 45mm / C 75mm 5356-1751-XXX	A 120mm / B 60mm / C 95mm 5356-1752-XXX	


Tray Options

Elbow Blocks	
Suits all sizes. 10cm (W) x 14cm (H) x 8.5cm (D) 5356-7351-000	

Adjustable Tray - Size 1	Adjustable Tray - Size 2	
A 650mm B 450mm C 370mm D (With Insert) 330mm / (Without Insert) 290mm 5230-7111-000	A 700mm B 530mm C 460mm D (With Insert) 410mm / (Without Insert) 360mm 5230-7112-000	

Foot Positioning

Sandals			
Size 1 170mm (L) x 80mm (W) 5356-6081-000	Size 2 220mm (L) x 95mm (W) 5356-6082-000	Size 3 275mm (L) x 105mm (W) 5356-6083-000	

Sandal Raiser Blocks			
Size 1 Kit of 4x15mm blocks 5356-6091-000	Size 2 Kit of 4x15mm blocks 5356-6092-000	Size 3 Kit of 4x15mm blocks 5356-6093-000	

Ankle Snugs		
Size 1 120mm - 225mm Circumference 120mm (H) 5356-6131-000	Size 2 165mm - 295mm Circumference 120mm (H) 5356-6132-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Multistander_20241031_NZ