

Jenx Standz



ABDUCTION STANDER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

October 2023

Client Details


Client/User Name:

Date:

NDIS/Funding:

User Measurements

| | |
|-----------------------|----------------------|
| A Hip width | <input type="text"/> |
| B Axilla height | <input type="text"/> |
| C Chest width | <input type="text"/> |
| D Overall user height | <input type="text"/> |
| E Heel to knee | <input type="text"/> |
| F Hip height | <input type="text"/> |
| User weight | <input type="text"/> |



Basic Packages

Standz Size 1, Yellow (Prone)

5240-0011-044

Includes:

- Gas Tilt,
- Angle Adjustable Footplates
- Size 2 Sandals
- Standard Knee Supports
- Pelvic Support Pads
- Trunk Support Pads
- Wipe-clean Yellow vinyl covers
- Tray and Tray extension wings

Standz Size 1, Yellow (Supine)

5240-0020-044

Includes:

- Gas Tilt
- Angle Adjustable Footplates
- Size 1 Sandals
- Knee Supports with Multi-axis technology and incorporated flexed knee pads
- Pelvic Support Pads
- Trunk Support Pads
- 1 Link Multigrip Head Support
- Wipe-clean Yellow vinyl covers
- Tray with extension wings and a tray holder

Standz Size 1, Yellow (Opt 2, Supine)

5240-0021-044

Includes:

- Gas Tilt
- Angle Adjustable Footplates
- Size 2 Sandals
- Knee Supports with Multi-axis technology and incorporated flexed knee pads
- Pelvic Support Pads
- Trunk Support Pads
- Flat shoulder support
- 2 Link Multigrip Head Support
- Wipe-clean Yellow vinyl covers
- Tray with extension wings and a tray holder

Product Details

Size 1 - Manual Tilt

| | |
|-----------------------------------|-----------------|
| Average Age | 1 - 9 Years |
| Weight Capacity | 45 kg |
| Hip Width | 17.5cm - 30.5cm |
| Hip Height | 45cm - 70cm |
| Chest Width | 17.5cm - 30.5cm |
| Chest Height | 58cm - 100cm |
| Max Height | 138cm |
| Tilt Range (90° = Upright) | 5° to 90° |
| Frame Width | 64.5cm |
| Frame Length | 89cm |
| Height of surface when horizontal | 81cm |
| Abduction Range | 30° Either Leg |

Key

| | |
|--------|--------------------------------------|
| Blue | Suitable for Prone and Supine |
| Green | Suitable for Supine only |
| Yellow | Suitable for Prone only |

Items outlined in blue are required and an option must be selected

Base unit



Standz Basic Package



Tool free adjustment standing frame includes base with castors, hydraulic assisted tilt, pelvic and trunk support pads with strap, abducting legs, foot positioning plates, tray mount.

5240-0005-044



Knee support options

| Standard Knee Support | | Multi-Axis Knee Support | |
|--|---|--|---|
| Dual wall knee supports with accessories available for extra function (see below). |  | Width, height, depth and angle adjustable, with soft arched foam pads behind the knee and support straps to enable knee flexion. |  |
| 5240-1750-044 | | 5240-1770-044 | |
| A: 11.5cm, B: 12.5cm, C: 9.5cm & Min: 6.5cm | | A: 13.5cm & B: 11.5cm | |

| Knee Blocks | | Knee Block Mounting Bracket | |
|---|---|---|---|
| Flexible PU knee cups, complete with padded knee straps. Adjustable for height, width, angle and depth. | | | |
| Knee Blocks Size 1 Internal Support Size: 7.5cm W x 10cm H |  | Mounting bracket required when using knee block option for knee support. |  |
| 5356-1751-011 | | | |
| Knee Blocks Size 2 Internal Support Size: 9.5cm W x 12cm H | | | |
| 5356-1752-011 | | 5240-1720-000 | |

Accessories for Standard Knee Support

| Standard Strap (Required for Supine use) | De-rotational Strap | Arched Foam Pads (Recommended for Supine use) | Qty |
|---|--|--|-----|
| Holds knee in position and prevents knee flexion. Can also be used in Prone. | Adjustable strap for reducing internal or external rotation. | Ensures a comfortable standing position when accommodating knee flexion. Supplied individually. 13.5cm W x 13.5 cm H x 5.5cm D | |
|  |  |  | |
| 5240-1766-044 | 5240-1786-044 | 5353-1775-000 | |

Foot support options

| Size 1 sandals | | Size 2 sandals | |
|-----------------------------|---|-----------------------------|---|
| Ensure solid foot placement |  | Ensure solid foot placement |  |
| 5356-6081-000 | | 5356-6082-000 | |
| 8cm x 17cm long | | 9.5cm x 22cm long | |

Head support option

Flat Pad

For children with reasonable head control. Can also be used as chest pad in prone.
5240-4110-044

19cm W x 10cm H



Multigrip - 1 Link

Single link wings for some degree of lateral support.
33cm W x 5cm D x 7.5cm H
5356-4601-000

Wipeable Cover (Optional)
5356-4671-000

Multigrip - 2 Link

2 Links in each wing for extra support and positioning.
47.5cm W x 5cm D x 12cm H
5356-4602-000

Wipeable Cover (Optional)
5356-4672-000

Multigrip - 3 Link

3 adjustable links in each wing for ultimate lateral headrest control and comfort.
62cm W x 5cm D x 12cm H
5356-4603-000

Wipeable Cover (Optional)
5356-4673-000



Shoulder support option

Flat Shoulder Support

Shoulder support for supine use only (too wide for prone)

5240-2750-044

31.5cm W x 9.5cm H



Multigrip body support (for prone and supine)

Adjustable contour shoulder support

5240-2721-044

41cm max width



Pelvic support options

Pelvic Support Board

Lower back support

5240-1901-044

19.5cm W x 15cm H



Back Pad (Requires Pelvic support board)






Central back support

5240-2402-044

19.5cm W x 15cm H



Accessories

| Anti Tamper Conversion Kit | | Multi-position tray (mounts to bracket in base package) | |
|--|---|--|---|
| Allows for all wing knob adjustments to be converted to a tool fixing – ideal for environments where easy adjustment can be problematic. |  | Includes extension wings with slot. |  |
| 5240-0170-000 | | 5240-7112-044 | |
| | | A: 58.5cm & B: 39.5cm C: (Prone - 25cm, Supine - 45cm) | |
| Tray holder (useful for Supine) | | Tray Elbow Blocks | |
| When in Supine, the tray needs to be removed to remove the child. Place the tray on the convenient tray holder to avoid tripping. |  | Aids in supporting the arms and promoting a functional arm and hand position. Only usable in conjunction with Tray Extension Wings. |  |
| 5240-7365-000 |  | 5356-7351-000 | |
| | | 10cm W x 14cm H | |

Medifab Notes

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Standz_20231025_NZ