

Jenx Standz 2



ABDUCTION STANDER

SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:

Date:

DOB:

User Measurements

A Hip width

B Axilla height


C Chest width

D Overall user height

E Heel to knee

F Hip height

User weight



Product Specifications

Average age

5 - 14 Years

Weight capacity

70 kg

Hip width

25 - 37 cm

Hip height

50 - 95 cm

Chest width

22 - 36 cm

Chest height

82 - 116 cm

Max height

160 cm

Tilt range (90° = upright)

0° - 90°

Frame width

76 cm

Frame length

107 cm

Height off surface when horizontal

90 cm

Abduction range

30° Either Leg

Items outlined in blue are recommended as part of standard package

Base unit

Standz Basic Package


- 4 lockable swivel castors
- Tilt mechanism with dedicated release lever on handle with safety lock
- 0° - 30° independent abduction movement of leg posts in 5 degree increments with easy to use locking plunger
- Individual footplates with angle adjustment for plantar and dorsiflexion
- Tilt range of 5° - 90° with easy to read angle gauge
- Pelvic and trunk support boards and cover pads
- Includes tray mounting bracket

5242-0005-054






Foot support options


Sandals				
Multi-adjustable foot supports to ensure secure placement and positioning of feet.	Size 1 8cm W x 17cm L 5356-6081-000	Size 2 9.5cm W x 22cm L 5356-6082-000	Size 3 10.5cm W x 27.5cm L 5356-6083-000	

Ankle Snugs			
Additional ankle support for users that need it, featuring unique magnetic buckle fastening.	Size 1 Inner Circumference: 12-22.5cm 5356-6131-000	Size 2 Inner Circumference: 16.5-29.5cm 5356-6132-000	

Knee support options

Multi-Axis Knee Supports			
Width, height, depth and angle adjustable supports with soft arched foam pads behind the knee and support straps. Ideal for accommodating knee flexion.	Size 1 5242-1770-054	Size 2 11.5cm W x 13.5cm H 5242-1771-054	

Knee Blocks		Knee Block Mounting Bracket	
Flexible PU knee cups, complete with padded knee straps. Adjustable for height, width, angle and depth.		 Mounting bracket required when using knee block option for knee support. 5240-1720-000	
Size 1 Internal Support Size: 7.5cm W x 10cm H 5356-1751-011			
Size 2 Internal Support Size: 9.5cm W x 12cm H 5356-1752-011			


Standard Knee Support	
Width, height, depth and angle adjustable supports. Provides width adjustable walls to correctly gutter the knee position. <i>Must be used with knee strap accessories.</i> A: 11.5cm B: 12.5cm C: 6.5 - 9.5cm 5242-1750-054	

Key Blue Prone and Supine Green Supine Only Yellow Prone Only

Knee support options (continued)

Standard Knee Straps			
<p>Holds knee in position and prevents knee flexion in supine.</p> <p><i>Only compatible with Standard Knee Support.</i></p>	<p>Size 1</p> <p>Vinyl Pad Length: 16.5 - 21cm</p> <p>5242-1766-054</p>	<p>Size 2</p> <p>Vinyl Pad Length: 24 - 36cm</p> <p>5242-1767-054</p>	

Rotation Control Knee Straps			
<p>Knee cuff and strap system suitable for adjusting internal and external knee rotation.</p> <p><i>Only compatible with Standard Knee Support.</i></p>	<p>Size 1</p> <p>Cuff Circumference: 38cm Max</p> <p>5242-1786-054</p>	<p>Size 2</p> <p>Cuff Circumference: 45cm Max</p> <p>5242-1787-054</p>	

Arched Foam Pads for Standard Knee Support		Qty
<p>(Recommended for Supine use).</p> <p>Ensures a comfortable standing position when accommodating knee flexion. Supplied individually.</p> <p>13.5cm W x 13.5 cm H x 5.5cm D</p> <p>5353-1775-000</p>		

Pelvic support options

Flip-away Lateral Support with Strap		Flexi Flip Away Lateral Support with Strap	
<p>Firm lateral support pads with incorporated support strap. Flip away to aid transfers.</p> <p>Support Pads: 9.5cm H x 12cm D</p> <p>5242-3371-054</p>		<p>Flexible lateral supports with incorporated support strap. Contour to the body to give firm yet comfortable positioning. Flip away to aid transfers.</p> <p>Support Pads: 11.5cm H x 18.5cm D</p> <p>5242-3372-054</p>	

Trunk support options

Flip-away Lateral Support with Strap		Flexi Flip Away Lateral Support with Strap	
<p>Firm lateral support pads with incorporated support strap. Flip away to aid transfers.</p> <p>Support Pads: 9.5cm H x 12cm D</p> <p>5242-3371-054</p>		<p>Flexible lateral supports with incorporated support strap. Contour to the body to give firm yet comfortable positioning. Flip away to aid transfers.</p> <p>Support Pads: 11.5cm H x 18.5cm D</p> <p>5242-3372-054</p>	

Key Blue Prone and Supine Green Supine Only Yellow Prone Only

Shoulder support options

Flat shoulder support	Shoulder support - Size 1	Shoulder support - Size 2
<p>A simple support pad for the shoulders. 31.5cm W x 9.5cm H 5242-2750-054</p>	<p>Shoulder protractors mounted to a flat shoulder pad. Offers angle adjustment and promotes shoulder girdle stability.</p>	
	<p>Back Pad: (22cm W x 9cm H) Protractor: (15.5cm W x 13cm H) 5242-2751-054</p>	<p>(Required to achieve maximum user height). Back Pad: (22cm W x 19cm H) Protractor: (15cm W x 21.5cm H) 5242-2752-054</p> 

Head support options









Oval Head Support	
<p>Simple headrest for children who require minimal head support. Comes mounted on an adjustable bracket. 19cm W x 10cm H 5356-4611-054</p>	

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
<p>Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000</p>	<p>2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000</p> 	<p>3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm W x 5cm D x 12cm H 5356-4603-000</p>	
<p>Wipeable Cover (Optional) 5356-4671-000</p>	<p>Wipeable Cover (Optional) 5356-4672-000</p>	<p>Wipeable Cover (Optional) 5356-4673-000</p>	

Optional Accessories

Pelvic Support Board	Pelvic Support Board Back Pad
<p>Offers a full contact surface to facilitate hip extension and neutral pelvic positioning. Pad Size: 19.5cm W x 15cm H 5242-1903-054</p>	<p>Offers support to the child's upper back. May be used in combination with a head rest to limit unwanted hyper extension in prone. Only usable with pelvic support board. Pad Size: 19.5cm W x 15cm H 5242-1902-054</p> 

Optional Accessories (continued)

Anti Tamper Conversion Kit		Multi-position tray (mounts to bracket in base package)	
<p>Allows for all wing knob adjustments to be converted to a tool fixing – ideal for environments where easy adjustment can be problematic.</p> <p>5240-0170-000</p>		<p>Can be used in both prone and supine. The tray is angle adjustable to ensure the correct position in either application.</p> <p>5242-7110-054</p>	
Multi-Position Tray with Bowl		Tray Elbow Blocks	
<p>Multi-position tray with included bowl which can be used for encouraging tactile play and hands to mid-line function. Comes with tray overlay when bowl use not required.</p> <p>5242-7113-054</p>		<p>Aid in supporting the arms and promoting a functional arm and hand position – only usable in conjunction with Tray Extension Wings.</p> <p>10cm W x 14cm H</p> <p>5356-7351-000</p>	
Tray Padded Wings		Full Tray Padding	
<p>Recommended for use in supine and just as effective in prone. Padded extension wings for comfort and protection for the elbow.</p> <p>15 cm W x 15cm L (3cm Padding)</p> <p>5242-7116-000</p>		<p>Pads the entire tray surface to give a soft surface to prevent any discomfort caused by sudden uncontrolled arm and hand movement.</p> <p>(3cm Padding)</p> <p>5242-7115-000</p>	
Tray Extension Wings			
<p>Recommended for use in supine and just as effective in prone. Extension wings offer a support surface for the child's elbows and prevents arms from slipping off the tray.</p> <p>17cm W x 16cm L</p> <p>5242-7355-000</p>			
Waistcoats			
<p>4 point detachable harness with zip. Used to provide additional trunk stability.</p>	<p>Size 2</p> <p>A: 5cm B: 32cm C: 14cm D: 15.5cm</p> <p>5242-6252-000</p>	<p>Size 3</p> <p>A: 6.5cm B: 37cm C: 17.5cm D: 18cm</p> <p>5242-6253-000</p>	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Standz Size 2_20241210