

Kidsflex




KIDSFLEX CAR SEAT

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Basic Packages


Size 1	Size 2	Size 3
<p>Includes:</p> <ul style="list-style-type: none"> • Colour • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20°) • Slim left turn swivel base • Seat belt guide on the left <p>*right turning available as special order</p>	<p>Includes:</p> <ul style="list-style-type: none"> • Colour • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20°) • Slim left turn swivel base • Seat belt guide on the left 	<p>Includes:</p> <ul style="list-style-type: none"> • Colour • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20°) • Slim left turn swivel base • Seat belt guide on the left
4170-0011-020	4170-0021-060	4170-0031-090

Specifications

Size	Total Height	Seat to Head	Back Height	Seat Width	Seat Depth	Lateral Height (Both Sizes)	Distance of Lateral Supports	Weight	Max User Weight
1	74 cm	67 cm	43-51 cm	31 cm	28-33 cm	26-40 cm	16-28 cm	8.2 kg	22 kg
2	80 cm	73 cm	49-57 cm	35 cm	34-39 cm	26-40 cm	19-32 cm	8.9 kg	36 kg
3	90 cm	82 cm	58-69 cm	40 cm	40-46 cm	34-48 cm	24-37 cm	12.3 kg	75 kg

ISOfix Adapter*


Clamp Bow

Select Size		To be used if ISOfix is not available.
Size 1 4170-6911-100		Size 1 4170-0261-000
Size 2/3 4170-6912-100		Size 2/3 4170-0262-000

*Only compatible with vehicles fitted with ISOFIX receivers

Foot Support

Shoulder Supports



Size	Size	1 Size
Size 1 4170-5241-000		4170-2750-000
Size 2/3 4170-5242-000		



Head Supports

Lateral Head Support (Recommended)		Extended Head Support	
Size 1 4170-4201-000		Size 1 4170-4301-000	
Size 2 4170-4202-000		Size 2 4170-4302-000	
Size 3 4170-4203-000		Size 3 4170-4303-000	

Lateral Supports

Small Lateral Supports (Recommended)		Large Lateral Supports	
Fixed 4170-3301-000		Fixed 4170-3302-000	
Swing-Away 4170-3401-000		Swing-Away 4170-3402-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Kidsflex_20240624_NZ