

BINGO Evo Stroller

Larger supportive seating system with many features, including a wide range of accessories and options to meet the needs of different children.

Modular seating & positioning stroller



Multi-directional suspension **protects the child** by absorbing ground bumps, and is effective off road or on rough terrain

Easily adjust to rear facing configuration for improved close interaction with caregiver or monitoring of safety

Stepless seat depth adjustment

Ergonomic, height-adjustable push handle

Sun canopy option available to allow parents and caregivers to still get out and about, knowing the child is protected from the sun

Adjustable seat unit
Seat depth and width, back support height and lower leg length

Note: Cover shown is an optional extra.

BINGO Switch:

At the heart of this innovative new design lies an interchangeable front end, a unique feature which eliminates the need to purchase multiple bases. Use the standard wheels for everyday journeys or quickly switch in the third wheel for off-road expeditions, delivering a safe, secure experience across all terrain.

Size 1: 2884-0001-000 **Size 2:** 2884-0002-000

Swivelling Front Wheel Attachment:

Size 1: 2884-9571-000 **Size 2:** 2884-9572-000



Stroller Packages Include:

Removable seat, matte black mobility base with swivel front black mag wheels and solid tyres, tie-downs, hip guides, basket, sun canopy/rain cover, adjustable suspension and 3 sided foot support.



TUBIE FUN COVERS

Personalise and protect your stroller with waterproof and washable covers from Tubie Fun!



ONE SEAT. TWO BASE OPTIONS

Indoor Cobra base option available. Height adjustable while the child is seated.



TILT

The special construction of the tilt mechanism allows a 45° tilt of the seat unit without compromising the overall stability.



COMPACT FOLDED SIZE



Easier to store in a car boot! The innovative folding mechanism can be released easily on the push bar.

BACK SUPPORT RECLINE

The operating elements of the tilt & back support recline are positioned centrally and allow for easy and secure handling.



COLOUR OPTIONS:



Turquoise Soft
(-349)



Berry Soft
(-347)



Black
(-346)

Size	Code	Seat Load Capacity	Overall Length	Overall Width	Seat Depth	Seat Width	Back Support Recline	Shoulder Height	Foot Support Height	Overall Weight	Indoor Base Seat Height	Indoor Base Footprint
1	2881-0015-000	50 kg	92 cm	63 cm	16-34 cm	18-30 cm	90-160°	28-46 cm	16-37 cm	15 kg	27-67.5*	61cm (W) x 71cm (L)
2	2882-0015-000			68 cm	22-40 cm	23-35 cm	90-160°	34-55 cm	16-40 cm	16.1 kg		
3	2883-0008-000			68 cm	30-45 cm	23-35 cm	90-135°	43-64 cm	16-40 cm	17.1 kg		

* Depending on size of castors

Accessories

**Occipital Head Support
Merino/Spacer Fabric**
2010-4110-099
Wipe Down
2010-4110-209



**Std Lateral Head Support,
Merino/Spacer Fabric**
Small: 2010-4210-099
Large: 2010-4211-099



**Standard Lateral Head
Support, Wipe Down**
Small: 2010-4210-209
Large: 2010-4211-209



**Ext Lateral Head Support,
Merino/Spacer Fabric**
Small: 2010-4310-099
Large: 2010-4311-099



**Extended Lateral Head
Support, Wipe Down**
Small: 2010-4310-209
Large: 2010-4311-209



Head Support Pads
1: 2880-4471-346
2: 2880-4472-346



Head Support Speakers
2880-8330-000
(Head support pads required)



Swivel Lock
2880-8440-000



Hand Strap
2005-6930-000



Medial Knee Support
1: 2880-1651-346
2: 2880-1652-346
3: 2880-1653-346




Std Lateral Trunk Support
1: 2880-3010-346
2: 2880-3020-346
3: 2880-3030-346



Spex Lateral Trunk Supports
Narrow Straight: 2880-3008-100
Broad: 2880-3008-310
Teardrop: 2880-3008-500



Cover Set - Llamas
Size 1: 2880-9561-121
Size 2: 2880-9562-121
5 Point Harness Cover:
2880-6511-121



Cover Set - Planets
Size 1: 2880-9561-036
Size 2: 2880-9562-036
5 Point Harness Cover:
2880-6511-036



Anterior Chest Support
1: 2880-6751-346
2: 2880-6752-346
3: 2880-6753-346



**Wonderseat Padded
Centrepoint Harness**
Size 0: 1409-5660-017
Size 1: 1409-6661-017



**5-Point Anterior
Chest Support**
Spex: 2006-6551-000



Hip Belt
Spex 2-Point:
1405-6220-027
Wonderseat 4-Point:
1406-5420-027



Groin Strap
1: 2880-6361-346
2-3: 2880-6362-346



Foot Support Spacer
4cm thick
1: 2880-5321-000
2: 2880-5322-000



Ankle Foot Fast
1: 2880-6101-000
2: 2880-6102-000



Therapy Tray
1: 2880-7111-000
2-3: 2880-7112-000
Infill: 2880-7180-000



Grab Rail
1: 2880-7251-346
2-3: 2880-7252-346



IV Pole
2880-8275-000
Extra Bracket
2880-8275-050




Oxygen Bottle Holder
2880-8253-000



Chassis Mount Tray
1: 2880-8058-000
2: 2880-8059-000



Sacral Support
2880-2510-060
Thigh Wedge
2006-0505-000
Trochanter Wedges
2006-0660-000



Cobra Indoor Base
1: 2880-8711-000
2&3: 2880-8721-000



Separate Seat Unit
1: 2881-0011-000
2: 2882-0012-000
3: 2883-0013-000



**Twin Bingo Stroller
Connection Mount**
2880-8630-000
Attendant Brake
2880-8452-000

