

Sharky



SPECIALISED TILT-IN-SPACE COMMUNE

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Product Details

Options Available	Size 1	Size 2
35cm Back Support Height	6081-0010-100	N / A
45cm Back Support Height	6081-0010-000	6082-0015-100
55cm Back Support Height	N / A	6082-0015-000
Seat Width (cm)	24, 27, 30	29, 32, 35
Seat Depth (cm)	27 - 36	27 - 40
Maximum Traversable Toilet Height (cm)	52	
Tilt Adjustment Angle (°)	-10 to +45	
Back Support Recline Angle (°)	72 - 120	
Leg Elevation Angle (°)	78 - 138	
Foot Support Angle (°)	80 - 100	
Arm Support Height (cm)	14.5 - 20.5	
Wheel Diameter (mm)	75	
Load Capacity (kg)	60	
Length (cm)	70	
Width (cm)	56	



Sharky Size 1



Sharky Size 2

Upgrade Kit to Sharky Size 2

Upgrade Sharky mobility base from size 1 to size 2 using this kit.

6085-8081-086



Mobility Base (required)

(Sizes are specific for each size Sharky)

Sharky mobility base with height adjustable compact folding base, -10 to 45° tilt in space and goes over toilet heights from 42 - 53cm.

For Sharky Size 1
6081-0015-000

For Sharky Size 2
6082-0010-000



Head Support

Suitable for use on either size Sharky

Sharky headrest, multi adjustable.
6085-4750-000



Lateral Supports

Both sizes suitable for use on either size Sharky

Sharky Lateral Pads.
11 x 6cm
6085-3310-000

Sharky Lateral Pads.
14 x 8cm
6085-3320-000



Hip Support (required)

Options for Size 1 Sharky

24cm

Sharky seat width reducing hip supports for seat width of 24cm.
6085-1110-000



27cm

Sharky seat width reducing hip supports for seat width of 27cm.
6085-1120-000



30cm

Sharky seat width reducing hip supports for seat width of 30cm.
6085-1130-000



Options for Size 2 Sharky

29cm

Sharky seat width reducing hip supports for seat width of 29cm.
6085-1210-000



32cm

Sharky seat width reducing hip supports for seat width of 32cm.
6085-1220-000



35cm

Sharky seat width reducing hip supports for seat width of 35cm.
6085-1230-000






Seat Options

Seat Infill Pad		Armrests	
Sharky removable seat infill pad for hygiene opening 6085-0401-000		Sharky removable and height adjustable armrests 6085-5500-000	
Medial Thigh Support (Both sizes suitable for use on either size Sharky)			
Size 1 8cm wide x 12cm deep 6085-1500-000	Size 2 10cm wide x 14cm deep 6085-1520-000		
Splash Guard		Commode Pan	
Sharky removable splash guard 6085-7675-000		Sharky removable commode pan with retainer 6085-8210-000	

Belt Options (recommended)

Pelvi-Loc (for use on either size Sharky)		Standard Hip Belt (for use on either size Sharky)	
Pelvi-Loc 3D for maximum pelvic security		Sharky hip belt with neoprene pads	
X-Small 1430-6340-024		Size 1 Pad Width 12.5cm 6085-1910-000	
Small 1430-6341-024		Size 2 Pad Width 16.5cm 6085-1920-000	

Foot Positioning

Foot Support System			
Sharky removable foot support system (These options are specific for each size Sharky).	Size 1 For Size 1 Sharky 6085-1810-000	Size 2 For Size 2 Sharky 6085-1820-000	
Calf Strap		Ankle Huggers (for use on either size Sharky)	
For Size 1 Sharky 6085-6160-000		Sharky neoprene ankle huggers	
For Size 2 Sharky 6085-6260-000		Size 1 6085-6101-000	
		Size 2 6085-6102-000	

Harnesses

Shoulder Harness (for use on either size Sharky)


Sharky neoprene shoulder harness.

Size 1

24.5cm high x 19cm wide
6085-6610-000

Size 2

31.5cm high x 21.5cm wide
6085-6620-000



Handles

Sharky Push Handles

6085-8600-000



Accessories

Toilet Adapter

Size 1

(for size 1 Sharky)
6085-0600-000

Size 2

(for size 2 Sharky)
6085-0620-000



Bath Subframe


Sharky bath subframe with seat angle adjustment. Minimum width of bath: 36cm

Size 1

(for size 1 Sharky)
6085-0500-000

Size 2

(for size 2 Sharky)
6085-0520-000



Tray Table

Sharky depth adjustable, removable tray table.
Note: Arm Supports required for use with Tray Table

Size 1

(for size 1 Sharky)
6085-7100-000

Size 2

(for size 2 Sharky)
6085-7200-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



SHARKY

DISTRIBUTED BY:

Script: Sharky_20240408