

# medifab

Shaping Better Lives



Juni Chair



Confidence,  
engagement,  
learning



Whether in the classroom, nursery or at home, better posture while seated creates a more positive experience for your child. It can help develop fine motor skills, increase confidence and the ability to engage and learn.

Juni chair comes in four colours and the contemporary design helps to promote social inclusion.

Juni's seat and backrest provide unique contouring to promote better posture, give higher levels of comfort and help your child feel more secure and stable in the chair.

It is our vibrant chair for children who require minimal levels of postural support when seated.

We understand that your child's needs change during each stage of growth and development. That's why the seat's height and depth, as well as armrest height, can easily be adjusted.



Base Package includes **seat unit** and **hip belt**.

# Features



## Colour Options

Juni comes in four colours (Red, Blue, Lime and Black) and the contemporary design helps to promote social inclusion.

## Tailor Juni to your child's needs

With exceptional adjustment of Juni, along with a great range of optional accessories, your child's Juni can be tailored to them, supporting them to achieve their goals and more!





## Contoured seat base and back support as standard

Gentle contouring provides support and gives stability for sitting, promoting engagement and focus on the task at hand.



## Low seat to floor height

Juni's low seat to floor height allows children to be sat with their feet on the floor, increasing support and allowing for more concentration on their play and learning.



## High weight rating

Juni has a maximum user weight of 30-60kg. Its sturdy frame helps to weight and prevent chair movements.

# Accessories



## Adjustable and Foldable Table

Adjustable height work surface, can be easily folded for storage, a great companion to the Juni chair.

This table comes in a beautiful 'natural' finish that complements the colour of the legs. (We don't currently offer the grey finish shown in the image).

**Size 1** 3209-7191-000      **Size 2** 3209-7192-000  
(Table height 31-45cm)      (Table height 40-60cm)



## Side Cushions

Reduce seat width using side supports. The thickness of the cushion equates to the amount reduced.

<b>For size 1 and 2 Juni:</b>	<b>For size 3 and 4 Juni:</b>
25mm 3209-0001-025	25mm 3209-0003-025
50mm 3209-0001-050	50mm 3209-0003-050
75mm 3209-0011-050	75mm 3209-0013-050



## Standard Skis

Helps to stabilise and prevent chair movements.

**Size 1** 3209-5021-000      **Size 3** 3209-5023-000  
**Size 2** 3209-5022-000      **Size 4** 3209-5024-000



## Sliding Footboard

Add footboard to mobile skis to support users legs

<b>Size 1</b> 3209-5071-000	<b>Size 3</b> 3209-5073-000
<b>Size 2</b> 3209-5072-000	<b>Size 4</b> 3209-5074-000



## Mobile Skis

Makes the Juni chair easier to move, ideal for pushing under tables!

**Size 1** 3209-5051-000      **Size 3** 3209-5053-000  
**Size 2** 3209-5052-000      **Size 4** 3209-5054-000



## Quick Release Tray

Adjustable depth tray table, providing a work surface and an upper limb support.

- Size 1 3201-7020-000
- Size 2 3201-7021-000
- Size 3 3201-7022-000
- Size 4 3201-7023-000



## Medial Thigh Support

Padded support for increased medial thigh support and seat contouring.

- For size 1,2 & 3 3209-1501-099
- For size 4 3209-1504-099



## 4 Point Harness

Colourful anterior chest support for improved posture.

- Size 1 Red 3201-6401-022
- Size 1 Lime 3201-6401-054
- Size 1 Blue 3201-6401-060
- Size 1 Black 3209-6351-099

- Size 2 Red 3201-6402-022
- Size 2 Lime 3201-6402-054
- Size 2 Blue 3201-6402-060
- Size 2 Black 3209-6352-099

- Size 3 Red 3201-6403-022
- Size 3 Lime 3201-6403-054
- Size 3 Blue 3201-6403-060
- Size 3 Black 3209-6353-099

- Size 4 Red 3201-6404-022
- Size 4 Lime 3201-6404-054
- Size 4 Blue 3201-6404-060
- Size 4 Black 3209-6354-099



## Backrest Extension

Extension to standard back support, increase back support height by 11-16cm (depending on size).

- Size 1 Red 3209-2101-022
- Size 1 Lime 3209-2101-054
- Size 1 Blue 3209-2101-060
- Size 1 Black 3209-2101-099

- Size 2 Red 3209-2102-022
- Size 2 Lime 3209-2102-054
- Size 2 Blue 3209-2102-060
- Size 2 Black 3209-2102-099

- Size 3 Red 3209-2103-022
- Size 3 Lime 3209-2103-054
- Size 3 Blue 3209-2103-060
- Size 3 Black 3209-2103-099

- Size 4 Red 3209-2104-022
- Size 4 Lime 3209-2104-054
- Size 4 Blue 3209-2104-060
- Size 4 Black 3209-2104-099

# Accessories

## Continued

### Height Adjustable Push Handle

Great for use with standard  
or mobile skis.

- Size 1 3209-8611-000
- Size 2 3209-8612-000
- Size 3 3209-8613-000
- Size 4 3209-8614-000



### Skis and Footboard

Standard skis with built in foot support,  
great for users who will push and could  
cause their chair to tip.

- Size 1 3209-5061-000
- Size 2 3209-5062-000
- Size 3 3209-5063-000
- Size 4 3209-5064-000



### Footsupport

Height and angle adjustable foot support.

- Size 1 3209-5001-000
- Size 2 3209-5002-000
- Size 3 3209-5003-000
- Size 4 3209-5004-000



### Outriggers

Optional accessory which increases  
seat to floor height by 5cm.

- Size 1 3209-2811-000
- Size 2 3209-2812-000
- Size 3 3209-2813-000





# Specifications



Code	Size 1	Size 2	Size 3	Size 4
Black	3209-1001-099	3209-1002-099	3209-1003-099	3209-1004-099
Blue	3209-1001-060	3209-1002-060	3209-1003-060	3209-1004-060
Red	3209-1001-022	3209-1002-022	3209-1003-022	3209-1004-022
Lime	3209-1001-022	3209-1002-054	3209-1003-054	3209-1004-054

Specifications	Size 1	Size 2	Size 3	Size 4
Seat Height (cm)	16-22	20-28	24-34	30-40
Seat Width (cm)	24	28	32.5	36.5
Seat Depth (cm)	18-24	22-30	26-36	32-42
Back Support Height (cm)	26	29	32.5	40
Load capacity of seat unit (kg)	30	40	50	60
Product weight (kg)	4	6	8	10



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)

DISTRIBUTED BY:

Flyer: Juni Chair\_20240719