



# Helio A7

ORDER FORM - AUSTRALIA(Unpriced)

### 1- Purchasing information

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_

Company Name:\* \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone  Email

Marked for: \_\_\_\_\_

### 3- Shipping information

Ship to:\* \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

Zip Code: \_\_\_\_\_

Phone: \_\_\_\_\_

### 2- ATP/Clinician information

ATP/Clinician:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone  Email

MODEL								AUD (taxes not included)
WC02	<input type="checkbox"/>	Ultralight Helio A7 wheelchair 7000 series aluminum (120 kg/265 lb weight capacity)						Base Frame
WC09	<input type="checkbox"/>	HD Heavy duty kit (159 kg/ 350 lb weight capacity)						Upcharge
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber						Upcharge
<b>Frame</b>								
<b>SEAT WIDTH<sup>3</sup> - cm (inches)</b>								
SW	<input type="checkbox"/>	35.56 (14") <input type="checkbox"/>	38.10 (15") <input type="checkbox"/>	40.64 (16") <input type="checkbox"/>	43.18 ( 17") <input type="checkbox"/>	45.72 ( 18") <sup>1</sup> <input type="checkbox"/>	48.26 (19") <sup>1</sup> <input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20") <sup>1</sup> <input type="checkbox"/> 53.34 (21") <sup>2</sup> <input type="checkbox"/> 55.88 (22") <sup>2</sup> <input type="checkbox"/>						Upcharge
<b>SEAT DEPTH - cm (inches)</b>								
SD	<input type="checkbox"/>	35.56 (14") <input type="checkbox"/>	38.10 (15") <input type="checkbox"/>	40.64 (16") <input type="checkbox"/>	43.18 ( 17") <input type="checkbox"/>	45.72 ( 18") <input type="checkbox"/>	48.26 (19") <input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20") <input type="checkbox"/>						Upcharge
<b>TRANSIT TIE DOWN</b>								
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)						NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Short belt 50-104 cm (20" - 41")						Upcharge
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Long belt 50-116 cm (20" - 46")						Upcharge
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>4</sup>						Upcharge

1- Available with HD option; 2- Only available with HD option [3. Width calculator](#) 4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

### A7 FRONT Seat-to-floor Height - Please select your choice

Forks / Position	3" caster			4" caster			5" caster			6" caster			7" caster			8" caster				
	Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length				
	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"		
4 in 10,2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"													
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"										
5 in 12,7 cm	P1																			
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"										
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"							
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"				
7 in 17,8 cm	P1									41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"					
	P2								41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"			
	P3								43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"	

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS							
CA01-03	<input type="checkbox"/>	Composite wheel w/ PU tire	3" x 1"			NC	
CA01-04	<input type="checkbox"/>		4" x 1"			NC	
CA01-05	<input type="checkbox"/>		5" x 1"			NC	
CA01-06	<input type="checkbox"/>		6" x 1"			NC	
CA01-07	<input type="checkbox"/>		7" x 1"			NC	
CA01-08	<input type="checkbox"/>		8" x 1"			NC	
CA10	<input type="checkbox"/>		Pneumatic casters	6" x 1 ¼"			Upcharge
CA11	<input type="checkbox"/>			8" x 1 ¼"			Upcharge
CA14	<input type="checkbox"/>	Newton UltraCaster Composite	4" x 1 ½ "			Upcharge	
CA15	<input type="checkbox"/>		5" x 1 ½ "			Upcharge	
CA16	<input type="checkbox"/>		6" x 1 ½ "			Upcharge	
CA17	<input type="checkbox"/>		4" x 1 ½ "			Upcharge	
CA18	<input type="checkbox"/>		Newton UltraCaster Aluminum	5" x 1 ½ "			Upcharge
CA19	<input type="checkbox"/>	6" x 1 ½ "				Upcharge	
CASTER OPTIONS							
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) <sup>1</sup>					Upcharge
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)					Upcharge

1 - Not Available with WC-19

### A7 - REAR Seat-to-floor Height

WHEEL SIZES	Height available		Specify desired height: _____
20" (451) wheel <input type="checkbox"/>	31.75cm/12.5 "	to 43.81cm/17.25"	
22" (501) wheel <input type="checkbox"/>	34.29cm/13.5 "	to 46.35cm/18.25"	
24" (540) wheel <input type="checkbox"/>	36.83cm/14.5 "	to 48.89cm/19.25"	
25" (559) wheel <input type="checkbox"/>	38.10cm/15 "	to 50.16cm/19.75"	
26" (590) wheel <input type="checkbox"/>	39.37cm/15.5 "	to 51.34cm/20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (½") of selected Seat to Floor height"

### Footrest

SWING-IN/SWING-OUT FOOTREST						
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at: *		NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at: *		NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest	39.29 - 52.07cm (15.5" - 20.5")	Set at: *		NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at: *		NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at: *		NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest	39.37 - 52.07cm (15.5" - 20.5")	Set at: *		NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest	34.29 - 40.64cm (13.5" - 16")	Set at: *		NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at: *		NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest	39.37 - 52.07cm (15.5" - 20.5")	Set at: *		NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) <sup>1</sup>	19.05 - 27.94cm (7.5" - 11")	Set at: *		NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension 2" <sup>1,2</sup>	24.13 - 33.02cm (9.5" - 13")	Set at: *		NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension 4" <sup>1,2</sup>	24.13 - 36.83cm (9.5" - 14.5")	Set at: *		NC
FOOTREST OPTIONS						
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) <sup>3</sup>	12.70 - 34.29cm (5" - 13.5")	Set at: *		Upcharge
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, length adjustable <sup>4</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	35.56 - 55.88 cm (14"-22")	Set at:		Upcharge
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount;

\* +/- 1.27 cm (1/2")

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite <sup>1</sup>				NC
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite				Upcharge
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum				Upcharge
FP04	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber with rubber inserts				Upcharge
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum				Upcharge
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				Upcharge
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop (not available with one piece footplate)				Upcharge
FP09	<input type="checkbox"/>	Calf strap				Upcharge
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				Upcharge

1- Not available with front mount

SWING AWAY RESIDUAL LIMB SUPPORT							
FR27	<input type="checkbox"/>	Residual limb support - flat 15.24x20.32cm (6"x8")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR28	<input type="checkbox"/>	Residual limb support - flat 20.32x24.4cm (8"x10")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR29	<input type="checkbox"/>	Residual limb support - flat 25.4x25.4cm (10"x10")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR30	<input type="checkbox"/>	Residual limb support - flat 25.4x35.56cm (10"x14")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR32	<input type="checkbox"/>	Residual limb support - contour 15.24x20.32cm (6"x8")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR33	<input type="checkbox"/>	Residual limb support - contour 20.32x24.4cm (8"x10")		Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge

### Backrest

FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE (Not available with push handles options)							
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	40.64cm (16")				NC
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	45.72cm (18")				NC
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	50.08cm (20")				NC
BC59	<input type="checkbox"/>	8° bend aluminum back cane	53.34 - 60.96cm (21" - 24")	increment of 1"	Set at <sup>2</sup> :		Upcharge

FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)							
BC11-912	<input type="checkbox"/>	Straight aluminum back cane	22.86 - 30.48cm (9" - 12")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC11-1215	<input type="checkbox"/>	Straight aluminum back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC11-1518	<input type="checkbox"/>	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC11-1821	<input type="checkbox"/>	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC12-1315	<input type="checkbox"/>	8° bend aluminum back cane	33.02 - 38.10cm (13.5" - 15")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC12-1518	<input type="checkbox"/>	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		NC
BC12-1821	<input type="checkbox"/>	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		NC

ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)								
BC13-1518	<input type="checkbox"/>	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC13-1821	<input type="checkbox"/>	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC14-1518	<input type="checkbox"/>	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC14-1821	<input type="checkbox"/>	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
		Select back cane angle:	85° <sup>1</sup> <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>

PUSH HANDLES OPTIONS (ONE option at a time)							
PH01	<input type="checkbox"/>	Push handles					STD
PH05	<input type="checkbox"/>	Fold down push handles <sup>3,5</sup>					Upcharge
PH06-920	<input type="checkbox"/>	clamp on height adjustable push handles 22.86 - 50.8cm (9" - 20") <sup>3,4,5</sup>					Upcharge
PH06-923	<input type="checkbox"/>	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") <sup>3,4,5</sup>					Upcharge
PH02	<input type="checkbox"/>	Omit push handles					NC

STROLLER HANDLE AND STABILIZER BAR (only available with push handle)							
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)					Upcharge
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)					Upcharge

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit; 4- Not compatible with swing away armrests, 8° bend back canes, 21-24" back canes and amputee axle plates; 5- Not compatible with stroller handle, the stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY							
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery					NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery					Upcharge
BH99	<input type="checkbox"/>	Omit Back upholstery					STD

SEAT CUSHION							
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion					Upcharge
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion					Upcharge
SC99	<input type="checkbox"/>	Omit seat cushion					STD

## Rear wheels

ONE ARM DRIVE <sup>1</sup> <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber</i>							Available	
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			Upcharge
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			Upcharge
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)						NC
REAR WHEELS							Available	
RW01	<input type="checkbox"/>	Mag <sup>4</sup>				20", 22" 24"		Upcharge
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>5</sup>				20", 22" 24"		NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		Upcharge
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	Upcharge
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>
RW07	<input type="checkbox"/>	Spinergy CLX <sup>7</sup>	Blk STD <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>
								Pink <sup>6</sup> <input type="checkbox"/>
								Grn <sup>6</sup> <input type="checkbox"/>
								Org <sup>6</sup> <input type="checkbox"/>
								24", 25"
								Upcharge
TIRES							Available	
TI01	<input type="checkbox"/>	Soft Urethane - 3.49cm (1 3/8")				20", 22" 24"		STD
TI02	<input type="checkbox"/>	Soft Urethane - knobby 3.49cm (1 3/8")				22", 24", 25", 26"		Upcharge
TI04	<input type="checkbox"/>	Soft Urethane - Shox 2.54cm (1")				22", 24", 25", 26"		Upcharge
TI05	<input type="checkbox"/>	Pneumatic - 3.49cm (1 3/8")				20", 22", 24", 25" 26"		NC
TI08	<input type="checkbox"/>	Pneumatic with Airless insert - 3.49cm (1 3/8") <sup>8</sup>				22", 24", 26"		Upcharge
TI06	<input type="checkbox"/>	Pneumatic HP - 2.54cm (1")				20", 22", 24", 25", 26"		Upcharge
TI07	<input type="checkbox"/>	Pneumatic HP - Puncture resistant, black 2.54cm (1")				24", 25", 26"		Upcharge

- 1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option  
 4- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires; 6- Colour spokes require an additional 10 days lead-time  
 7- Only available with 1" tires; 8- Only available with Mag, Newton One and Newton Gravity rear wheels

## Handrims

HANDRIM							Available	
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>				20", 22", 24", 25"		Upcharge
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>				20", 22", 24", 25", 26"		Upcharge
HR03	<input type="checkbox"/>	High friction coated <sup>1,2</sup>				22", 24", 25", 26"		Upcharge
HR07 to HR18	<input type="checkbox"/>	Newton Air Grip <sup>2</sup>	<i>See page 6 for the color selection and price</i>			20", 22", 24", 25", 26"		
HR28	<input type="checkbox"/>	Surge <sup>1,2</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		Upcharge
HR05	<input type="checkbox"/>	Surge LT <sup>1,2</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		Upcharge
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2</sup>				22", 24", 25", 26"		Upcharge
HR62	<input type="checkbox"/>	Tour H <sup>1,2</sup> (1 5/16" Oval handrim with high friction strip)				22", 24", 25", 26"		Upcharge
HR63	<input type="checkbox"/>	Nova H <sup>1,2</sup> (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"		Upcharge
HR64	<input type="checkbox"/>	Optimum H <sup>1,2</sup> (1 1/8" heart-shaped handrim)				22", 24", 25", 26"		Upcharge
HR99	<input type="checkbox"/>	Omit handrim						NC

- 1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

## Handrims option

ASSEMBLY POSITION							Available	
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab					STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab					NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup>	Short tab 2 positions					Upcharge
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup>	Short tab 2 positions STD <input type="checkbox"/>	Short tab 1 position <input type="checkbox"/>				Upcharge
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							Available	
HR23	<input type="checkbox"/>	Standard grip				22", 24", 25", 26"		Upcharge
HR27	<input type="checkbox"/>	Super grip				22", 24", 25", 26"		Upcharge
HR099	<input type="checkbox"/>	Omit thumb piece						NC


- 3- Only available with Aluminum and Air Grip handrims;

AXLE							Available	
AX01	<input type="checkbox"/>	Quick release axle <sup>1</sup>						STD
AX02	<input type="checkbox"/>	Permanent axle						NC
AX03	<input type="checkbox"/>	Quad release axle <sup>1</sup>						Upcharge
AXLE PLATE							Available	
AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>						\$196

- 1- Not required with one arm drive 2- May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER							Available	
CM00	<input type="checkbox"/>	0° camber						STD
CM03	<input type="checkbox"/>	3° camber						NC
CM06	<input type="checkbox"/>	6° camber (not available with amputee axle plate)						Upcharge

**CENTER OF GRAVITY**



**Choose position of rear wheel axle plate:**

Measure distance X from front of back post to center of axle receiver.

Most stable						Most responsive					
<b>2,5</b>	3,2	3,8	4,4	<b>5,1</b>	5,7	6,3	7	<b>7,6</b>	8,3	8,9	9,5
<b>1"</b>	1/4"	1/2"	3/4"	<b>2"</b>	1/4"	1/2"	3/4"	<b>3"</b>	1/4"	1/2"	3/4"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**If not indicated, default setting is at 3,8 cm (1.5")**

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6.3 cm (2.5").

Possible conflict between a 26" rear wheel and a center of gravity over 6.3 cm (2.5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

**WHEEL LOCK**

WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock			STD
WL16		Aluminum push to lock with extension			Upcharge
WL02	<input type="checkbox"/>	Pull to lock			Upcharge
WL04	<input type="checkbox"/>	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			Upcharge
WL13 WL14	<input type="checkbox"/>	6" removable extension handle			Upcharge
WL05	<input type="checkbox"/>	Scissors Lock			Upcharge
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)			Upcharge
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>1</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>	Upcharge
WL99	<input type="checkbox"/>	Omit wheel locks			NC

1- For more information about incompatibilities and restriction, refer to the table at page 8 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_intl/support-and-education/how-to-documents](https://www.motioncomposites.com/en_intl/support-and-education/how-to-documents)

**Armrests**

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) <sup>1,4</sup>						
AR01-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	Upcharge
AR01-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	Upcharge
AR02-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	Upcharge
AR02-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	Upcharge
HEIGHT ADJUSTABLE "T" ARMREST <sup>1</sup>						
AR13-610	<input type="checkbox"/>	25.4cm (10") Desk length armpad <sup>5</sup>	15.25 - 25.40cm (6" - 10")			Upcharge
AR03-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")			Upcharge
AR03-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")			Upcharge
AR14-610	<input type="checkbox"/>	35.56cm (14") Full length armpad <sup>5</sup>	15.25 - 25.40cm (6" - 10")			Upcharge
AR04-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")			Upcharge
AR04-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")			Upcharge
HEIGHT ADJUSTABLE SWING AWAY ARMREST <sup>2,4</sup>						
AR05	<input type="checkbox"/>	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")			Upcharge
AR23	<input type="checkbox"/>	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")			Upcharge
AR06	<input type="checkbox"/>	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")			Upcharge
AR17	<input type="checkbox"/>	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")			Upcharge
ARMREST OPTIONS						
AR12	<input type="checkbox"/>	35.56cm (14") Full Length Armpad - Gel Ovation (pair) <sup>4</sup>				Upcharge
AR99	<input type="checkbox"/>	Omit armrest				NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; 3- Not available with swing away and flip back armrest;

4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
SG01	<input type="checkbox"/>	Plastic side guards (removable)				Upcharge
SG02	<input type="checkbox"/>	Carbon fiber side guards (fixed)				Upcharge
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)				Upcharge
SG06	<input type="checkbox"/>	Plastic fender side guards - straight bracket (removable)				Upcharge
SG08	<input type="checkbox"/>	Plastic fender side guards - offset bracket (removable) <sup>1</sup>				Upcharge
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)				Upcharge
SG99	<input type="checkbox"/>	Omit side guard				NC
SIDE GUARDS OPTION						
SG07	<input type="checkbox"/>	Side guards clip (pair) <sup>2</sup>				NC

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

POSITIONING BELT			
SB01-S	<input type="checkbox"/>	Velcro adjustable - <b>Standard 88.9cm (35")</b>	Upcharge
SB01-L	<input type="checkbox"/>	Velcro adjustable - <b>Long 152.40cm (60")</b>	Upcharge
SB02-S	<input type="checkbox"/>	Auto buckle - <b>Standard 88.9cm (35")</b>	Upcharge
SB02-L	<input type="checkbox"/>	Auto buckle - <b>Long 152.40cm (60")</b>	Upcharge
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Long 154.94cm (61")</b>	Upcharge
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Small 129.54cm (51")</b>	Upcharge
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Long 154.94cm (61")</b>	Upcharge
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Small 129.54cm (51")</b>	Upcharge
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Long 154.94cm (61")</b>	Upcharge
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	Upcharge
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Medium 73.66 (18")</b>	Upcharge
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Long 93.98cm (24")</b>	Upcharge
SB99	<input type="checkbox"/>	Omit positioning belt	STD

HEADREST <sup>1</sup>			
AC01-S	<input type="checkbox"/>	Profiled Headrest - <b>Small 8.8x15.2cm (3.5"x6")</b>	Upcharge
AC01-M	<input type="checkbox"/>	Profiled Headrest - <b>Medium 11.4x17.8cm (4.5"x7")</b>	Upcharge
AC01-L	<input type="checkbox"/>	Profiled Headrest - <b>Large 14x2.3cm(5.5"x8")</b>	Upcharge
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar <sup>1</sup>	Upcharge

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	Upcharge
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	Upcharge
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>	Upcharge
AC06	<input type="checkbox"/>	Cane and crutch holder      Left <input type="checkbox"/> Right <input type="checkbox"/>	Upcharge
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - <i>Installed to the left</i>	Upcharge
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup> Available: 20", 22", 24", 26"	Upcharge
AC09	<input type="checkbox"/>	Footrest neoprene impact guard ( <i>pair</i> ) - Black	Upcharge
AC10	<input type="checkbox"/>	Seat pouch	Upcharge
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	Upcharge
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	Upcharge
AC23	<input type="checkbox"/>	Newton Backpack	Upcharge
AC12	<input type="checkbox"/>	Velcro Closing Strap	Upcharge
AC13	<input type="checkbox"/>	Tool kit	Upcharge
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt      W/S <input type="checkbox"/> W/M <input type="checkbox"/> W/L <input type="checkbox"/> M/S <input type="checkbox"/> M/M <input type="checkbox"/> M/L <input type="checkbox"/> M/XL <input type="checkbox"/> M/XXL <input type="checkbox"/> M/XXXL <input type="checkbox"/>	Upcharge

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4")


2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerdy; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down

Colors																				
Options <sup>1</sup>	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Upcharge
Frame <sup>2</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	Custom Color <a href="#">Specify PPG Paintit color:</a>																	additional 10 days lead-time		Upcharge

1- Finish between Airgrip and frame might be different; 2- with black hydroformed aluminum crossbrace

### Fork and Stem bolt combinations



### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :