



# Helio C2

ORDER FORM - AUSTRALIA(UNPRICED)

### 1- Purchasing information

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_

Company Name:\* \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone  Email

Marked for: \_\_\_\_\_

### 3- Shipping information

Ship to:\* \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

Zip Code: \_\_\_\_\_

Phone: \_\_\_\_\_

### 2- ATP/Clinician information

ATP/Clinician:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone  Email  \* Required

MODEL								AUD (Taxes not included)
WC05	<input type="checkbox"/>	UltraLight Folding Carbon Fiber Helio C2 Wheelchair (120 kg - 265 lb weight capacity)						BASE FRAME
WC09	<input type="checkbox"/>	HD Heavy Duty package (159 kg - 350 lb weight capacity)						Upcharge
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber						Upcharge
Frame								
SEAT WIDTH <sup>3</sup> - cm (inches)								
SW	<input type="checkbox"/>	35.56 (14")	38.10 (15")	40.64 (16")	43.18 (17")	45.72 (18") <sup>1</sup>	48.26 (19") <sup>1</sup>	NC
	<input type="checkbox"/>			50.80 (20") <sup>1</sup>	53.34 (21") <sup>2</sup>	55.88 (22") <sup>2</sup>		Upcharge
SEAT DEPTH								
SD	<input type="checkbox"/>	35.56 (14")	38.10 (15")	40.64 (16")	43.18 (17")	45.72 (18")	48.26 (19")	NC
	<input type="checkbox"/>			50.80 (20")				Upcharge
TRANSIT TIE DOWN								
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)						NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup>			Short belt 50-104 cm (20" - 41")			Upcharge
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup>			Long belt 50-116 cm (20" - 46")			Upcharge
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>4</sup>						Upcharge

1- Available with HD option; 2- Only available with HD option [3. Width calculator](#) 4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

NEWTON ULTRALIGHT PACKAGE			
PK01	<input type="checkbox"/>	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels	Save \$260 Upcharge

### C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position	7.62cm (3") caster			10.16cm (4") caster			12.7cm (5") caster			15.24cm (6") caster			17.78cm (7") caster			20.32cm (8") caster			
	Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			
	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	
4" 10.2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"												
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
5" 12.7 cm	P1																		
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"						
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"			
7" 17.8 cm	P1										41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"			
	P2						41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"				
	P3						43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"

Height available +- 0.65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS							
CA01-03	<input type="checkbox"/>	Composite wheel w/ PU tire	3" x 1"			NC	
CA01-04	<input type="checkbox"/>		4" x 1"			NC	
CA01-05	<input type="checkbox"/>		5" x 1"			NC	
CA01-06	<input type="checkbox"/>		6" x 1"			NC	
CA01-07	<input type="checkbox"/>		7" x 1"			NC	
CA01-08	<input type="checkbox"/>		8" x 1"			NC	
CA10	<input type="checkbox"/>		Pneumatic casters	6" x 1 ¼"			Upcharge
CA11	<input type="checkbox"/>			8" x 1 ¼"			Upcharge
CA14	<input type="checkbox"/>	Newton UltraCaster Composite	4" x 1 ½"			Upcharge	
CA15	<input type="checkbox"/>		5" x 1 ½"			Upcharge	
CA16	<input type="checkbox"/>		6" x 1 ½"			Upcharge	
CA17	<input type="checkbox"/>		4" x 1 ½"			Upcharge	
CA18	<input type="checkbox"/>		Newton UltraCaster Aluminum	5" x 1 ½"			Upcharge
CA19	<input type="checkbox"/>	6" x 1 ½"				Upcharge	
CASTER OPTIONS							
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) <sup>1</sup>					Upcharge
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)					Upcharge

1- Not available with WC-19

### C2 - REAR Seat-to-floor Height

WHEEL SIZES	Height available		Specify desired height: _____
20" (451) wheel <input type="checkbox"/>	31.75cm / 12.5"	to 43.81cm / 17.25"	
22" (501) wheel <input type="checkbox"/>	34.29cm / 13.5"	to 46.35cm / 18.25"	
24" (540) wheel <input type="checkbox"/>	36.83cm / 14.5"	to 48.89cm / 19.25"	
25" (559) wheel <input type="checkbox"/>	38.10cm / 15"	to 50.16cm / 19.75"	
26" (590) wheel <input type="checkbox"/>	39.37cm / 15.5"	to 51.34cm / 20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1cm" of selected seat to floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")			Set at: *	NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *	NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 39.29 - 52.07cm (15.5" - 20.5")			Set at: *	NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")			Set at: *	NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *	NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")			Set at: *	NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 34.29 - 40.64cm (13.5" - 16")			Set at: *	NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *	NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")			Set at: *	NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) <sup>1</sup>	19.05 - 27.94cm (7.5" - 11")		Set at: *	NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension 2" <sup>1,2</sup>	24.13 - 33.02cm (9.5" - 13")		Set at: *	NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension 4" <sup>1,2</sup>	24.13 - 36.83cm (9.5" - 14.5")		Set at: *	NC
FOOTREST OPTIONS <sup>3</sup>						
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) <sup>3</sup>	12.70 - 34.29cm (5" - 13.5")		Set at: *	Upcharge
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, length adjustable <sup>4</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	35.56 - 55.88 cm (14"-22")		Set at:	Upcharge
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- Available only with two piece adjust. footplates; 2- May interfere with front casters; 3- Not available with front mount;

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite <sup>1</sup>				NC
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite				Upcharge
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum				Upcharge
FP04	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber with rubber inserts				Upcharge
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum				Upcharge
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				Upcharge
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop (not available with one piece footplate)				Upcharge
FP09	<input type="checkbox"/>	Calf strap				Upcharge
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				Upcharge

1- Not available with front mount

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 15.24x20.32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR28	<input type="checkbox"/>	Residual limb support - flat 20.32x24.4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR29	<input type="checkbox"/>	Residual limb support - flat 25.4x25.4cm (10"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR30	<input type="checkbox"/>	Residual limb support - flat 25.4x35.56cm (10"x14")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR32	<input type="checkbox"/>	Residual limb support - contour 15.24x20.32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge
FR33	<input type="checkbox"/>	Residual limb support - contour 20.32x24.4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>		Upcharge

Backrest								
<b>FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)</b>								
BC06-912		Straight Carbon fiber back cane	22.86 - 30.48cm (9" - 12")	<b>*Not available until Mid-2023</b>			NC	
BC06-1215		Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	<b>*Not available until Mid-2023</b>			NC	
BC06-1518		Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	<b>*Not available until Mid-2023</b>			NC	
BC06-1821		Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	<b>*Not available until Mid-2023</b>			NC	
BC08-1315		8° Bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	<b>*Not available until Mid-2023</b>			Upcharge	
BC08-1518		8° Bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	<b>*Not available until Mid-2023</b>			Upcharge	
BC08-1821		8° Bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	<b>*Not available until Mid-2023</b>			Upcharge	
<b>* Replacement options below: Aluminum fixed angle and height adjustable or straight carbon fiber angle and height adjustable</b>								
<b>FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)</b>								
BC11-912	<input type="checkbox"/>	Straight Aluminum back cane	22.86 - 30.48cm (9" - 12")	increment of 1/2"	Set at <sup>2</sup> :		NC	
BC11-1215	<input type="checkbox"/>	Straight Aluminum back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at <sup>2</sup> :		NC	
BC11-1518	<input type="checkbox"/>	Straight Aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		NC	
BC11-1821	<input type="checkbox"/>	Straight Aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		NC	
BC12-1315	<input type="checkbox"/>	8° Bend Aluminum back cane	34.29 - 38.1cm (13.5" - 15")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC12-1518	<input type="checkbox"/>	8° Bend Aluminum back cane	38.1 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC12-1821	<input type="checkbox"/>	8° Bend Aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
<b>ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)</b>								
BC09-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC09-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC09-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at <sup>2</sup> :		Upcharge	
BC10-1315		8° bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	<b>*Not available until Mid-2023</b>			Upcharge	
BC10-1518		8° bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	<b>*Not available until Mid-2023</b>			Upcharge	
BC10-1821		8° bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	<b>*Not available until Mid-2023</b>			Upcharge	
		Select back cane angle:	85° <sup>1</sup> <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>
<b>PUSH HANDLES OPTIONS (ONE option at a time)</b>								
PH01	<input type="checkbox"/>	Push handles					STD	
PH05	<input type="checkbox"/>	Fold down push handles <sup>3,5</sup>					Upcharge	
PH06-920	<input type="checkbox"/>	clamp on height adjustable push handles 22.86 - 50.8cm (9" - 20") <sup>3,4,5</sup>					Upcharge	
PH06-923	<input type="checkbox"/>	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") <sup>3,4,5</sup>					Upcharge	
PH02	<input type="checkbox"/>	Omit push handles					NC	
<b>STROLLER HANDLE AND STABILIZER BAR (only available with push handle)</b>								
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)					Upcharge	
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)					Upcharge	
<sup>1</sup> - Not available with 8 degree bend cane; <sup>2</sup> - If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; <sup>3</sup> - 120kg (265 lb) weight limit <sup>4</sup> - Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; <sup>5</sup> - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock								
<b>BACK UPHOLSTERY</b>								
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery					NC	
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery					Upcharge	
BH99	<input type="checkbox"/>	Omit Back upholstery					STD	
<b>SEAT CUSHION</b>								
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion					Upcharge	
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion					Upcharge	
SC99	<input type="checkbox"/>	Omit seat cushion					STD	

Rear wheels														
<b>ONE ARM DRIVE</b> <sup>1</sup> <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber</i>														
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			Upcharge						
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			Upcharge						
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)						NC						
<b>REAR WHEELS</b>														
Available														
RW01	<input type="checkbox"/>	Mag <sup>4</sup>				20", 22" 24"		Upcharge						
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>5</sup>				20", 22" 24"		NC						
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		Upcharge						
RW04	<input type="checkbox"/>	Spineray Spox	Blk STD <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	22", 24", 25", 26"						
RW05	<input type="checkbox"/>	Spineray LX	Blk STD <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	Pink <sup>6</sup> <input type="checkbox"/>	Grn <sup>6</sup> <input type="checkbox"/>	Org <sup>6</sup> <input type="checkbox"/>	22", 24", 25", 26"		Upcharge
RW07	<input type="checkbox"/>	Spineray CLX <sup>7</sup>	Blk STD <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	Pink <sup>6</sup> <input type="checkbox"/>	Grn <sup>6</sup> <input type="checkbox"/>	Org <sup>6</sup> <input type="checkbox"/>	24", 25"		Upcharge
<b>TIRES</b>														
Available														
TI01	<input type="checkbox"/>	Soft Urethane - 3.49cm (1 3/8")				20", 22" 24"		STD						
TI02	<input type="checkbox"/>	Soft Urethane - knobby 3.49cm (1 3/8")				22", 24", 25", 26"		Upcharge						
TI04	<input type="checkbox"/>	Soft Urethane - Shox 2.54cm (1")				22", 24", 25", 26"		Upcharge						
TI05	<input type="checkbox"/>	Pneumatic - 3.49cm (1 3/8")				20", 22", 24", 25", 26"		NC						
TI08	<input type="checkbox"/>	Pneumatic with Airless insert - 3.49cm (1 3/8") <sup>8</sup>				22", 24", 26"		Upcharge						
TI06	<input type="checkbox"/>	Pneumatic HP - 2.54cm (1")				20", 22", 24", 25", 26"		Upcharge						
TI07	<input type="checkbox"/>	Pneumatic HP - Puncture resistant, black 2.54cm (1")				24", 25", 26"		Upcharge						

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option

4- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires; 6- Colour spokes require an additional 10 days lead-time

7- Only available with 1" tires; 8- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims								
<b>HANDRIM</b>								
Available								
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>				20", 22", 24", 25"		Upcharge
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>				20", 22", 24", 25", 26"		Upcharge
HR03	<input type="checkbox"/>	High friction coated <sup>1,2</sup>				22", 24", 25", 26"		Upcharge
HR07 to HR18	<input type="checkbox"/>	Newton Air Grip <sup>2</sup>	<i>See on page 6 for color choice and price</i>			20", 22", 24", 25", 26"		
HR28	<input type="checkbox"/>	Surge <sup>1,2</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)				**Limited availability	22", 24", 25", 26"	Upcharge
HR05	<input type="checkbox"/>	Surge LT <sup>1,2</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)				**Limited availability	22", 24", 25", 26"	Upcharge
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2</sup>				**Limited availability	22", 24", 25", 26"	Upcharge
HR62	<input type="checkbox"/>	Tour H <sup>1,2</sup> (1 5/16" Oval handrim with high friction strip)					22", 24", 25", 26"	Upcharge
HR63	<input type="checkbox"/>	Nova H <sup>1,2</sup> (1 1/4" Oval handrim with high friction strip)					22", 24", 25", 26"	Upcharge
HR64	<input type="checkbox"/>	Optimum H <sup>1,2</sup> (1 1/8" heart-shaped handrim)					22", 24", 25", 26"	Upcharge
HR99	<input type="checkbox"/>	Omit handrim						NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option								
<b>ASSEMBLY POSITION</b>								
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab					STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab					NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup>	Short tab 2 positions					Upcharge
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup>	Short tab 2 positions STD <input type="checkbox"/>	Short tab 1 position <input type="checkbox"/>				Upcharge
<b>THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS</b>								
HR23	<input type="checkbox"/>	Standard grip				**Limited availability		Upcharge
HR27	<input type="checkbox"/>	Super grip				**Limited availability		Upcharge
HR099	<input type="checkbox"/>	Omit thumb piece						NC

3- Only available with Aluminum and Air Grip handrims;

AXLE								
<b>AXLE</b> <sup>1</sup>								
AX01	<input type="checkbox"/>	Quick release axle <sup>1</sup>						STD
AX02	<input type="checkbox"/>	Permanent axle						NC
AX03	<input type="checkbox"/>	Quad release axle <sup>1</sup>						Upcharge
<b>AXLE PLATE</b>								
AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>						Upcharge

1- Not required with one arm drive 2- May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER								
CM00	<input type="checkbox"/>	0° camber						STD
CM03	<input type="checkbox"/>	3° camber						NC
CM06	<input type="checkbox"/>	6° camber (not available with amputee axle plate)						Upcharge



POSITIONING BELT			
SB01-S	<input type="checkbox"/>	Velcro adjustable - <b>Standard 88.9cm (35")</b>	Upcharge
SB01-L	<input type="checkbox"/>	Velcro adjustable - <b>Long 152.40cm (60")</b>	Upcharge
SB02-S	<input type="checkbox"/>	Auto buckle - <b>Standard 88.9cm (35")</b>	Upcharge
SB02-L	<input type="checkbox"/>	Auto buckle - <b>Long 152.40cm (60")</b>	Upcharge
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Long 154.94cm (61")</b>	Upcharge
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Small 129.54cm (51")</b>	Upcharge
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Long 154.94cm (61")</b>	Upcharge
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Small 129.54cm (51")</b>	Upcharge
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Medium 142.24cm (56")</b>	Upcharge
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Long 154.94cm (61")</b>	Upcharge
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	Upcharge
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Medium 73.66 (18")</b>	Upcharge
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Long 93.98cm (24")</b>	Upcharge
SB99	<input type="checkbox"/>	Omit positioning belt	STD

HEADREST <sup>1</sup>			
AC01-S	<input type="checkbox"/>	Profiled Headrest - <b>Small 8.8x15.2cm (3.5"X6")</b>	Upcharge
AC01-M	<input type="checkbox"/>	Profiled Headrest - <b>Medium 11.4x17.8cm (4.5"X7")</b>	Upcharge
AC01-L	<input type="checkbox"/>	Profiled Headrest - <b>Large 14x2.3cm(5.5"X8")</b>	Upcharge
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar <sup>1</sup>	Upcharge

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	Upcharge
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	Upcharge
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>	Upcharge
AC06	<input type="checkbox"/>	Cane and crutch holder      Left <input type="checkbox"/> Right <input type="checkbox"/>	Upcharge
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - installed on the left	Upcharge
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup> Available: 20", 22", 24", 26"	Upcharge
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black	Upcharge
AC10	<input type="checkbox"/>	Seat pouch	Upcharge
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	Upcharge
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	Upcharge
AC23	<input type="checkbox"/>	Newton Backpack	Upcharge
AC12	<input type="checkbox"/>	Velcro Closing Strap	Upcharge
AC13	<input type="checkbox"/>	Tool kit	Upcharge
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt      W/S <input type="checkbox"/> W/M <input type="checkbox"/> W/L <input type="checkbox"/> M/S <input type="checkbox"/> M/M <input type="checkbox"/> M/L <input type="checkbox"/> M/XL <input type="checkbox"/> M/XXL <input type="checkbox"/> M/XXXL <input type="checkbox"/>	Upcharge

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4")


2 - Cannot be combined with Anti Tipplers

3 - Not available with Mag wheels and Spinegy; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																					
Options <sup>1</sup>	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Upcharge
Frame <sup>2</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	<input type="checkbox"/>	Custom Color	<a href="#">Specify PPG Paintit color:</a>																additional 10 days lead-time	Upcharge	

1- Finish between Airgrip and frame might be different 2- With carbon fiber crossbrace;

### Fork and Stem bolt combinations



### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :