# spex for BINGO

### FOR HOGGI STROLLERS

NEW ZEALAND SCRIPTING AND ORDERING GUIDE

## Client Details

Client/User name:


Date:


DOB:

## Measurements

Select Size		Size 1		Size 2	
Order Code		2885-0031-032		2885-0032-032	
	A Hip width	17-28cm (6.5-11")		22-33cm (8.5-13")	
	<b>B</b> Seat depth	18-30cm (7-12")		23-38cm (9-15")	
	<b>C</b> Back height	36-44cm (14-17")		40-48cm (16-19")	
	<b>D</b> Lower leg length*	25-40cm (10-16")		25-45cm (10-18")	
	E Trunk lateral width	12-28cm (5-11")		15-30cm (6-12")	
	<b>F</b> Trunk lateral height	15-34cm (6-13.5")		17-40cm (6.5-16")	
	<b>G</b> Arm support height	15-22cm (6-8.5")		15-22cm (6-8.5")	
	H Internal arm support width	20-32cm (8-13")		28-38cm (11-15")	
	I Head support height	45-65cm (18-25.5")		50-70cm (20-28")	
Seating	J Cushion total width	33cm (13") Size 1.1	33cm (13") Size 1.2	38cm (15") Size 2.1	38cm (15") Size 2.2
specs	K Cushion total depth	33cm		38cm	
	L Internal cushion width	23cm (9") Size 1.1	28cm (11") Size 1.2	28cm (11") Size 2.1	33cm (13") Size 2.2
	M Centre of thigh contours	15cm (6") Size 1.1	17cm (6.5") Size 1.2	18cm (7") Size 2.1	20cm (8") Size 2.2
	N Thigh contour depth	10cm (4") Size 1.1	12cm (5") Size 1.2	12cm (5") Size 2.1	15cm (6") Size 2.2
	O Cushion side height	14cm (5.5")		14cm (5.5")	
	P Thigh contour height	9cm (3.5")		9cm (3.5")	
	<b>Q</b> Abductor contour height	5cm (2")		5cm (2")	
	<b>R</b> Back cushion width	33cm (13")		38cm (15")	
	<b>S</b> Trunk cushion width	26cm (10")		30cm (12")	
	Frame seat depth capacity	18-30cm (7-12")		20-38cm (8-15")	
	Frame back height capacity	36-46cm (14-18")		41-55cm (16-22")	
	Backrest recline angle	90-160°			
Stroller	Height of push bar		75-110c	cm (29.5-43")	
specs	Overall frame width	63cm (25")		68cm (27")	
	Max load of seat unit		50kg (110lb) (body w	eight including other load)	
	Base frame capacity	60kg (132lb)			

Notes





### The BINGO Spex package includes:

- Matte black frame with tie downs
- 3kg fabric basket
- Tray table
- Sun canopy



FABRIC BASKET



TRAY TABLE



## Colour Selection

(Express delivery on colours stocked in our Devonport Warehouse).

Jet	Royal	Chilli	Mango	Bumblebee
		Express Delivery		
-099	-077	-022	-030	-044
Ocean	Leaf	Calypso	Granite	Ink / Wipe Down

. Express Delivery



Calypso





## Cushion Selection

Size 1.1	Size 1.2	Size 2.1	Size 2.2
Fits BINGO Size 1		Fits BINGO Size 2	
33cm (13") total width	33cm (13") total width	38cm (15") total width	38cm (15") total width
23cm (9") internal width	28cm (11") internal width	28cm (11") internal width	33cm (13") internal width
(Recommended for	(Recommended for	(Recommended for	(Recommended for
18-23cm seat depth)	23-30cm seat depth)	23-33cm seat depth)	28-38cm seat depth)
2885-0610-XXX	2885-0611-XXX	2885-0612-XXX	2885-0613-XXX
Spare Cover:	Spare Cover:	Spare Cover:	Spare Cover:
2885-0660-XXX	2885-0661-XXX	2885-0662-XXX	2885-0663-XXX

## Back Supports

### Spex Contouring Back System

For Size 1 BINGO 2885-2012-XXX

For Size 2 BINGO 2885-2014-XXX

Spare Cover: 2885-2062-XXX Spare Cover: 2885-2064-XXX



## Head Supports

Standard Lateral	Extended Lateral	Mount and Stern Hardware	
	With Stylo hardware & BINGO attachment		
For moderate head positioning needs. 2885-4241-XXX Spare Cover: 2885-4261-XXX	For complex head positioning needs. 2885-4341-XXX Spare Cover: 2885-4361-XXX	For use with other Spex head pads 2885-4510-000	

## Arm Supports

Pelvic Supports

### Arm Supports

With height-adjustable, quick-release hardware and neoprene cover.

To fit Size 1 2885-5621-291

To fit Size 2 2885-5622-291



### Pelvic Supports

Mounts onto arm support bracket (pair) behind cushion trough (arm supports must be selected).

To fit Size 1 2885-1111-000

To fit Size 2 2885-1112-000



## Lateral Supports


### Narrow lateral supports

Soft padding/neoprene cover (pair) with Axial Technology. 8.5cm high x 11.5cm deep 2885-3410-291



#### **Broad lateral supports**

Soft padding/neoprene cover (pair) with Axial Technology. 10.5cm high x 13cm deep 2885-3431-291



### **OPTIONAL SELECTION ITEMS**

## Shoulder Harness

## X-Small Centrepoint

Small Centrepoint

Pad height 26cm, total pad width 22cm 1409-6660-017

### Pad height 36cm, total pad width 32cm 1409-6661-017



## Hip Belt

### 2 Point Hip Belt 4 Point Hip Belt 2-point, padded dual centre pull, side release buckle. 4-point, padded dual centre pull, side release buckle. X-Small (Pad width 30cm) 1405-6220-027 4-point, padded dual centre pull, side release buckle. Small (Pad width 36cm) 1405-6221-027 Small (Pad width 36cm) 1406-6421-027

## Feet Positioning

### Foot positioning sandals

### Height kit / raiser blocks

**Size 1** (10/11") 2885-5321-000

Footplate Spacer Block\*

**Size 2** (12/13") 2885-5322-000



May require height kit/s. 2885-6081-000



For foot positioning sandals (kit of 4). 2885-6091-000



\*Can not be used with foot positioning sandals

## Critical Care Options

IV Pole		xtra Bracket Oxygen Bottle Holder		
With mobility bracket for attachment to mobility base. Fits both Size 1 and Size 2. 2880-8275-000	to secon other in	mounting of IV pole nd base (Cobra or door base). 3275-050	With mounting bracket for attachment to mobility base. 120mm diameter for 'C Type' bottle. Fits both Size 1 and Size 2. 2880-8253-000	
Chassis Mount Trav				

**Chassis Mount Tray** For high dependency equipment - suction/oxygen etc.

Size 1 (for BINGO Size 1) 2880-8058-000

Size 2 (for BINGO Size 2) 2880-8059-000



## Swivel Lock

### For both Size 1 and Size 2

To hold castors straight when brisk walking. 2880-8440-000



## Hand Strap

Wrist safety tether strap

2005-6930-000



## Tray Therapy Tray Infill

### **Therapy Tray Infill**

For smaller trunk cut out. 2880-7180-000



## Indoor Base

Cobra high-low indoor frame with 75mm castors

Size 1 (for BINGO Size 1) 2880-8711-000

Size 2 (for BINGO Size 2) 2880-8721-000



## Separate Seat Unit

Use with Cobra Base

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

Size 1 2885-0011-032 **Size 2** 2885-0012-032



Manufactured by Spex Ltd solutions@spexseating.com spexseating.com

