

# Carrot 3000




**CAR SEAT**

Distributed by  
**medifab**

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	



# Seat Specifications



Age Range	User Weight Range	Seat Depth	Shoulder Height	Seat Weight	Testing Data
3-15 years	13.6-49kg	25-40 cm	32-47 cm 45-60 cm (when shoulder protector wings are added)	8.7kg	Tested to ECE R44/04 <div>ECE R44/04 Universal 15-36 kg S (E1) 04 30 1283 MARLBUSH CARROT II</div>

# Basic Packages

Includes:		
<ul style="list-style-type: none"><li>Carrot seat unit</li><li>5-point harness</li><li>Headrest wings</li><li>Tee-pad</li><li>Free recline backrest</li></ul>	4116-0010-259	

# Behavioural Packages



Includes:		
<ul style="list-style-type: none"><li>Carrot seat unit</li><li>5-point harness</li><li>Headrest wings</li><li>Tee-pad</li><li>Free recline backrest</li></ul>	Small Package with Anti-Escape Buckle 4116-0010-026	
	Small Package with Pen-Key Buckle 4116-0010-025	

Swivel Base Attachment	Swivel Base Protection Plate		
<p>ISOfix essential. (Standard adjustment of mandatory footprop is 34–40cm). 4116-0180-104</p>		<p>Foam plate to protect car seat cover, mounts under swivel base 4116-0190-000</p>	




# Seat Depth Growth

For growth, please select one of the following 3 rodded options, and associated extensions below, if required:  
(Rodded components cannot be used in combination with each other).

## Option 1

Rodded 50mm Seat Depth Extension		50mm Growth Pad
4116-0281-259	+	4116-0285-259
		

## Option 2

Rodded Waterfall Seat Depth Extension		50mm Growth Pad*	100mm Growth Pad*
4116-0283-259	+	4116-0285-259	4116-0289-259
			

## Option 3

Rodded Foot Support		50mm Growth Pad*	100mm Growth Pad*
<p><b>Small</b> Lower leg (D): 25-30 cm 4116-5241-250</p> <p><b>Medium</b> Lower leg (D): 30-40 cm 4116-5242-250</p> <p><b>Large</b> Lower leg (D): 40-50 cm 4116-5243-250</p>	+	4116-0285-259	4116-0289-259
			
		Adaptor	
		<p>For use of foot support with swivel base 4116-5120-000</p>	




\*50mm and 100mm growth pads cannot be used in combination with each other.

# Back Support Height Growth





50mm Height Extension	100mm Height Extension	Shoulder Protector Wings
50mm back support height extension infill pad. 4116-2481-259 	50mm back support height extension infill pad. 4116-2482-259 	(For children over 45cm shoulder height). 4116-2753-259 

\*50mm and 100mm infill pads cannot be used in combination with each other




# Installation Options

Tether Strap	ISOfix lower straps	Airplane Fixing Strap
4116-6953-017 	No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6910-000 	For securing to aircraft seat. 4116-6960-017 



# Cover Options

Wipe-down			
Wipe-down seat base. 4116-0476-209		Tee-pad, Standard. 4116-6553-209	
Wipe-down backrest. (Main unit only). 4116-2476-209		Tee-pad, Large. 4116-6558-209	
Removable Incontinence Cover		Air Mesh Cushion Pad	
Attaches using Velcro. Unable to be attached to wipe down covers. 4116-0475-209		Attaches using Velcro. Unable to be attached to wipe down covers. 4116-8880-099	



# Accessories (continued)

Tray / Forearm Support		Tilt Wedge	Upgraded Lateral Supports	
4116-5534-259		Provides 20° tilt. 4116-0291-259		Firm padding for added trunk control. 4116-3250-090 

Seat Belt Holder		Shoulder Strap Retainer	
Designed to keep seat belt away from the neck. 4116-6790-017		<b>Small</b> 4116-6721-000	
		<b>Medium</b> 4116-6722-000	
		<b>Large</b> 4116-6723-000	

Medial Thigh Support		Hinge Locking Belt	
Detachable & Washable. Attaches to chest support and crotch strap. 4116-1510-000		For easy carrying. 4116-8170-000	

Hip Belt		Accessory Cross Strap	
<b>Small</b> 4116-6531-027		<b>Anti-Escape</b> 4604-3004-026	
<b>Medium</b> 4116-6532-027		<b>Pen Release</b> 4604-3004-025	
<b>Large</b> 4116-6533-027		<b>Side Release</b> 4604-3112-019	

Deeper Head Support Wings		Larger Tee Pad	
4116-4653-090		4116-6542-049	
		Cover 4116-6558-259	
		Black Spextex Cover 4116-6558-209	



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Carrot\_20231215\_NZ