

# Juni Chair



## INDOOR SEATING

### SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Size & Colour Options

	Size 1	Size 2	Size 3	Size 4
<b>Black</b>	3209-1001-099	3209-1002-099	3209-1003-099	3209-1004-099
<b>Blue</b>	3209-1001-060	3209-1002-060	3209-1003-060	3209-1004-060
<b>Red</b>	3209-1001-022	3209-1002-022	3209-1003-022	3209-1004-022
<b>Lime</b>	3209-1001-054	3209-1002-054	3209-1003-054	3209-1004-054


Seat Unit Measurements				
Seat height (cm)	16-22	20-28	24-34	30-40
Seat width (cm)	24	28	32.5	36.5
Seat depth (cm)	18-24	22-30	26-36	32-42
Back support height (cm)	26	29	32.5	40
Load capacity (kg)	30	40	50	60
Product weight (kg)	4	6	8	10



### Base package includes:



- Seat Unit
- Hip Belt




# Accessories

Side Cushions - For Size 1 & 2 Juni	Side Cushions - For Size 3 & 4 Juni	
Reduce seat width using side supports.		
To reduce width by 2.5cm 3209-0001-025	To reduce width by 2.5cm 3209-0003-025	
To reduce width by 5cm 3209-0001-050	To reduce width by 5cm 3209-0003-050	
To reduce width by 7.5cm 3209-0011-050	To reduce width by 7.5cm 3209-0013-050	

Standard Skis		Mobile Skis	
Makes moving Junji chair in a straight line easier (under tables).		Makes the Junji chair easier to move, ideal for pushing under tables!	
<b>Size 1</b> 3209-5021-000		<b>Size 1</b> 3209-5051-000	
<b>Size 2</b> 3209-5022-000		<b>Size 2</b> 3209-5052-000	
<b>Size 3</b> 3209-5023-000		<b>Size 3</b> 3209-5053-000	
<b>Size 4</b> 3209-5024-000		<b>Size 4</b> 3209-5054-000	

Sliding Foot Support (for use with mobile skis)		Adjustable & Foldable Table	
Add foot support to mobile skis to support users legs.		Adjustable height work surface for children which can be easily folded for storage, a great companion to the Junji chair.	
<b>Size 1</b> 3209-5071-000		<b>Size 1</b> Table height 31-45cm 3209-7191-000	
<b>Size 2</b> 3209-5072-000			
<b>Size 3</b> 3209-5073-000		<b>Size 2</b> Table height 40-60cm 3209-7192-000	
<b>Size 4</b> 3209-5074-000			

Lateral Supports - For Size 1 & 2 Juni	Lateral Supports - For Size 3 & 4 Juni	
Added support to help with sitting,		
3209-3001-000	3209-3003-000	

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY: