

Duro Stroller



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Product Details

	Size 1	Size 2
Turquoise	2780-8801-349	2780-8802-349
Black	2780-8801-346	2780-8802-346

Seat Specifications		
Seat width	23-34cm	28-39cm
Seat depth	26-40cm	31-45cm
Seat angle	16 or 20°	16 or 20°
Back height	45-63cm	61-83cm
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length	20-44cm	20-44cm
Folded size (L x(W) x H)	84 x 63 x 32cm	84 x 68 x 38cm

Base Specifications		
Max weight	75kg	75kg
Front wheel diameter	20cm	20cm
Rear wheel diameter	30cm	30cm



Stroller package includes:

- Hoggi Duro Rehab-Buggy
- Hip Support, Aluminium
- Seat and Backrest Upholsteries
- Small Spex Hip Belt
- Canopy
- PU Wheels with Black Rim
- Swivelling Front Wheels
- Supporting Brackets
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support with Heel Stop
- Accessories Bag (Max weight: 3kg)
- Tie Down Kit (According to ISO 7176-19)
- Retainer for Accessories
- User Manual










Items outlined in blue are required as standard package

Upholstery Colour Options

Turquoise		Black	
Seat Upholstery (required)			
Size 1 2780-8801-349		Size 1 2780-8801-346	
Size 2 2780-8802-349		Size 2 2780-8802-346	

Head Supports

Head Support Pads	
Size 1 11cm Deep x 16cm High 2880-4471-346	
Size 2 13cm Deep x 18cm High 2880-4472-346	

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209

Lateral Supports


Standard Lateral Supports		
Height and angle adjustable. All sizes are compatible with both Duro sizes.		
Size 1 8cm(H) x 10cm(L) 2880-3010-346	Size 2 10cm(H) x 12cm(L) 2880-3020-346	Size 3 12cm(H) x 13cm(L) 2880-3030-346



Lateral Supports (continued)

Spex Teardrop Lateral Supports		Spex Broad Lateral Supports	
<p>Size 1 13.5cm(H) x 7.5cm(L) 2880-3008-500</p>	<p>Size 2 15.5cm(H) x 10.5cm(L) 2006-3120-099</p>		<p>10.5cm(H) x 13cm(L) 2880-3008-310</p> 



Harness Options

Spex Centrepoint Harness	5-Point Harness
<p>X-Small Recommended for size 1 Duro Total pad dimensions: 26cm(H) x 22cm(W) 1409-6660-017</p>	<p>With front pull adjusters and neck padding 2880-6511-346</p> 
<p>Small Recommended for size 1 & 2 Duro Total pad dimensions: 36cm(H) x 32cm(W) 1409-6661-017</p>	
<p>Medium Recommended for size 2 Duro Total pad dimensions: 36cm(H) x 32cm(W) 1409-6662-017</p>	


Hip Belt Options

Spex Push Button	Pen Popper
<p>Small Pad width 36cm 1405-6211-028</p>	<p>Small Pad width 36cm 1405-6211-025</p>
<p>Medium Pad width 46cm 1405-6212-028</p>	<p>Medium Pad width 46cm 1405-6212-025</p>


Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			
To fit Size 1 2880-7111-000	To fit Size 2 2880-7112-000		


Grab Rail	
Cannot be used with therapy tray	
To fit Size 1 2880-7251-346	To fit Size 2 2880-7252-346



Hand Strap

Wrist safety tether strap	
2005-6930-000	

Swivel Lock

To hold castors straight when brisk walking	
2880-8440-000	

medifab

Distributed by Medifab Ltd

22 Factory Drive
Rolleston 7675
New Zealand

0800 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY: