

Manual Wheelchair Comparison Chart



SPECIFICATIONS	NITRUM	7 ^R	QUICKIE 2 LITE	QUICKIE 2	M6	IRIS
SEAT WIDTH	12 - 20"	12 - 20"	14 - 20"	12 - 22"	22 - 30"	14 - 22"
SEAT DEPTH	12 - 20"	12 - 20"	14 - 20"	12 - 20"	18 - 22"	14 - 22"
FRONT SEAT HEIGHT	16 - 21"	16 - 21"	14.5 - 21.5"	13.5 - 22"	17 - 20"	12.5 - 19"
REAR SEAT HEIGHT	13 - 20"	13 - 20"	13 - 19.5"	13 - 22"	17 - 20"	12.5 - 19"
FRAME STYLE	Z-Frame 7020 Aluminium	Z-Frame 7000 Aluminium	Folding frame with fixed or Swing-away front	Folding frame with fixed or Swing-away front	Folding Heavy Duty	Tilt In Space
BACK HEIGHT	Height Adjustable: 10" - 20" Fixed Height: 8" - 19"	8 - 20"	9 - 21"	9-21"	17 - 20"	15 - 24"
BACK SUPPORT ANGLE	-18° - Aft +10°	-27° - +8°	-14° - +5° (Optional)	-15° - +3° (Optional)	N/A	-5° - +30°
FRAME / HANGER ANGLE	75° / 80° / 85°	75° / 80° / 85°	70° / 80°	60°/ 70°/ 80°/ 90°	N/A	60° / 70° / 80° / 90° / Elevating
INSET	0" / 1" / 2"	0" / 1" / 2"	0 / 1" each side	N/A	N/A	N/A
CAMBER	0° / 3° / 6°	0° / 3° / 6°	(Adjustable) Set at 3°	(Adjustable) Set at 3°	0°	0°
MAX. USER WEIGHT	125kg	120kg	120kg	136kg Std / 158kg HD	295kg	136kg / 181kgs HD
PRODUCT WEIGHT	From 8.3kg	From 8.5kg	From 10.6kg	From 11.5kg	From 24kg	From 17.6kg
TRANSIT APPROVED	Yes*	Yes*	Yes*	Yes*	No	Yes
FRAME COLOUR	30 colours	29 colours	24 colours	Black Opal Standard, 24 colours	24 colours	Black Opal Standard, 25 colours
WARRANTY	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime
FOLDING	No	No	Yes	Yes	Yes	No

* When fitted with Transit Option Conforms to ANDI/RESNA WC-19 Standards. Please refer to the Quickie & Zippie WC-19 Transit Securement Supplement

Also refer Paediatric Manual Wheelchair Comparison Chart for paediatric wheelchair solutions.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.