

Simple Stuff Works



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

March 2024

Client details

Client/User name:			
DOB:		Date:	
Therapist:			

Notes

Further Information

CHOOSE A KIT TO USE ON **FOAM MATTRESS**

Kit Suggestion

Baby Kit	Child Kit	Teen/Adult Kit
<p>0-5 Years</p> <p>Package includes:</p> <ul style="list-style-type: none">• Cot Stabilising Mesh Pack• Cot Temperate Fibre Topper• Fitted Cotton Sheets – Cot Size• Small Laterals• Small Wedges• Small Horseshoe pillow• Simple Stabiliser Size 1• Travel Bag <p>7150-7001-C00</p>	<p>5-12 Years</p> <p>Package includes:</p> <ul style="list-style-type: none">• Single Bed Stabilising Mesh Pack• Single Bed Temperate Fibre Topper• Fitted Cotton Sheets – Single Bed Size• Medium Lateral Supports• Small Wedges• Small Lateral Supports• Simple Stabiliser Size 2• Small Terry Towelling Pillow• Travel Bag <p>7150-7002-S00</p>	<p>12+ Years</p> <p>Package includes:</p> <ul style="list-style-type: none">• Single Bed Stabilising Mesh Pack• Single Bed Temperate Fibre Topper• Fitted Cotton Sheets – Single Bed Size• Medium Lateral Supports• Large Lateral Supports• Medium Wedges• Large Terry Towelling Pillow• Travel Bag <p>7150-7003-S00</p>

OR REFER TO **PAGE 3** TO **CONFIGURE YOUR OWN**

CONFIGURE YOUR OWN
FOLLOW THE PATHS DEPENDING ON CLIENT'S MATTRESS



Choose your set for foam mattress

Cot Bed	Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x cot temperate toppers • 2 x cot fitted sheets <p>7150-7021-C00</p>	<p>Set Includes:</p> <ul style="list-style-type: none"> • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x single temperate toppers • 2 x single fitted sheets <p>7150-7022-S00</p>	

Choose your set for air mattress

Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> • 2 x single temperate toppers • 2 x single fitted sheets <p>7150-7024-S00</p>	




Soft wedges for foam mattress

Soft wedges for air mattress

Wedges with padded covers (pair)		Wedges with padded covers and integrated mesh (pair)	
<p>Small 7150-7351-000</p>	<p>Medium 7150-7352-000</p>	<p>Small 7150-7326-000</p>	<p>Medium 7150-7327-000</p>




Lateral supports for foam mattress

Lateral Supports with padded covers (pair)

Small 7150-7201-000	Medium 7150-7202-000	Large 7150-7203-000
		

Lateral supports for air mattress

Lateral Supports with padded covers and incorporated mesh (pair)

Small 7150-7211-000	Medium 7150-7212-000	Large 7150-7213-000
		

Neck support pillow

Promotes neutral neck alignment


(Suitable for air & foam mattresses).

300mm x 350mm x 100mm

7150-4050-000

Spare terry cotton cover

7150-4055-000




Terry Towling Pillow

Assists with head positioning

<p>Small A: 480mm B: 350mm 7150-7501-000</p>	<p>Large A: 680mm B: 450mm 7150-7502-000</p>	
<p>Spare terry cotton cover 7150-7511-000</p>	<p>Spare terry cotton cover 7150-7512-000</p>	


Roll Pillows

Assists with trunk and/or thigh positioning

<p>Small A: 150mm B: 400mm</p>	<p>Large A: 180mm B: 500mm</p>	
<p>With breathable inner & terry cotton outer 7150-7851-000</p>	<p>With breathable inner & terry cotton outer 7150-7852-000</p>	
<p>With wipeable incontinence inner & terry cotton outer 7150-7855-000</p>	<p>With wipeable incontinence inner & terry cotton outer 7150-7856-000</p>	
<p>Spare terry cotton cover 7150-7857-000</p>	<p>Spare terry cotton cover 7150-7858-000</p>	


Horseshoe Pillows

Multi-purpose positioning pillow

<p>Small A: 190mm B: 460mm</p>	<p>Large A: 200mm B: 500mm</p>	
<p>With breathable inner & terry cotton outer 7150-7711-000</p>	<p>With breathable inner & terry cotton outer 7150-7712-000</p>	
<p>With wipeable incontinence inner & terry cotton outer 7150-7721-000</p>	<p>With wipeable incontinence inner & terry cotton outer 7150-7722-000</p>	
<p>Spare terry cotton cover 7150-7731-000</p>	<p>Spare terry cotton cover 7150-7732-000</p>	





Simple Stabiliser

Supports hips and thighs in neutral supine alignment

Size 1	Size 2	Size 3	Size 4
H: 11cm D: 11cm W: 42cm To fit approx leg circumference of 25cm. Suggested age range 2 – 4yrs.	H: 19cm D: 35cm W: 46cm To fit approx leg circumference of 35cm. Suggested age range 2 – 10yrs.	H: 25cm D: 40cm W: 65cm To fit approx leg circumference of 55cm. Suggested age range 5 – 13yrs.	H: 28cm D: 40cm W: 65cm To fit approx leg circumference of 60cm. Suggested age range 13 – adult.
With breathable inner & terry cotton outer 7150-7341-000	With breathable inner & terry cotton outer 7150-7342-000	With breathable inner & terry cotton outer 7150-7343-000	With breathable inner & terry cotton outer 7150-7344-000
With wipeable incontinence inner & terry cotton outer 7150-7345-000	With wipeable incontinence inner & terry cotton outer 7150-7346-000	With wipeable incontinence inner & terry cotton outer 7150-7347-000	With wipeable incontinence inner & terry cotton outer 7150-7348-000
			
Spare terry cotton cover 7150-7221-000	Spare terry cotton cover 7150-7222-000	Spare terry cotton cover 7150-7223-000	Spare terry cotton cover 7150-7224-000
Spare wipeable inner cover 7150-7225-000	Spare wipeable inner cover 7150-7226-000	Spare wipeable inner cover 7150-7227-000	Spare wipeable inner cover 7150-7228-000






Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3	
A: 270mm B: 380mm C: 280mm D: 80mm E: 130mm F: 80mm	A: 400mm B: 560mm C: 340mm D: 100mm E: 200mm F: 90mm	A: 400mm B: 610mm C: 400mm D: 110mm E: 150mm F: 70mm	
With breathable inner & terry cotton outer 7150-7971-000	With breathable inner & terry cotton outer 7150-7972-000	With breathable inner & terry cotton outer 7150-7973-000	
With wipeable incontinence inner & terry cotton outer 7150-7981-000	With wipeable incontinence inner & terry cotton outer 7150-7982-000	With wipeable incontinence inner & terry cotton outer 7150-7983-000	
			
Spare terry cotton cover 7150-7991-000	Spare terry cotton cover 7150-7992-000	Spare terry cotton cover 7150-7993-000	

Foot Supports

Provides positioning & skin protection for the feet (pair)

Size 1	Size 2	Size 3	
A: 21cm (max) 7150-5351-000	A: 25cm (max) 7150-5352-000	A: 30cm (max) 7150-5353-000	
			



Optional Extras

Stabilising Mesh Set (suitable for foam mattresses only)

Used to stabilise equipment on bed during sleep.
Includes 2x mesh sheets; 1 bottom, 1 top

Cot Bed Size
7150-7011-C00

Single Bed Size
7150-7012-S00





Temperate Topper

Used to allow air to circulate and take moisture away from the body.

Cot Bed Size
7150-7051-C00

Single Bed Size
7150-7052-S00




Fitted Sheets

To regulate body temperature and ensure that the client is comfortable.

Cot Bed Size
7150-7121-C00

Single Bed Size
7150-7122-S00



SSW Sleep System Carry Bag

7060-8180-000



Note: For co-sleeping needs, please get in touch with our friendly team.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: SimpleStuffWorks_20240312