

# Q700

## POWERCHAIR SCRIPTING AND ORDERING GUIDE

Distributed by  
**medfab**



### CLIENT DETAILS

Client/User Name	
Assessor	
Funding Body	
User Weight	
Date	
Notes	

### Q700 BASE PACKAGE

- R-Net Expandable Electronics
- Lights
- 14" Drive Wheels
- 9" Castors
- Transit Tie Downs
- 4 Pole Motors
- 70Ah Gel Batteries
- Battery Charger

### DRIVE WHEEL OPTION

Front	Mid	Rear
-------	-----	------

### SEAT SELECTION

Seat Type	
Sedeo Lite Seat Top	(Width: 15" - 20" Depth: 16" - 22")
Sedeo Pro Short Seat Top	(Width: 16" - 22" Depth: 14" - 18")
Sedeo Pro Long Seat Top	(Width: 16" - 22" Depth: 17" - 22")

### SEAT SIZE SELECTION

Seat Width:	Seat Depth:
-------------	-------------

### SEAT TO FLOOR HEIGHT SELECTION

16.5	17	17.5	18	18.5
19	Other:			

### POWER SEAT FUNCTIONS

Elevate 12"	Tilt 50°	Power Recline 130°
Power Recline 150° with Anti-Shear		

### LEGREST OPTIONS

Legrest Type	
Swing Away 70°	Swing Away Manual Elevating
Swing Away Power Elevating	Fixed Centre Mount
Powered Center Mount	

Stump Support	
Left	Right

### ARM SUPPORT OPTIONS

Arm Support Assembly	
Sedeo Lite Cantilever Flip Back Armrests (N/A with PWR recline)	Single Post Flip Back Armrests
Sedeo Pro Cantilever Flip Back Armrests	Dual Locking Armrest Post for Cantilever Flip Back

Armrest Pads	
Sedeo Lite Cantilever Armrest Pads 380mm (N/A with PWR recline)	Sedeo Pro Contour Arm Pads 300mm
Sedeo Pro Contour Arm Pads 400mm	Comfort Gel 2" Armrest Pads
Comfort Foam 2" Armrest Pads	Gel Arm Troughs Medium
Gel Arm Troughs Large	Spex Comfi Wide Armrests

### JOYSTICK OPTIONS

Joystick	
Rnet Joystick Controller Standard LED	Rnet Joystick Controller LCD Screen CJSM1
Rnet Joystick Controller Advanced LED Screen CJSM2 with Assignable Buttons, Bluetooth & IR	Attendant Control

Joystick Mount	
Swing Away Interchangeable	Round Style Swing Away Joystick Mount - Left
Midline Joystick Mount	Round Style Swing Away Joystick Mount - Right

Joystick Knob Option		
Goal Post	T-Bar	Mushroom
Solid Ball	Stick	Large Softball

### ACCESSORIES

Armrest Bag	Ram Tablet Mount
Side Guards	Ram Cup Holder
Tray Table	Rear View Mirror
Transfer Bar - Pair	Control 5 Switch Box - Button Controls
Ram Phone Mount	Control 5 Switch Box - Toggle Controls
Upgrade motors to Suretrac	Upgrade to 8km Hi Torque motors
Upgrade to 13kmh motors (N/A with FWD)	

Colour Accent Kit		
White	Orange	Green
Red	Blue	Black