

Jenx Junior+



**SEATING SYSTEM
FOR MILD, MODERATE
AND COMPLEX NEEDS**

Distributed by
medifab


SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Notes

User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	



Specifications

	Size 1	Size 2	Size 3
Average age	18 Months–5 Years	3–8 Years	5–12 Years
Weight capacity (kg)	50	50	50
Seat width (cm)	15*–37.5	15*–37.5	15*–37.5
Seat depth (cm)	21–30.5	24–35.5	26–44
Chest support width (cm)	18.5–28	18.5–28	21–30
Seat-top of backpad height (cm)	34.5–51**	34.5–51**	41–57.5**
Footplate angle (°)	-35–5	-35–5	-35–5
Footplate height (cm)	19.5–27	23.5–35	25–40.5
Backrest recline (°)	-15–15	-15–15	-15–15
Tilt range (°)	30 (Posterior)–10 (Anterior)	30 (Posterior)–10 (Anterior)	30 (Posterior)–10 (Anterior)


* 15cm achieved only when using the hip reducing pads on Page 3.

** 51cm back height in Size 1 + 2 and 57.5cm in Size 3 only achieved with back post extension on Page 4.

Colour Options

Red	Lime	Blue
		
Fabric Cover -022	Fabric Cover -054	Fabric Cover -066
Wipeable Cover -W22	Wipeable Cover -W54	Wipeable Cover -W66


Medifab recommends that items outlined in blue are selected as standard package.

 **Therapists Normally Choose**




Base

Y-Base

Height adjustment of 34 - 76cm, lockable swivel castors, tilt-in-space; 30° posterior, 10° anterior.
3615-8740-000




Seat Unit

Size 1	Size 2	Size 3
Age range: 18 months–4 years. Max user weight: 50kg. 3610-1001-000	Age range: 3–8 years. Max user weight: 50kg. 3610-1002-000	Age range: 6–12 years. Max user weight: 50kg. 3610-1003-000
		
<p>Contoured PU seat with removable abduction pad and memory foam seat pad, adjustable and flip-away PU hip pads, angle-adjustable footplate, 4-point hip belt, height-adjustable PU back pad, an independently adjustable PU lumbar/sacral pad, arm supports, tray table and push handle.</p>		



Covers Required:


Seat Cover Fabric 3612-0421-XXX	Seat Cover Fabric 3612-0423-XXX	Seat Cover Fabric 3612-0423-XXX
Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2423-XXX
Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX
Cover set (wipe clean only) 3612-0521-WXX	Cover set (wipe clean only) 3612-0522-WXX	Cover set (wipe clean only) 3612-0523-WXX

Hip Width Reducing Pads (Pair)


2.5cm	
Reduces Hip Width by 2.5cm on each side. 3612-1040-000	

Head Supports



Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm (W) x 5cm (D) x 12cm (H) 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm (W) x 5cm (D) x 12cm (H) 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Oval Head Support	
19cm (W) x 10cm (H) 5356-4611-XXX	

Shoulder Protractors (Pair)

Size 0	Size 1	Size 2	
Can be used on all three size seat units. Pad size 13cm (H) x 13.5cm (D) 3612-2750-000	Can be used on all three size seat units. Pad size 21cm (H) x 14.5cm (D) 3612-2751-000	Can be used on all three size seat units. Pad size 23.5cm (H) x 17.5cm (D) 3612-2752-000	
Covers (required) 3612-2760-XXX	Covers (required) 3612-2761-XXX	Covers (required) 3612-2762-XXX	

Trunk Supports (Pair)

Flexi Thoracics			
Size 1 18cm (L) x 8cm (H) 3615-3251-000	Size 1 Long  22cm (L) x 8cm (H) 3615-3252-000	Size 2 22cm (L) x 11cm (H) 3615-3253-000	
Covers (required) 3615-3271-XXX	Covers (required) 3615-3272-XXX	Covers (required) 3615-3273-XXX	

Waistcoat Harness

Size 1	Size 2	Size 3
A 5cm B 28cm C 12cm D 13.5cm 3615-6701-000	A 5cm B 32cm C 14cm D 15.5cm 3615-6702-000	A 6.5cm B 37cm C 17.5cm D 18cm 3615-6703-000




Back Post Extension

7.5cm	
7.5cm Back Post Extension for extra backrest height 3610-2480-000	


Abduction Block


Size 0	Size 1	
Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm 3615-1650-000	Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm 3615-1651-000	

Leg Positioning


Leg Control Pads	
Control leg position independently of pelvic position, flip completely out of the way to facilitate smooth and easy transfers. 12.5cm (H) x 7.5cm (D) 3615-1861-000	
Covers (required) 3615-3171-XXX	

Sandals


Size 1	Size 2	Size 3	
17cm (L) x 8cm (W) 5356-6081-000	22cm (L) x 9.5cm (W) 5356-6082-000	27.5cm (L) x 10.5cm (W) 5356-6083-000	

Sandal Raiser Blocks			
Size 1 Kit of 4x1.5cm blocks. 5356-6091-000	Size 2 Kit of 4x1.5cm blocks. 5356-6092-000	Size 3 Kit of 4x1.5cm blocks. 5356-6093-000	


Ankle Snugs

Size 1	Size 2 & 3	
Inner Circumference: 12cm - 22.5cm 5356-6131-000	Inner Circumference: 16.5cm - 29.5cm 5356-6132-000	

Foot Plate Pad

Foot Plate Pad	
Footplate pad not compatible with sandals. 4cm thick. 3615-5301-000	

Accessories for Arm Support & Tray

Tray Infill		
Great for attaching suction toys.		
Size 1 3610-7180-000	Size 2 3610-7182-000	

Half Tray Padding	Full Tray Padding
Comfort for forearms and elbows, whilst still enjoying a hard tray surface.	Added comfort for hands.
Size 1 3612-7441-000	Size 1 3612-7451-000
Size 2 3612-7442-000	Size 2 3612-7452-000

Grab Bar	Activity Frame
Used for children who need to fix one hand to increase functionality. 54cm (W) x 9.5cm (D) x 7cm (H) 3615-7200-000	Great for encouraging head control and hand to eye co-ordination. 54cm (W) x 16cm (D) x 35cm (H) 3615-7475-000

Elbow Blocks	Arm Support Covers
Used to encourage arms into a mid-line position and/or to block excessive retraction. 10cm (W) x 14cm (H) 3615-7350-000	3615-5671-XXX

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: JenxJunior+_20241030