

Standing Frame Comparison Chart



Specifications	Prone only								Supine or prone								Supine only		
	JENX MONKEY	JENX PRONE STANDER		DANDY			SMART-STAND HIP		JENX MULTISTANDER		JENX STANDZ	ELLA STANDER				JENX SUPINE			
		Size 1	Size 2	Size 1	Size 2	Size 3	100	200	Size 1	Size 2		Size 1	Size 2	Size 3	Size 4	Size 1	Size 2	Size 3	
Description	Gentle introduction to standing with easy to adjust positioning.	Dedicated to assisted prone standing with the child supported at the front		Popular standing frame that gives tool-less adjustability for a wide range of basic needs.			Ultimate 'body-hugging' prone standing frame.		Highly supportive standing frame in prone, supine, or upright.		Abducted standing systems offering prone and supine abducted standing.	Versatile standing configurations for symmetrical and asymmetrical positions.				Robust and easy to use standing solution.			
APPROX AGE OF USER*	9mths-4yrs	3-7yrs	6-14yrs	1-4yrs	4-8yrs	6-15yrs	1-5yrs	4-12yrs	9mths-6yrs	3-13yrs	1-9yrs	1-3yrs	2-8yrs	6-12yrs	10-18yrs	9mths-4yrs	3-11yrs	10-Adult	
USER WEIGHT (kg)	25	55kg	70kg	35	45	55	40	60	35	60	45	30	45	55	75	25	50	100	
MAX USER HEIGHT (cm)	-	-	-	85	125	145	110	150	116	175	138	60-95	90-130	120-155	130-180	105	145	190	
FOOT TO AXILLA (cm)	(Top of chest pad) 51-75	66-100	82-127	39-57	60-86	76-110	50-90cm	67-113	Prone 49-86 Supine 57-86	Prone 70-117 Supine 77-117	58-100	40-60	56-88	77-106	89-113	50-107	90-147	154-164	
CHEST WIDTH (cm)	15-23	19-29	23-35	17-20	25-30	25-35	21 (min.)	26 (min.)	15-26	17-32	18-30	14-22	16-29	19-32	26-37	13-31	13-31	22-44	
HIP WIDTH (cm)	15-23	19-29	23-35	17-20	25-30	25-35	21 (min.)	26 (min.)	12-26	17-32	18-30	14-22	16-29	19-32	26-37	13-31	13-31	22-44	
SUPINE TILT	N/A	N/A		N/A			N/A		90-10°		90°-0°		0-90°				0-80°		
PRONE TILT	20°	0-40°		19°	14°	11°	0-80°	0-70°	10-90°		0-90°		0-90°				N/A		
MULTI ADJUST KNEE SUPPORT	✗	✗		✓ (Joined together)			✓		✓		✓		✓				✓		
MULTI ADJUST FOOT SUPPORTS	✓	✓		✓			✓		✓		✓		✓				✓		
TRAY TABLE	✓	✓		✓			✓		✓		✓		✓				✓		
HEAD SUPPORT	✗	✗		✗			Custom		✓ (Supine)		✓ (Supine)		✓ (Supine)				✓		
POSTURAL SUPPORT LEVEL	Mild/Mod.	Mild/Mod.		Mild			Mild/Mod.		Mild/Mod.		Mild/Mod./Complex		Mod./Complex				Mod./Complex		
FRAME FOOTPRINT (cm)	57x88	90x65	99.5x65	63x67	68x80	72x87	66x50	86x80	56x71	65x95	63.5x90	66x53	78x58	104x71	106x76	57x66	65x95	65x95	
LEG ABDUCTION RANGE**	~10°	~10°		~10°			~12°		~10°		30°		25°	20°	25°	20°	Neutral		
OPERATION	Gas Assist	Gas Assist		Gas Assist			Gas Assist		Gas Assist (size 2 electric option)		Gas Assist		Gas or Electric				Battery/Electric		
KNEE CONTRACTURE ACCOMMODATE	Nil	Nil		Mild			Mild		Supine: Extensive Prone: Mild		Supine: Extensive Prone: Mild		Supine: Very Extensive Prone: Mild				Mild	Mild	Extensive

* Indicator only, for accurate product selection client measurements are the best guide.
** Measured from vertical (0°) leg alignment to abducted angle

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.