

Spex Classic Back Support with WheelAir

Spex with WheelAir: Custom Comfort Meets Active Cooling

Discover the ideal combination of advanced posture support and active cooling with the Spex with WheelAir back support. Engineered for wheelchair users who value comfort, control, and adaptability, this system provides on-the-spot customisation with the reliable Spex Contouring Cube System while WheelAir's six cooling settings keep you cool and dry all day. Lightweight, low-noise, and USB-C chargeable, it's support and comfort—elevated.

FEATURES & BENEFITS

- 24 Hours Battery Life: Keeps you cool all day
- Six Cooling Settings: Adjust to your comfort level
- Wireless Remote: Seamless control with one touch
- Low Noise: Quiet operation to keep you at ease
- USB-C Chargeable: Convenient charging on the go



SIZE CHART

Wheelchair Height (Inches)	Wheelchair Width (Inches)										
	11	12	13	14	15	16	17	18	19	20	21
12											
13											
14											
15											
16											
17											
18											
19											
20											
21											
22											
23											

- Blue squares are standard sizes.
- Non-standard or custom sizes are available on request.

PRODUCT DETAILS

Width	Capacity	Product Weight
10"-15"	200lbs (90kg)	(16x16)
16"-20"	300lbs (136kg)	1.7kg

- Note: This is the weight of the back support only and does not include the hardware.

PRODUCT PERFORMANCE

Postural Control				4
Lateral Stability				5
Postural Asymmetries				5
Changing Needs				5
Pressure Care				5

- Note: This is a guide only and ratings are dependent on correct use of product