

Starlight-NXT Car Seat

Fully adjustable special needs seat, offering excellent posture support options and an inbuilt tilt function. Starlight-NXT is suitable for young children from 9kg, and will grow with the child right up to 36kg.

Excellent postural support



Inbuilt adjustment to the rear of the seat allows the backrest height to be increased for growth

Uses **vehicle seatbelt** to secure the child


Large abdominal pad **protects child** at harness buckle

5 point comfort harness with shoulder and abdominal padding, and centre pull adjuster

Crash tested to ECE R44-04

Package includes

- Backrest infill padding
- Soft shoulder padding
- Slim lateral trunk supports
- 10° tilt
- Base board accessory
- Swivel base system (left turning)
- Seat stabilisation clamp bow
- Grey upholstery




SLIM SWIVEL BASE
Allows more space and headroom inside the car when turning out.

Left Turning Swivel Package: 4141-0011-090
Right Turning Swivel Package: 4141-0011-R90

ACCESSORIES

ISOfix Adapter for Stability in Vehicle
4141-6910-100



Seat Protection Mat Protects Vehicle Upholstery
4141-0190-000



Adjustable Footrest
4141-5240-000



Padded Foot Roll
4141-5340-000



Forearm Support Tray
4141-5530-000



Abduction Pommel
4141-1510-000



Cushion to Reduce Seat Depth to 15 cm and Back Height to 25 cm
4141-6570-000



Extra 10° Tilt System (total 20° tilt)
4141-0241-015





Size	Code	User Weight	Seat Width	Seat Depth*	Shoulder Height	Car Seat Weight
Standard	4141-0011-090	9-36 kg	30 cm	30 cm	27-37 cm	8.5 kg
With XL Kit added	4141-0140-000		30 cm	34 cm	44-53 cm	

*Reducer cushion available