

# Daessy Mounting System



## DAESSY MOUNTING SYSTEM

Distributed by  
**medifab**

SCRIPTING AND ORDERING GUIDE

December 2024

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Rigid Mount

Simple, sturdy, cost effective and quick and easy to remove from a wheelchair. It can be lifted and removed or swung to the side for easy chair access. Variations can be made in the length and shape of the tubes, and in the other parts, to adapt the Rigid Mount to any wheelchair.



### For Manual Wheelchairs

1611-9083-000  
(TC90) Tube Connector 90° Elbow

**Main Component**

1611-9101-000  
(UFC1000IP) Frame Clamp Inner Piece 1" Round

*Refer to Spare Parts for other frame clamps options.  
Requires full access to the round tube*

1611-9064-000  
(S-22x3) Standard S-Bend Tube 22" length x 3" offset

*Refer to Spare Parts for other S-Bend Tube sizes  
or custom tubes*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

### For Power Wheelchairs

1611-9083-000  
(TC90) Tube Connector 90° Elbow

**Main Component**

1601-0020-M6SS  
Universal T-Nut/+Nut for Power Wheelchair Seat Rails.  
Fits Quickie Q series, Invacare TDX, Quantum

*Refer to Spare Parts for other TNUT options*

1611-0049-000  
(MH2/IPA) Multi-hole 2" Inner Piece Adapter

*Refer to Spare Parts for other seat rail attachment options*

1611-9064-000  
(S-22x3) Standard S-Bend Tube 22" length x 3" offset

*Refer to Spare Parts for other S-Bend Tube sizes  
or custom tubes*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

### For Fixed Wheelchairs

1611-0122-000  
(UFCOP) Frame Clamp Outer Piece

1611-0034-000  
(IC1) Index Clamp

*UFCOP and IC1 are required.*

### For Tilting Wheelchairs

1611-0071-000  
(RFCR) Removable Frame Clamp Receiver

1611-0076-000  
(ROP) Removable Outer Piece

*RFCR and ROP are required*

1611-0128-000  
(USBF) Folding Quick Release Base  
*Refer to Spare Parts for other Quick Release Bases*

1611-9008-000  
(UDS-22-DC) Rehadapt  
*Therapist choice*

# Folding Mount

The Folding Mount can quickly and easily be removed from the wheelchair. It folds from in front, down to the side of the wheelchair. It consists of two lengths of stainless-steel tube joined by a folding mechanism. The horizontal tube can be lifted through an angle of 90° and rotated down 180° to place the horizontal tube and vertical tube parallel along the side of the wheelchair.



## For Manual Wheelchairs

1611-0078-000  
(RTH+2RTHM) Folding Mechanism with tube mounts

### Main Component

1611-9101-000  
(UFC1000IP) Frame Clamp Inner Piece 1" Round

*Refer to Spare Parts for other frame clamps options.  
Requires full access to the round tube*

1611-9064-000  
(S-22x3) Standard S-Bend Tube 22" length x 3" offset

*Refer to Spare Parts for other S-Bend Tube sizes  
or custom tubes*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

## For Power Wheelchairs

1611-0078-000  
(RTH+2RTHM) Folding Mechanism with tube mounts

### Main Component

1601-0020-M6SS  
Universal T-Nut/+Nut for Power Wheelchair Seat Rails.  
Fits Quickie Q series, Invacare TDX, Quantum

*Refer to Spare Parts for other TNUT options*

1611-0049-000  
(MH2/IPA) Multi-hole 2" Inner Piece Adapter

*Refer to Spare Parts for other seat rail attachment options*

1611-9064-000  
(S-22x3) Standard S-Bend Tube 22" length x 3" offset

*Refer to Spare Parts for other S-Bend Tube sizes  
or custom tubes*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

## For Fixed Wheelchairs

1611-0122-000  
(UFCOP) Frame Clamp Outer Piece

1611-0034-000  
(IC1) Index Clamp

*UFCOP and IC1 are required.*

1611-0128-000  
(USBF) Folding Quick Release Base  
*Refer to Spare Parts for other Quick Release Bases*

## For Tilting Wheelchairs

1611-0071-000  
(RFCR) Removable Frame Clamp Receiver

1611-0076-000  
(ROP) Removable Outer Piece

*RFCR and ROP are required*

1611-9008-000  
(UDS-22-DC) Rehadapt  
*Therapist choice*

# Swing Aside Mount

For greater independence the user may release the lock pin to swing the mount away. It's easily removed from wheelchair and consists of two lengths of stainless-steel tube joined to form a right-angle structure, that is supported on the wheelchair by a frame clamp assembly, which is permanently attached to either the left or right side of the wheelchair. The design of this mount allows use on a tilt seat system.



Right Angle Tube Option

Straight Vertical S-Tube Option

## For Manual Wheelchairs

1611-9101-000  
(UFC1000IP) Frame Clamp Inner Piece 1" Round

*Refer to Spare Parts for other frame clamps options.  
Requires full access to the round tube*

## For Power Wheelchairs

1601-0020-M6SS  
Universal T-Nut/+Nut for Power Wheelchair Seat Rails.  
Fits Quickie Q series, Invacare TDX, Quantum

*Refer to Spare Parts for other TNUT options*

## Right Angle Tube Option

1611-0045-000  
(LSAOP) Locking Swing Away Outer Piece

1611-0035-000  
(ICLSA) Index Clamp for Locking Swing Away

*LSAOP and ICLSA are required.  
Not recommended for wheelchairs that has Tilt Function*

1611-0148-000  
Release Cable Assembly with 90° Tube

1611-0049-000  
(MH2/IPA) Multi-hole 2" Inner Piece Adapter

*Refer to Spare Parts for other seat rail attachment options*

## Straight Vertical S-Tube Option

1611-0046-000  
(LSARH/TM/16) Complete Swing-Aside Arm consisting of a Rotate Head, Tube Mount, 16" Tube and Release Cable for the Swing-Aside Mount

1611-0071-000  
(RFCR) Removable Frame Clamp Receiver

1611-0076-000  
(ROP) Removable Outer Piece.

*RFCR and ROP are required.  
Recommended for wheelchairs that has Tilt Function*

1611-0054-000  
(O3L) Offset Link Offsets the Outer Piece 3" away from the Frame Clamp Inner Piece  
*Optional Part*

1611-9064-000  
(S-22x3) Standard S-Bend Tube 22" length x 3" offset  
*Refer to Spare Parts for other S-Bend Tube sizes or custom tubes*

1611-9065-010  
Daessy Custom Tube.  
*Please advise Tube Size – For Eyegaze device*

S-Bend Tube:  
L = Overall length  
O = Offset S-bend  
S = Straight before bend



1611-0128-000  
(USBF) Folding Quick Release Base  
*Refer to Spare Parts for other Quick Release Bases*

1611-9008-000  
(UDS-22-DC) Rehadapt  
*Therapist choice*

# Rear Folding or Lockable Rear Folding Mount

Daessy Rear Folding Mount or Lockable Rear Folding Mount folds from front to behind, protecting the device and remains attached to the wheelchair. The Lock Mechanism allows this Mount assembly to be used with a tilting seat system. This set up is directional and the attachment side of the wheelchair must be specified at time of order to receive the correct Lock Mechanism.



## For Manual Wheelchairs

1611-0026-000  
(ELRH/RTHTM) Elbow Rotate Head and Tube Mount  
*Main component of Rear or Lockable Folding Mount*

1611-9101-000  
(UFC1000IP) Frame Clamp Inner Piece 1" Round  
*Refer to Spare Parts for other frame clamps options.  
Requires full access to the round tube*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

1611-0094-000  
(STR-22) Standard Straight Tube 22" Length 7/8" OD

## For Power Wheelchairs

1611-0026-000  
(ELRH/RTHTM) Elbow Rotate Head and Tube Mount  
*Main component of Rear or Lockable Folding Mount*

1601-0020-M6SS  
Universal T-Nut/+Nut for Power Wheelchair Seat Rails.  
Fits Quickie Q series, Invacare TDX, Quantum  
*Refer to Spare Parts for other TNUT options*

1611-0049-000  
(MH2/IPA) Multi-hole 2" Inner Piece Adapter  
*Refer to Spare Parts for other seat rail attachment options*

1611-0093-000  
(STR-16) Standard Straight Tube 16" Length 7/8" OD

1611-0094-000  
(STR-22) Standard Straight Tube 22" Length 7/8" OD

## For Removable Option

1611-0071-000  
(RFCR) Removable Frame Clamp Receiver  
1611-0070-000  
(RFCI) Removable Frame Clamp Insert  
*RFCR and RFCI are required*

## Non-Lockable Option

1611-0069-000  
(RFA+RTHTM) Rear Folding Adapter and Tube Mount  
*Not recommended for wheelchairs that has Tilt Function*

## Lockable Option

1611-0044-000  
(LRFA+RTHTM+O3L) Right  
1611-0159-000  
(LRFA+RTHTM+O3L) Left

1611-0054-000  
(O3L) Offset Link Offsets the Outer Piece 3" away from the Frame Clamp Inner Piece  
*Optional Part*

1611-0128-000  
(USBF) Folding Quick Release Base  
*Refer to Spare Parts for other Quick Release Bases*

1611-9008-000  
(UDS-22-DC) Rehadapt  
*Therapist choice*

# Positioner Mount

Achieve the ideal position and angle through 6 points of adjustment. Easily removed from the wheelchair, it can be used on any style of wheelchair or alternate location. The Removable Frame Clamp Receiver has a locking mechanism that can be engaged to secure the mount for use on a tilting seat system.



## Standard Mounting Assembly

### For Manual Wheelchairs

1611-0320-000

(DPM9) Positioner Mount Default

*Complete Assembly*

**Components:** UFCxxxxIP, AHOP, RFCR, MTDPM9, Stainless Steel Tubes & TUSB

### For Power Wheelchairs

1611-0330-000

(DPM9-NFC) Positioner Mount excluding a frame clamp assembly

**Components:** AHOP, MTDPM9, Stainless Steel Tubes & TUSB

1601-0020-M6SS

Universal T-Nut/+Nut for Power Wheelchair Seat Rails.

Fits Quickie Q series, Invacare TDX, Quantum

*Refer to Spare Parts for other TNUT options*

1611-0049-000

(MH2/IPA) Multi-hole 2" Inner Piece Adapter

*Refer to Spare Parts for other seat rail attachment options*

1611-0054-000

(O3L) Offset Link Offsets the Outer Piece 3" away from the Frame Clamp Inner Piece

*Optional Part*

1611-0128-000

(USBF) Folding Quick Release Base

*Refer to Spare Parts for other Quick Release Bases*

1611-9008-000

(UDS-22-DC) Rehadapt

*Therapist choice*

# Desk and Table Mounts

## Desk Mount

1611-0020-000  
(DM)

The original Daessy Desk Mount featuring an angle adjustable U-Tube and choice of Quick-Release Base. The Desk Mount is most commonly used in conjunction with one of the Daessy wheelchair Mounts and allows a device to be positioned on a desk or table at precisely the same orientation as it is on the wheelchair. The Desk Mount consists of a base and U-shaped tube that holds the attachment fitting for the device on its adapter.



## Clamp Desk Mount

1611-0022-000  
(DMC)

This variation on the Desk Mount uses a Manfrotto Super Clamp to clamp onto the edge of any suitable surface between 15mm and 50mm thick. Similar to the standard Desk Mount, the U-Tube can be adjusted for angle and position.

It can be supplied with choice of Quick-Release Base to customer preference.



## Vertical Desk Stand with Quick Handles

1611-0238-000  
(VDS-2H)

The Daessy Vertical Desk Stand provides a solution for larger devices and taller placement to use on a desk or table surface. The Articulating Quick Release Base, included by default, allows infinite adjustment of the angle and rotation of the attached device or computer. The Articulating Quick Release Base position can be set anywhere along the vertical tube for height adjustment. The Vertical Desk Stand disassembles for compact storage or transportation.

### Standard configurations:


- VDS Vertical Desk Stand (AQRL-87)
- VDS-RH Vertical Desk Stand with Articulating Quick Release Base quick handle (AQRL-87-RH)
- VDS-2H Vertical Desk Stand with Height Adjustment and Articulating Quick Release Base quick handles (AQRL-87-2H)



## X Base Folding Desk Stand

1611-0132-000  
(XDS)

The X-Base Folding Desk Stand provides a sturdy vertical desk stand support which can be folded down easily and compactly. The Articulating Quick Release Base, included by default, allows infinite adjustment of the angle and rotation of the attached device or computer. The Articulating Quick Release Base position can be set anywhere along the vertical tube for height adjustment. The X-Base Folding Desk Stand easily folds for compact storage or transportation. Adjustable, solid construction and portable.




# Rolling Mounts


The Daessy Rolling Mount provides an independent floor stand solution for any AAC device or laptop. Fitted with the Daessy quick release base the Rolling Mount is an excellent companion to the Daessy wheelchair mounts allowing easy transfer of an AAC device or laptop from wheelchair to the Rolling Mount for use while in bed or an alternate location.

Daessy Rolling Mounts are now available in a variety of sizes and configurations.

- Independent stand support
- Height-adjustable
- Locking caster wheels

Rolling Mounts (Floor Stands) Includes TUSB quick release base, unless otherwise noted.

Offset Base	
1611-0227-000 (ROM24) 24" Base Std Height	
1611-0229-000 (ROM32) 32" Base Std Height	
1611-0072-000 (ROM32T) 32" Base Tall Height	
1611-0179-000 (ROM36) 36" Base Std Height	
1611-0074-000 (ROM36T) 36" Base Tall Height	

Standard Base	
1611-0177-000 (RM24) 24" Standard Base	
1611-0178-000 (RM32) 32" Standard Base	
1611-0241-000 (RM32T) 32" Standard Base Tall	

Articulated Joint Horizontal Tube (AJH)	
1611-0073-000 (ROM32T-AJH-AQL) 32T with AQRL-87	1611-0242-000 (ROM32T-NQR) Offset 32" Base, Exl TUSB
1611-0075-000 (ROM36T-AJH-AQL) 36T with AQRL-87	1611-0243-000 (ROM32T-AJH-NQR) Offset 32" Base Articulating Horizontal



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY: