


Hero-NXT



Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Basic Packages

Includes:	Package Sizes
Available in Black/Red 10° Tilt-in-space Swivel base 360° Rear/Forward facing position Five-point harness Backrest Height Adjustable Adjustable Lateral side Supports	One Size 4181-0010-020



Specifications

Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Crown Height	Distance Laterals	Total Height	Total Width	Total Depth	Capacity
280mm	210mm	250mm	210mm - 340mm	410mm - 530mm	120mm - 210mm	660mm - 730mm	430mm	630mm	18kg

*measurements with standard cushion

Tray Table

Padded
One Size 4181-5530-000



medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Hero-NXT_20221221_AU