


# Hero-NXT




# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
<b>User Measurements</b>	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



# Basic Packages


Includes:	Package Sizes
<p>Available in Black/Red                      10° Tilt-in-space                      Swivel base 360°                      Rear/Forward facing position                      Five-point harness                      Backrest Height Adjustable                      Adjustable Lateral side Supports</p>	<p><b>One Size</b>                      4181-0010-020</p> 

# Specifications

Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Crown Height	Distance Laterals	Total Height	Total Width	Total Depth	Capacity
280mm	210mm	250mm	210mm - 340mm	410mm - 530mm	120mm - 210mm	660mm - 730mm	430mm	630mm	18kg

\*measurements with standard cushion

# Tray Table

<b>Padded</b>
<p><b>One Size</b>                      4181-5530-000</p> 

# Lateral Supports

# Abduction Block

Pair (Recommended)



For Hero-NXT



One Size  
4181-3310-000



One Size  
4181-1510-000



# Seat Belt Extension

5-Point\*



One Size  
4081-6530-000



\*available as special order only

# Seat Protection Mat

For Hero-NXT



One Size  
4181-0190-000



# Notes

---

---

---

---

---

---

---

---

---

---

# medifab

Distributed by Medifab Ltd

22 Detroit Drive  
Rolleston 7675  
New Zealand

+64 3 307 9790

[solutions@medifab.com](mailto:solutions@medifab.com)

[www.medifab.com](http://www.medifab.com)



DISTRIBUTED BY:

Script: Hero-NXT\_20221221\_NZ