


# Kidsflex



# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



# Basic Packages

Size 1	Size 2	Size 3
<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>• Colour</li> <li>• 5 Point Harness</li> <li>• Recline Function</li> <li>• Tilt angle adjustment (0 to 20 degs)</li> <li>• Slim left turn swivel base</li> <li>• Seat belt guide on the left</li> </ul> <p>*right turning available as special order</p>	<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>• Colour</li> <li>• 5 Point Harness</li> <li>• Recline Function</li> <li>• Tilt angle adjustment (0 to 20 degs)</li> <li>• Slim left turn swivel base</li> <li>• Seat belt guide on the left</li> </ul>	<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>• Colour</li> <li>• 5 Point Harness</li> <li>• Recline Function</li> <li>• Tilt angle adjustment (0 to 20 degs)</li> <li>• Slim left turn swivel base</li> <li>• Seat belt guide on the left</li> </ul>
4170-0011-020 <input type="radio"/>	4170-0021-060 <input type="radio"/>	4170-0031-090 <input type="radio"/>


# Specifications

Size	Total Height	Seat to Head	Back Height	Seat Width	Seat Depth	Lateral Sizes	Distance of Lateral Supports	Weight	Max User Weight
2	80 cm	73 cm	41-57 cm	35 cm	34-39 cm	26-36 cm	19-32 cm	8.9 kg	36 kg
3	90 cm	82 cm	47 - 65cm	40 cm	40-46 cm	38-48 cm	24-37 cm	12.3 kg	75 kg

Required: Stability Option Must be Selected

# ISOfix Adapter\*

# Clamp Bow

Select Size		To be used if ISOfix is not available.		
<b>Size 1</b> 4170-6911-100	<input type="radio"/>		<b>Size 1</b> 4170-0261-000	<input type="radio"/>
<b>Size 2/3</b> 4170-6912-100	<input type="radio"/>		<b>Size 2/3</b> 4170-0262-000	<input type="radio"/>

\*Only compatible with vehicles fitted with ISOFIX receivers

# Foot Support

# Shoulder Supports

Size	1 Size
<b>Size 1</b> 4170-5241-000	
<b>Size 2/3</b> 4170-5242-000	
<b>4170-2750-000</b>	

# Head Supports

Lateral Head Support (Recommended)	Extended Head Support
<b>Size 1</b> 4170-4201-000	<b>Size 1</b> 4170-4301-000
<b>Size 2</b> 4170-4202-000	<b>Size 2</b> 4170-4302-000
<b>Size 3</b> 4170-4203-000	<b>Size 3</b> 4170-4303-000

# Lateral Supports

Small Lateral Supports (Recommended)	Large Lateral Supports
<b>Fixed</b> 4170-3301-000	<b>Fixed</b> 4170-3302-000
<b>Swing-Away</b> 4170-3401-000	<b>Swing-Away</b> 4170-3402-000



# medifab

Distributed by Medifab Ltd

22 Detroit Drive  
Rolleston 7675  
New Zealand

+64 3 307 9790

[solutions@medifab.com](mailto:solutions@medifab.com)

[www.medifab.com](http://www.medifab.com)



DISTRIBUTED BY:

Script: Kidsflex\_20230126.NZ