

Stroller Comparison Chart



SPECIFICATIONS	BINGO MINI	BINGO EVO			BINGO SPEX		PIXI		DURO		KOOL STRIDE
		Size 1	Size 2	Size 3	Size 1	Size 2	Size 1	Size 2	Size 1	Size 2	
DESCRIPTION	Ultimate early intervention seating with smart outdoors stroller	Versatile and modular buggy system for basic needs, up to critical care users.			A smart-looking stroller with superior postural support, that enables parents to take their child outside and explore!		Light and nimble buggy with comfortable rigid seating. Child can step on footplate to self-transfer.		Configurable stroller for everyday challenges, with a great range of support and positioning accessories.		Go off road and beyond!
APPROX AGE OF USER*	4mths-4yrs	1-5yrs	3-12yrs	8-15yrs	2-5yrs	3-7yrs	3-7yrs	6-14yrs	3-12yrs	8-15yrs	8-21yrs
MAX LOAD (kg)	25	50			50		35	50	75		70
SEAT DEPTH (cm)	13-25	18-30	23-35	23-35	17-28	22-33	28	33	23-34	28-39	43
SEAT WIDTH (cm)	12-27	16-34	22-40	30-45	18-30	23-38	25 or 30	35 or 40	26-40	31-45	30
SHOULDER HEIGHT (cm)	25.5-39.5	28-46	34-55	42-63	36-44	40-48	54 or 60	64 or 70	32-40	42-62	60
FOOTPLATE HEIGHT RANGE (cm)	16-34	19-40	19-40	19-40	25-40	25-45	13-35	13-39	20-44	20-44	45
FRONT/REAR FACING	Both	Both			Both		Front		Front		Front
INDOOR FRAME OPTIONS	✓	✓			✓		✗		✗		✗
APPROX PRODUCT WEIGHT (kg)	17.5	15	16.1	18	29	30	12.5	15	19.2	21	15
TILT IN SPACE	45°	45°			45°		✗		Static 16° or 20°		✗
RECLINE	80°-180° in 10° steps	90°-160°	90°-135°		90°-160°		90°-125°		78°-120°		✗
FOOTPLATE ELEVATION	✓	✓			✓		✗		✓		✗
GROWTH ADJUST SEAT	✓	✓			✓		5cm growth bracket, seat depth & back height		✓		✗
SWIVEL FRONT WHEELS	✓	✓			✓		✓		✓		✗
ACTIVITY TRAY	✓	✓			✓		✓		✓		✗
SUN CANOPY	✓	Yes			✓		✓		✓		✓
RAIN SHIELD	✓	✓			✓		✓		✓		✓
BASKET STORAGE	✓	✓			✓		✓		✓		✓
POSTURAL SUPPORT LEVEL	Mod./Complex	Mod./Complex			Complex		Moderate		Moderate		Mild
TYRE SOLID / PNEUMATIC	Solid	Solid			Solid		Combination		Solid		Pneumatic
TRANSPORT APPROVED	✓	✓			✓		✓		✓		✗

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.