

# Pixi Stroller



**RIGID &  
LIGHTWEIGHT**

Distributed by  
**medifab**

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE


# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Measurements	Size 1	Size 2
Pixi Order Code	2563-0003-191	2566-0003-191
Max. User Weight	35kg	50kg
Max. User Height	100cm (120cm with extension)	125cm (150cm with extension)
A Seat Depth	25cm (30cm with extension)	35cm (40cm with extension)
B Seat Width	28cm	33cm
C Shoulder Height	36-45cm	45-55cm
D Footplate Height	13-35cm	13-39cm
E Overall Width	58.5cm	63.5cm
F Push Handle Height	57-110cm	62-120cm
Stroller Weight (approx)	12.5kg	15kg
Folded Dimensions	L: 73cm W: 59cm H: 40cm	L: 85cm W: 66cm H: 51cm
Back Angle	90-125°	90-123°
Basket Weight Limited	5kg	5kg

**Stroller Package includes:**




- Foldable Frame
- Sun Roof
- Grey Upholstery
- Flip-up Foot Plate
- 2 Point Hip Belt











# Lateral Head Support

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
<b>Merino/Spacer Fabric</b> 2010-4110-099		<b>Small (17cm Internal width)</b>		<b>Small (17cm Internal width)</b>	
		<b>Merino/Spacer Fabric</b> 2010-4210-099		<b>Merino/Spacer Fabric</b> 2010-4310-099	
		<b>Wipe Down Cover</b> 2010-4210-209		<b>Wipe Down Cover</b> 2010-4310-209	
<b>Wipe Down Cover</b> 2010-4110-209		<b>Large (21cm Internal width)</b>		<b>Large (21cm Internal width)</b>	
		<b>Merino/Spacer Fabric</b> 2010-4211-099		<b>Merino/Spacer Fabric</b> 2010-4311-099	
		<b>Wipe Down Cover</b> 2010-4211-209		<b>Wipe Down Cover</b> 2010-4311-209	

# Growth Adjustments

Infill Spacer Block		Seat Depth Extension Bracket		Foot Plate Spacer	
To decrease seat depth. Specify thickness: mm		To increase seat depth by 5cm		<b>Size 1</b> To decrease footplate height 2560-5301-299	
2560-2800-032		<b>Size 1</b> 2560-0281-199			
		<b>Size 2</b> 2560-0282-199			


# Harnesses & Hip Belts

Pixi 5-Point Harness		Spex 5-Point Harness	
With front pull adjusters and camlocks 2560-6511-000		On tee pad for babies and very small children 2006-6551-000	
Anti-Escape Tee Pad Harness		Chest Strap with Lateral Pads	
Fitted with anti-escape buckles to prevent user from escaping stroller. 2560-6556-026		<b>Size 1</b> 2560-6801-090	
		<b>Size 2</b> 2560-6802-090	
2-Point Padded Centre Pull Hip Belt			
<b>Centre-Release Buckle</b> 1405-6211-028	<b>Anti-Escape Buckle</b> 1405-6211-026	<b>Pen Popper Buckle</b> 1405-6211-025	
			

# Abduction Pommel

To keep user's legs in the correct position

2560-1501-090



# Tray Table

Activity and Feeding Tray

**Size 1**

2560-7120-013

**Size 2**

2560-7130-013



# Grab Rail

Attaches to the Pixi frame

**Size 1**

2560-7253-299

**Size 2**

2560-7254-299



# Weather Protection

Rain Shield Cover

For protection from the rain  
2560-7601-010



UV Shade Cover

For protection from the sun  
2560-7701-010



# Storage Options

Mesh Basket with Plastic Bottom (weight limit 5kg)

**Size 1**

2560-8011-017

**Size 2**

2560-8012-017



# Tie-Downs

ISO 7176-19

2560-8501-021



# medifab

Distributed by Medifab Ltd

22 Detroit Drive  
Rolleston 7675  
New Zealand

0800 543 343

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Pixi\_20230705\_NZ