


Starlight-NXT




Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Foot Support Height	
E Chest Width	
F Axilla Height	
User Weight	




Basic Packages

Includes:	Package Style
<ul style="list-style-type: none"> Adjustable head support Back support infill padding Slim lateral supports 20° tilt Base board accessory Swivel base system - (Left turning) to enable turning of the car seat towards the doorway ISOFIX seat attachment Grey/Black upholstery 	<p>Left Turning 4141-0021-090</p> 

Specifications

	Seat Width (pelvis)	Seat Width (knee)	Seat Depth	Back Height	Crown Height	Lateral Height Adjustment	Lateral Width Adjustment	Product Weight	Weight Capacity
Standard	30cm	35cm	26cm	32–42cm	51–62cm	17–27cm	14–19cm	8.5kg	9–36kg
XL	30cm	35cm	32cm	44–53cm	67–76cm	25–38cm	14–19cm	8.5kg	9–36kg

XL Upgrade Kit

Includes:	
<ul style="list-style-type: none"> Extended Harness Longer back support padding Additional mounting plates <p>(Harness upgrade only possible when child is over 15kg)</p>	<p>4141-0140-000</p> 

medifab

Distributed by Medifab Ltd

26 Pardoe Street,
East Devonport,
Tasmania, 7310

+61 3 6427 0105

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Starlight-NXT_20230417_AU