

Symmetrisleep



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client details

Client/User name:	
Date:	
DOB:	

Contents:

- Symmetrisleep Positioning System
- Symmetrisleep Thermal Regulating Options
- Lowzone Posture Cushions

Velco Receptor Sheet



Stretchy velcro-like fitted sheet for securing of positioning brackets and other supports

Cot 7010-8005-CT0 <input type="radio"/>	3/4 7010-8000-T00 <input type="radio"/>	Single 7010-8005-SN0 <input type="radio"/>	Single King 7010-8005-KS0 <input type="radio"/>	Double 7010-8005-D00 <input type="radio"/>
Length 100 cm Width 75 cm	Length 160 cm Width 80 cm	Length 190 cm Width 100 cm	Length 205 cm Width 106 cm	Length 205 cm Width 150 cm

Positioning Brackets

Brackets that attach to Velco Receptor Sheet

Small 7011-7240-00B <input type="radio"/>	Medium 7011-7241-00B <input type="radio"/>	Large 7011-7242-00B <input type="radio"/>	Extra Large 7011-7230-001 <input type="radio"/>	
Colour: Red	Colour: Yellow	Colour: Blue	Colour: Green	
Width 10 cm Height 8 cm Length 13cm	Width 16 cm Height 13 cm Length 13 cm	Width 15 cm Height 25 cm Length 25 cm	Width 15 cm Height 28 cm Length 25 cm	


Coolover Bracket Pads

Padding to be installed on plastic brackets


Small 7012-7165-001 <input type="radio"/>	Medium 7012-7175-001 <input type="radio"/>	Large 7012-7185-001 <input type="radio"/>	Extra Large 7012-7195-001 <input type="radio"/>	
Length 13 cm Width 10 cm Height 8 cm	Length 13 cm Width 16 cm Height 13 cm	Length 25 cm Width 15 cm Height 25 cm	Length 25 cm Width 15 cm Height 28 cm	

Monoblok Knee Supports

Monoblok Supine Knee Support


Small 7013-7552-000 <input type="radio"/>	Medium 7013-7553-000 <input type="radio"/>	Large 7013-7554-000 <input type="radio"/>	
Length 10 cm Width 36 cm	Length 20 cm Width 40 cm	Length 25 cm Width 60 cm	

Monoblok Rolling Knee Support


Small 7013-7562-000 <input type="radio"/>	Medium 7013-7563-000 <input type="radio"/>	Large 7013-7564-000 <input type="radio"/>	
Length 18 cm Width 10 cm	Length 20 cm Width 21 cm	Length 25 cm Width 27 cm	

Monoblok Knee Support Accessories

Coolover Legwrap Liners (suits either of Monoblok supports)

	Small (8cm Diameter)	Medium (8cm Diameter)	Large (16cm Diameter)	
5mm Code	7013-7401-CVR <input type="radio"/>	7013-7421-CVR <input type="radio"/>	7013-7421-CVR <input type="radio"/>	
10mm Code	7013-7411-CVR <input type="radio"/>	7013-7412-CVR <input type="radio"/>	7013-7413-CVR <input type="radio"/>	

Monoblok Ramp Wedges


Medium 7013-7573-000 <input type="radio"/>	Large 7013-7575-000 <input type="radio"/>	
Length 40 cm Width 60 cm	Length 50 cm Width 70 cm	

Monoblok Safety Straps to secure to mattress

Single Strap 7013-7585-000 <input type="radio"/>	Double Strap 7013-7586-000 <input type="radio"/>	
---	---	--


Side Lying Leg Support

Stabilises the pelvis in a side-lying position

Small 7013-7445-000 <input type="radio"/>	Medium 7013-7450-000 <input type="radio"/>	Large 7013-7455-000 <input type="radio"/>	
Depth 36 cm Width 36 cm	Depth 60 cm Width 60 cm	Depth 74 cm Width 74 cm	

Airmantle

Mattress topper for draping over brackets and supports for air circulation/thermal comfort

Single 7014-8010-S00 <input type="radio"/>	
Length: 190cm Width: 100cm Depth: 2cm	

Coolover Sheet

Temperature regulation sheets



Cot 7015-7120-C00 <input type="radio"/>	3/4 7015-7120-T00 <input type="radio"/>	Single 7015-7120-S00 <input type="radio"/>	Single King 7015-7120-SK0 <input type="radio"/>	Double/Queen 7015-7120-D00 <input type="radio"/>	Flat Sheet 7015-7120-SFS <input type="radio"/>
Length 130 cm Width 60 cm	Length 160 cm Width 80 cm	Length 190 cm Width 106 cm	Length 210 cm Width 105 cm	Length 215 cm Width 150 cm	Length 320 cm Width 210 cm

Symmetrisleep Stickers

Provides the carers with customised user instruction sheets












- Sheets allow you to choose a lying position and then select the parts you need, to assist with planning and recording sleep system set-up.
- Add the selected model to the chart, and accessories can be placed underneath as well as on top, to mirror the sleep system set-up of the client.
- Space to add name and instructions.

7010-0061-000



Lowzone Posture Cushions

Wipe-down PU fabric pillows with chipped memory foam filling accessible via zip opening

A: Square 7016-7500-001 <input type="checkbox"/>	C: Med Oblong 7016-7580-001 <input type="checkbox"/>	D: Lrg Oblong 7016-7610-001 <input type="checkbox"/>	ES: Sml Horseshoe 7016-7720-001 <input type="checkbox"/>
Length 30 cm Width 30 cm	Length 45 cm Width 30 cm	Length 75 cm Width 45 cm	Length 37 cm Width 25 cm
			
Coolover Cover 7016-7500-CVR <input type="checkbox"/>	Coolover Cover 7016-7580-CVR <input type="checkbox"/>	Coolover Cover 7016-7610-CVR <input type="checkbox"/>	Coolover Cover 7016-7720-CVR <input type="checkbox"/>
Bamboo Cover 7016-7500-BAM <input type="checkbox"/>	Bamboo Cover 7016-7580-BAM <input type="checkbox"/>	Bamboo Cover 7016-7610-BAM <input type="checkbox"/>	Bamboo Cover 7016-7720-BAM <input type="checkbox"/>
EM: Med Horseshoe 7016-7710-001 <input type="checkbox"/>	EL: Lrg Horseshoe 7016-7700-001 <input type="checkbox"/>	EX: Horseshoe Extended 7016-7702-001 <input type="checkbox"/>	F: Sml Roll 7016-7800-001 <input type="checkbox"/>
Length 65 cm Width 30.5 cm	Length 80 cm Width 35 cm	Length 65 cm Width 30 cm	Length 30 cm Width 8 cm
			
Coolover Cover 7016-7710-CVR <input type="checkbox"/>	Coolover Cover 7016-7700-CVR <input type="checkbox"/>	Coolover Cover 7016-7702-CVR <input type="checkbox"/>	Coolover Cover 7016-7800-CVR <input type="checkbox"/>
Bamboo Cover 7016-7710-BAM <input type="checkbox"/>	Bamboo Cover 7016-7700-BAM <input type="checkbox"/>	Bamboo Cover 7016-7702-BAM <input type="checkbox"/>	Bamboo Cover 7016-7800-BAM <input type="checkbox"/>
G: Med Roll 7016-7850-001 <input type="checkbox"/>	G+: Med Roll Extended 7016-7860-001 <input type="checkbox"/>	Z Pillow 7016-7810-001 <input type="checkbox"/>	
Length 75 cm Width 45 cm	Length 37 cm Width 25 cm	Length 65 cm Width 30.5 cm	
			
Coolover Cover 7016-7850-CVR <input type="checkbox"/>	Coolover Cover 7016-7860-CVR <input type="checkbox"/>	Coolover Cover 7016-7810-CVR <input type="checkbox"/>	
Bamboo Cover 7016-7850-BAM <input type="checkbox"/>	Bamboo Cover 7016-7860-BAM <input type="checkbox"/>	Bamboo Cover 7016-7810-BAM <input type="checkbox"/>	

medifab

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Symmetrisleep 20221221_NZ