

CareWave Details

- Body positioning system which offers lying and seating support to infants and adults.
- Various cushion shapes are available depending on the user's need for positioning.
- Unique micro-particle padding provides flexibility and a moulding effect, essential to support asymmetric postures, fill body gaps and for comfort around the user's body and limbs. Available in two sizes.
- Outer cover of CareWave cushions: Waterproof and breathable Polyamid bi-elastic fabric (47%) coated with polyurethane (53%) Coating with antibacterial treatment.
- Filling: 100% Polystyrene micro-particles.
- Fire standard: Passes the requirements of EN 597 - 1&2 and EN ISO 12952 - 1 to 4.



COLOUR IDENTIFICATION

GREEN CAREWAVE = Small sizing
 BLUE CAREWAVE = Large sizing



Cushion Covers

- Cushion covers can be added to most sleep system items. Medifab suggest that at least one cover is purchased with each cushion. Benefits include:
 - Quick and easy to launder.
 - Breathability, allows skin to release moisture without undue discomfort or pressure marking
 - Warm to touch, adds comfort.
 - Absorbent - wicks away perspiration, saliva etc.



Cleaning Instructions

- Cushion Maintenance:
 Wipe with usual detergent.



CareWave Kit



POSTURAL CARE MADE EASY

- For supine, prone or side lying.
- Cylindrical bolsters and abduction wedge are Velcro® attached for additional postural control.
- Easily transported.
- Easy to set up.

CAREKIT INCLUDES

- Micro-particle positioning mattress.
- 2 Cylindrical bolster cushions.
- 1 Small half-moon cushion.
- 1 Abduction cushion.
- 1 Small Universal cushion.
- Transport bag.

	Code	Mattress Length
Small	7007-0001-066	120 cm
Large	7007-0002-066	160 cm



CareWave Delta

- Create lateral stability either alone, with one on each side, or an extra pillow to create a neutral position in the knee & hip joints.
- Supports and prevents the patient from falling forward from a decubitus position.
- Can be used in a semi fowler position, by placing under the legs.

	Small	Large	XX Large
Code	7007-0552-045	7007-0551-066	7007-0554-066
Width	40 cm	55 cm	180 cm
Depth	20 cm	30 cm	30 cm
Height	20 cm	30 cm	30cm
Cover Code	7007-0552-270	7007-0551-270	7007-0554-270



CareWave Universal

- Use to avoid rubbing between knees, calves or heels.
- Placed under the feet, it helps to reduce pressure on heels.
- Can also be used simply as a pillow.

	Small	Large
Code	7007-1252-045	7007-1251-066
Length	45 cm	60 cm
Width	30 cm	40 cm
Cover Code	7007-1252-046	7007-1251-067



CareWave Cylindrical

- Can be used between legs, under heels and ankles in order to reduce pressure.
- Prevents the user from falling forward from a decubitus position.
- Provides sufficient support to achieve a corrected position for many users.
- Use for corrected side lying positioning for sufficient angles to keep the airway open and to ensure a neutral position in the involved joints.

	Small	Large
Code	7007-1502-045	7007-1501-066
Length	45 cm	70 cm
Diameter	14 cm	17 cm
Cover Code	7007-1502-270	7007-1501-270



CareWave Abduction

- Use when the patient is lying down, placed in crotch area.
- Can also be placed under the arm.

	Small	Large	XX Large
Code	7007-1552-045	7007-1551-066	7007-1554-066
Width	12 cm	15 cm	90 cm
Depth	23 cm	30 cm	95 cm
Height	15 cm	20 cm	30 cm
Cover Code	7007-1552-270	7007-1551-270	



CareWave Semi-Fowler

- Gives support and holds the user comfortably in a semi-fowler position.
- Place the cushion under the legs, the abduction module forward between the knees.

	Small	Large
Code	7007-1852-045	7007-1851-066
Width	50 cm	70 cm
Depth	35 cm	50 cm
Height	12 cm	14 cm
Cover Code	7007-1852-270	7007-1851-270



CareWave Decubitus

- Supports shoulders, back, hips and knees of the user in a decubitus position. Reduces rubbing.
- Provides excellent trunk stability while the hip and knee angle can be varied.
- Lower part of the "S" can be drawn up between the knees so that the right median and good body symmetry can be achieved.
- Can also be shaped into a U-shaped cushion by adjusting the fill and rotating one end to thus provide a neutral position for the knee and hip joints, while the pad is pushed under the seat to fix the pelvis.
- The ends also provide support for the elbows.

	Small	Large	XX Large
Code	7007-3002-045	7007-3001-066	7007-3004-066
Width	25 cm	30 cm	36 cm
Length	140 cm	180 cm	200 cm
Cover Code	7007-3002-270	7007-3001-270	



CareWave Decubitus Slim

- Designed for positioning of lower limbs to prevent them from crossing. Also prevents lateral sliding.

	Small	Large
Code	7007-3012-045	7007-3011-066
Width	9 cm	12 cm
Length	224 cm	300 cm
Cover Code	7007-3012-270	7007-3011-270



CareWave Back Cushion

- Maintains the back in a semi-sitting position when user is sitting in bed or in a lateral decubitus position. Helps reduce pressure on ischium and heels.
- Can also be used in a sitting position in a chair.

	Small	Large
Code	7007-2452-045	7007-2451-066
Length	200 cm	200 cm
Width	40 cm	60 cm
Cover Code	7007-2452-270	7007-2451-270



CareWave Body Bolster

- New addition to the CareWave range, which allows full support of the body from head to feet.
- An integrated pillow facilitates installation and user comfort and successfully maintains user position.
- When used in lateral position, the Body Bolster protects against pressures and avoids friction and pressure points on arms, knees and calves.

	Small	Large
Code	7007-3502-045	7007-3501-066
Length	135 cm	210 cm
Width	20 cm	35 cm
Cover Code	7007-3502-270	7007-3501-270



CareWave Half-Moon

- It gives support and holds the patient comfortably in a lying or sitting position, and can be used to stabilize the upper body. The filling can be moved within the cushion so you can arrange it to adjust the head and arm height.
- If one wishes to counteract unwanted lateral movement of the body, the ends of the Half-Moon can be pushed under the hip area. Body weight will stabilize the pelvis and prevent unwanted movement. The idea is to eliminate movement that will take the body out of the desired symmetrical position.

	Small	Large
Code	7007-4052-045	7007-4051-066
Length	140 cm	180 cm
Width	25 cm	30 cm
Cover Code	7007-4052-270	7007-4051-270



CareWave Ring Cushion

- Comfortably supports the user's head. If support is required especially on one side, the special filling can be shaken to one side of the cushion and positioned slightly skewed to one side if needed.
- Can be positioned with the opening towards the patient's visual direction which allows for a stable position for the head and neck, backwards and sideways.
- If the patient has the ability to make a lot of head movement the pillow can move out of place so the patient must be able to move the pillow or get help with this.



	Small	Large
Code	7007-4082-045	7007-4081-066
Width	35 cm	40 cm
Cover Code	7007-4082-270	7007-4081-270

CareWave Heel Cushion

- Fixes to the mattress at the heel area with an adjustable strap.
- Supports ankles and relieves pressure on the heel zone.

	One Size
Code	7007-5201-066
Length	90 cm
Width	18 cm



CareWave Heel Bootie

- Designed to protect and prevent the feet/heels from injury and bedsores and to provide comfort whatever the patient's position.

	Small	Large
Code	7007-5302-045	7007-5301-067
Length of foot surface	25 cm	30 cm



Bow Cushion

- Can be used between the knees as an abduction cushion.
- Use under heels to reduce pressure.

	Small	Large
Code	7007-1602-045	7007-1601-066
Length	28 cm	34 cm
Width	19 cm	24 cm
Height	6 cm	10 cm
Cover Code	7007-1602-270	7007-1601-270



Stabilegs

- Stabilegs uses the same filling as the CareWave cushions.
- Positioning device has removable and adjustable modules.
- Abduction module (removable/adjustable) to limit friction and prevent leg-crossing in cases of spasticity.
- Adduction module (removable/adjustable) ensures alignment of lower limbs and pelvis stability.
- Heel relief and prevention equinus for prevention of bedsores in the heel area. End modules (removable/adjustable) are necessary to ensure pressure relief between the bed footboard and Stabilegs.



	One Size
Code	7007-1831-066
Length	58 cm
Width	55 cm
Height	12 cm



Fitleg Non-Slip

- Positioning device for proper orthopaedic alignment of lower limbs and full or partial relief on the heels.
- Prevents knee hyper-extension.
- Full relief on heels with the Non-Slip Module.
- High resiliency foam for stable positioning & optimized product life.
- Two distinct modules allow for customised patient positioning.



	Fitleg Non-Slip	Fitleg Module
Code	7007-1821-066	
Length	55 cm	50 cm
Width	58 cm	21 cm
Height	12 cm	10.5 cm



Fitleg

- Positioning device for proper orthopaedic alignment of lower limbs which prevents front sliding, friction and shear. Anatomic shape provides maximum comfort and stable positioning.
- Fitleg Evo has abduction module to limit friction around knee articulations and prevent leg-crossing in case of spasticity.
- Prevention of knee hyper-extension due to the anatomic shape.
- Viscoelastic foam combined with HR foam base to provide maximum comfort.
- Reduces pressure on heels area with support on calves and thigh.



	Fitleg Small	Fitleg Large	Fitleg Evo
Code	7007-1802-066	7007-1801-066	7007-1811-066
Length	58 cm	58 cm	58 cm
Width	40 cm	55 cm	55 cm
Height	12 cm	12 cm	12 cm

Ear Protection Cushion

- Comfortably supports the patient's head and relieves contact on the ear area in order to prevent bedsores.
- Comes standard with two cotton covers.

	One Size
Code	7007-4751-066
Length	30 cm
Width	30 cm



Floor Pod

- Floor positioner which can be used for supine and prone positioning and postural drainage.
- This positioning cushion adapts to the user's body shape and comfortably keeps them in a half-seated position.
- Filling is not free-flow like a bean bag; unique CareWave filling supports the user in a stable position.

	Small	Large
Code	7007-0402-066	7007-0401-066
Diameter	110 cm	130 cm
Cover Code	7007-0402-270	7007-0401-270



Hemi-Care Cushion

- Supports and holds the upper limbs in a stable position. Both parts of the Hemi Arm are fixed together with Velcro strips.
- Maintains the arm in a position which cannot slide and reduces pressure on elbow area so that there is no sensation of tear strength in the shoulder. For users with kyphosis, the Hemi-Care cushion helps maintain the back straight.

	Small	Large
Code	7007-5522-045	7007-5521-066
Length	45 cm	48 cm long
Width	30 cm	35 cm

