

# CentroGlide Bariatric<sup>250</sup>



## Australia's bariatric power chair.

We developed the CentroGlide Bariatric to cater for weights of up to 250kg with the same easy drive capabilities as the standard CentroGlide. Designed and manufactured by Glide in Australia, the CentroGlide is the result of decades of experience, user feedback and refinement.

### Completely customisable.

As a fully-scripted wheelchair, CentroGlide Bariatric is custom-built to the occupant's size and comfort requirements. The options are almost endless.

### Mobility and independence.

The powerful motor gives the user the freedom to travel across grass, uneven footpaths, and ramps. CentroGlide has great manoeuvrability in tight spaces, and stability over kerbs. It comes with 50° seat tilt as standard, with options for vertical lift and anterior reach to allow the users to reach high shelves or transfer from the chair in a standing position.

## Designed for comfort and ease of use.

We have put decades of user feedback and experience to develop all the features that make the CentroGlide Bariatric a fit for purpose, durable and robust chair.



### Customised seating.

Many seating options available. For ultimate comfort, choose your customised cushion, backrest and headrest.



### Active suspension.

Our patented 6-wheel APT suspension improves traction and aids ascending and descending kerbs safely without affecting the stability and angle of the seat.



### Wheels.

CentroGlide has 12" drive wheels in either standard width and wide format. The front castor are 8".



### Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



### Electronics.

Fully programmable R-net joystick control system. Compatible with many other 3rd party speciality electronics.



### Power seating.

As standard, the seat has a 50° tilt with other options available including power recline and power leg rest elevation.



### Battery range.

The powerful battery gives you a range of up to 40km on a single charge.



# CentroGlide sets the benchmark.

## Features.

- Power base with removable seating module
- High quality German 4 pole high torque motors
- Vehicle tie down points
- Complex electronics available
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) suspension
- 90mm ground clearance
- Max ramp incline/decline 12° (as per AS3695 testing)
- Drive wheels 12"
- Electronic brakes with manual park
- R-net LCD controls as standard

## Safety.

- Full compliance AS 3695 250kg
- For Occupant Vehicle Transport 159kg  
(Crash tested to AS 3696:19 159kg)

## Arrange a test drive.

Chat to a Glide Wheelchair specialist and experience the extraordinary features of the CentroGlide Bariatric in person. Find your nearest distributor at [glide.com.au/nearest](http://glide.com.au/nearest)

## Summary Specifications.



### Recommended User Weight

Up to 250kg AS 3695 full compliance



### Overall Width

620mm (wheel base)



### Back Height

1030mm Top of push handle (custom height available)



### Brakes

Electronic and manual park



### Overall Length

900mm (foot plate removed)



### Seat Depths

40 to 58cm



### Arm Rests

Flip up (adjustable)



### Battery

Up to 40km range



### Seat Widths

49 to 73cm



### Seat Height

Std 48cm top of rail 42cm (lower seat option)



### Leg Rests

Swing away with aluminium footplates detachable std. Centre mount and power elevating also available



### Wheelchair Weight

Starting at 85kg (without batteries and upholstery)

# medifab

[www.medifab.com](http://www.medifab.com)

+64 3 307 9790 | [solutions@medifab.com](mailto:solutions@medifab.com)

22 Detroit Drive, Rolleston 7675