

## The versatile power wheelchair.

CentroGlide is our leading model. It's a robust mid-wheel-drive power wheelchair, designed and built by us in Australia to provide remarkable stability and trouble-free service for years.

### Completely customisable.

CentroGlide is a fully-scripted wheelchair, built to the occupant's precise requirements.

### Independence and freedom of movement.

CentroGlide has been designed to allow remarkable manoeuvrability in tight spaces, and stability over kerbs. It comes with 50° seat tilt as standard, with options for vertical lift and anterior reach to allow the users to reach high shelves or transfer from the chair in a standing position.

### Power through the day.

CentroGlide is built with the functionality and reliability to give users total mobility, indoors or out.

## Designed for comfort and ease of use.

We have put decades of user feedback and experience to develop all the features that make CentroGlide a popular choice.



### Customised seating.

For ultimate comfort many seating options are available. Choose customised cushion, backrest and headrest.



### Electronics.

Fully programmable R-net joystick control system. Compatible with many other 3rd party specialty electronics.



### Active suspension.

Our patented 6-wheel APT suspension improves traction and aids ascending and descending kerbs safely without affecting the stability and angle of the seat.



### Power seating.

As standard, the seat has a 50° tilt. Options include power seat lift and anterior reach.



### Wheels.

CentroGlide has 14" drive wheels, 8" and 7" castor wheels as standard. Other sizes are available.



### Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



### Battery range.

The powerful battery gives you a range of up to 40km on a single charge.

# CentroGlide sets the benchmark.

## Features.

- Power base with removable seating module
- High quality German 2 pole with optional 4 pole high torque motors
- Vehicle tie down points
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) kerb climbing ability
- 80mm ground clearance
- Max ramp incline/decline 10° (as per AS3695 testing)
- 4.1/3.50-6 (14") drive wheels
- Electronic brakes with manual park
- R-net LCD Std controls

## Safety.

- Full compliance AS 3695 175kg
- Crash tested to AS 3696:19 to max user weight 175kg

## Arrange a test drive.

Chat to a Glide Wheelchair specialist and experience the extraordinary features of the GlideCentro in person. Find your nearest distributor at [glide.com.au/nearest](http://glide.com.au/nearest)

## Specifications.



**Recommended User Weight**  
Up to 175kg AS 3695 full compliance



**Overall Length**  
900mm (foot plate removed)



**Seat Widths**  
31 to 61cm



**Overall Width**  
620mm (wheel base)



**Seat Depths**  
31 to 58cm



**Seat Height**  
46cm (top of rail using 0-46 set up)



**Back Height**  
1030mm Top of push handle (custom height available)



**Arm Rests**  
Flip up (adjustable)



**Battery**  
Up to 40km range



**Leg Rests**  
Swing away detachable



**Frame Colour**  
Wide choice



**Brakes**  
Electronic and manual park