

CentroGlide XT

GLIDE

The power to go off-pavement.

CentroGlide XT is a special version of our leading power wheelchair, with features to give users the freedom to navigate grass playing fields, dirt tracks or uneven outdoor surfaces with confidence.

XT stands for 'Extreme Terrain'.

Our flagship CentroGlide is already an industry-leading mid-wheel-drive power wheelchair with remarkable stability and the ability to traverse kerbs and median strips with ease. However, the CentroGlide XT adds even more capability with low pressure wide knobby drive tyres and 9" front castors.

How can something so tough be so comfortable?

All our wheelchairs are scripted and made-to-measure right here in Australia. Your Glide chair is literally built for your comfort, with almost endless customisation to ensure all-day comfort, freedom of movement, and independence.

Built to keep up with your lifestyle.

CentroGlide XT is ready for adventure, and built with the functionality and reliability to give users total mobility, indoors or out.

Designed for comfort and ease of use.



Customised seating.

For ultimate comfort many seating options are available. Choose customised cushion, backrest and headrest.



Active suspension.

Our patented 6-wheel APT suspension improves traction and aids ascending and descending kerbs safely without affecting the stability and angle of the seat.



Wheels.

CentroGlide XT has 14" wide drive wheels with special low-pressure off-road tread tyres and 9" castor wheels.



Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



Electronics

Fully programmable R-net joystick control system. Compatible with many other 3rd party speciality electronics.



Power seating.

As standard, the seat has a 50° tilt. Options include power seat lift and anterior reach.



Battery range.

The powerful battery gives you a range of up to 40km on a single charge.



CentroGlide XT sets the benchmark.

Features.

- Power base with removable seating module
- High quality German 4 pole high torque motors
- Vehicle tie down points
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) kerb climbing ability
- 80mm ground clearance
- Max ramp incline/decline 10°
(as per AS3695 testing)
- 14" drive wheels with off-road low pressure tyres
- Electronic brakes with manual park
- R-net LCD Std controls (dynamic optional)

Safety.

- Full compliance AS 3695 175kg
- Crash tested to AS 3696:19
to max user weight 175kg

Arrange a test drive.

Chat to a Glide wheelchair specialist and experience the extraordinary features of the GlideCentro XT in person. Find your nearest distributor at glide.com.au/nearest

Specifications



**Recommended
User Weight**
Up to 175kg AS 3695
full compliance



Overall Length
980mm (foot plate
removed)



Seat Widths
31 to 61cm



Overall Width
740mm (wheel base)



Seat Depths
31 to 58cm



Seat Height
46cm (top of rail using
0-46 set up)



Back Height
1030mm Top of push
handle (custom height
available)



Arm Rests
Flip up (adjustable)



Battery
Up to 40km range



Leg Rests
Centre mount



Frame Colour
Wide choice



Brakes
Electronic and
manual park