

Grillo Gait Trainer




**VERSATILE POSTERIOR
OR ANTERIOR WALKING AID**

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Key

BLUE	Suitable for Anterior and Posterior
GREEN	Suitable for Anterior only
YELLOW	Suitable for Posterior only

 Therapists Normally Choose


Size Options

Measurements	Size Mini (0)	Size Small (1-2)	Size Medium (3-4)	Size Large (5)
Anterior Code	5750-0011-000 <input type="radio"/>	5750-0051-12B <input type="radio"/>	5750-0053-34B <input type="radio"/>	5750-0055-05B <input type="radio"/>
Posterior Code	5750-0140-000 <input type="radio"/>	5750-0051-12B + 5750-1000-000 <input type="radio"/>	5750-0053-34B + 5750-1000-000 <input type="radio"/>	5750-0055-05B + 5750-1000-000 <input type="radio"/>
Pelvic Support Height	42-61cm	49-84cm	54-100cm	+
Trunk Support Height	51-75cm	69-100cm	76-125cm	5750-1000-000
Support Circumference	53-79cm	60-86cm	74-120cm	84-140cm
Frame Width	57cm	65cm	76cm	80cm
Frame Length	73cm	75cm	91cm	111cm
User Height Range	75-100cm	90-135cm	110-165cm	150-190cm
Max Load	35kg	45kg	90kg	110kg
Total Weight of Anterior figuration	13.5kg	16.2kg	19.5kg	23kg
Total Weight of Posterior figuration	13.5kg	15.7kg	17.7kg	22kg

ANTERIOR PACKAGE	POSTERIOR PACKAGE
-------------------------	--------------------------

Anterior Package Contains:

- Pelvic support
- Trunk support (not inc. in Mini)
- D Handle (not inc. in Mini)
- Transparent tray (Mini only)
- Grab rail (Mini only)
- Pelvic harness
- Front direction locks
- Drag brakes
- Anti-reverse locks
- Foot activated brakes



Mini
Sml/Med/Lrg

Posterior Package Contains:

- Trunk support
- Pelvic support
- Pelvic harness
- Adjustable handle bars
- Direction locks
- Drag brakes
- Anti-reverse locks
- Foot activated brakes



Sml/Med/Lrg
Mini

Head Support

Multiple Adjustable Head Support

For Posterior Only

5750-4600-000



Arm Support Options

Arm Gutter Supports

Mini Posterior
5750-5500-000

Small Anterior
5750-5501-1-2

Small Posterior
5750-5501-1-2 + 5750-5550-000

Medium/Large Anterior
5750-5503-3-4

Medium/Large Posterior
5750-5503-3-4 + 5750-5550-000



Arm Straps

Mini
5750-6860-000

Small
5750-6862-000

Medium/Large
5750-6864-000



Arm Retainers

5640-6865-000



D-Handles

Included in Anterior package

Small
5750-5701-1-2

Medium
5750-5703-3-4

Large
5750-5705-005



Multi-Adjustable Grip Handles

For Mini only

5750-5700-000



Ergonomic Handles

Posterior
5750-5735-000

Anterior
5750-5730-000



Hand Grips

Small
5750-5721-1-2

Medium/Large
5750-5723-3-4



Key **BLUE** Suitable Anterior and Posterior **GREEN** Suitable for Anterior **YELLOW** Suitable for Posterior

Trunk Support

Included in the standard Anterior and Posterior package

Mini <input type="radio"/>	Small <input type="radio"/>	Medium <input type="radio"/>	Large <input type="radio"/>
5750-2850-000	5750-2851-1-2	5750-2853-3-4	5750-2855-005



Lumbar Support

Adjustable lower back support

Small <input type="radio"/>	5750-2501-1-2
Medium <input type="radio"/>	5750-2503-3-4



Black Belt

For Size 0

Requires tray for use

5750-6820-000



Seat Options



Padded Saddle Seat

Fits Mini/Small/Medium <input type="radio"/>	5750-0211-1-2
Fits Medium <input type="radio"/>	5750-0213-3-4
Fits Large <input type="radio"/>	5750-0215-005



Height Adjustable Folding Seat

Small <input type="radio"/>	5750-0201-1-2
Medium <input type="radio"/>	5750-0203-3-4



Ergonomic Saddle Seat

Mini <input type="radio"/>	5750-0250-000
Small <input type="radio"/>	5750-0251-1-2
Medium <input type="radio"/>	5750-0252-3-4



Key **BLUE** Suitable Anterior and Posterior **GREEN** Suitable for Anterior **YELLOW** Suitable for Posterior


Leg Dividers

Padded Leg Divider

For Mini only

A 53cm B 20cm C Ø5cm

5750-5400-000



Leg Divider with Moving Leg Cuffs

Small



A 60cm B 20cm C Ø5cm

5750-5401-1-2

Medium



A 64cm B 28cm C Ø5cm

5750-5403-3-4

Large



A 79cm B 28cm C Ø5cm

5750-5405-005



Additional Accessories

Weight Bars

Small

5750-8431-1-2



Medium

5750-8433-3-4



Assistant Push Handle

Mini/Small/Medium

5750-8600-000



medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Grillo_20200730