

Hero Car Seat

The Hero child car seat is a specifically modified car seat to accommodate the postural requirements of children with additional needs. It is suitable for installation in most commercially available vehicles that have an ISOFIX connection. The car seat is designed for use on the back seat of a car, and with the ISOFIX anchorage system providing a secure attachment to the vehicle, it does not have to be fastened with the car seat belts, guaranteeing easy handling.

Swivelling car seat for forward or rear-facing

Swivel by 90° so that the child can be seated in the car seat without straining the parent or caregiver's back

Large seat recline

Height adjustable back support and head supports

5-point-harness with upholstery, height adjustable

Belt mechanics installed into the front of the seat

360° swivel base with adjustable footwell stilt

ACCESSORIES

Security Tray
4181-5530-000



FUNCTIONAL BENEFITS

- Switch the Hero around for either forward or rear facing. The Hero can be used both forward and rear-facing direction.
- To support children with limited upper body control, the Hero has a tilt function to adjust the seat angle.
- The Hero is equipped with its own 5-point chest and shoulder strap that securely holds your child in place.

Package includes

- Available in Black/Red
- Rear/Forward facing position
- Back Support Height Adjustable
- 10° tilt with swivel base 360°
- Five-point harness
- Lateral Supports as standard

Practical and functional



In direction of travel; adjustable back height



Rotatable by 90°, 180° or 360° for ease of transfers



In rearward position



Lateral Supports fitted as standard

| Code | Seat Width (knee) | Seat Width (pelvis) | Seat Depth | Back Height | Crown Height | Distance Laterals | Total Height | Total Width | Total Depth | Capacity |
|---------------|-------------------|---------------------|------------|-------------|--------------|-------------------|--------------|-------------|-------------|----------|
| 4181-0010-020 | 280 mm | 210 mm | 250 mm | 210 mm | 410 mm | 120 mm | 660 mm | 430 mm | 630 mm | 18 kg |
| | | | | 340 mm | 530 mm | 210 mm | 730 mm | | | |

*Measurements with standard cushion