Ipai LGT Car Seat

The Ipai LGT lives up to its name, characterised by a light weight easy-to-use design and provides mild body support and care for children with additional needs. Comes standard equipped with ISOFIX, thoracic supports and a 5-point harness that is padded at the shoulders. As part of the supportive harness a large, stable tee pad guarantees firm support across the child's pelvis.

Light weight, comfortable and safe!



Package includes

- Available in Black
- Five-point harness
- Backrest height adjustable



Backrest Height



ISO Fix



Side Impact Protection System

ACCESSORIES



FUNCTIONAL BENEFITS

- Simple height adjustment function allows the backrest 20cm of height range providing greater comfort and safety for the child.
- Easy release ISOFIX attachment makes lpai LGT a breeze to install in any ISOFIX compatible vehicle. Please note, the vehicle seat belt must be used in conjunction with the car seat.
- Have the peace of mind that your child is kept safe with the Side Impact Protection System that can be unscrewed to secure the car seat against the side of the vehicle door, improving safety in the event of side impact to the vehicle.



Code	Seat Width (rear)	Seat Width (front)	Seat Depth	Back Height	Head Support Height	Lateral Support Height	Distance between Laterals	Back Support Recline Angle	Total Height	Overall Width	Total Depth	Total Weight	Capacity
4081-0011-060	31 cm	34 cm	32 cm	38-58 cm	59-79 cm	20-33 cm	14-28 cm	85°-162°	67-85 cm	48 cm	39 cm	7 kg	36 kg