



THE SUPPORT YOU NEED TO BE ACTIVE

The J3 Carbon back provides great stability and comfort. Through the right pelvis support it offers firm stabilisation of the upper trunk and encourages an optimal spinal curve. This enhances the stability you need for a very active live, reduces back pain and provides a great comfort. Plus, the force that's generated when propelling your wheelchair is transferred more effectively and efficiently through a solid back into the frame. That means you're generating greater movement of your wheelchair with less effort.

IT'S SO LIGHT

Carbon is lightweight and strong – therefore the ideal material for this backrest. Offering the same strength as aluminium (but 60 % lighter), the J3 Carbon back reduces the weight that a wheelchair user has to push, day-by-day. A low thoratic back weighs just 690 grams! You will feel the difference!



CHOOSE YOUR BACK

Different needs demand adequate solutions. The four back heights (from 17 to 38 cm) focus on the active user with a need for low to moderate posterior trunk support. Mild lateral trunk support is also given by the 5 cm contour depth.



→ Active contour

Mild lateral trunk support
through 5 cm contour depth




→ Different support levels
Always the right support
(four heights from 17 to 38 cm)



WHERE ARE MY KEYS?

You now have somewhere to store those easy to misplace items with the new accessory bag. Not only will it keep all of your valuables in one place, but it also smartly covers the gap between your back and cushion. You'll wonder what you ever did without it.







Max. user weight:	125 kg
Backrest width:	32 / 36 / 41 / 46 cm
Backrest height:	17 / 24 / 30 / 38 cm
Contour depth:	5 cm
Back angle:	± 20°
Backrest weight:	From 690 g
Hardware	Fixed and quick-release version
Accessories	Accessorie bag and lumbar supports

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



Sunrise Medical Australia Pty Ltd 6 Healey Circuit, Huntingwood, NSW 2148 Ph: (02) 9678 6600

Fax: (02) 9678 6655

Email: enquiries@sunrisemedical.com.au

www.sunrisemedical.com.au