

Corner Seat & Nursery Table




CORNER SEAT & NURSERY TABLE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Measurements	User Measurements
A: Hip Width	
B: Seat Depth	
C: Seat to top of shoulder	
D: Seat to top of head	
E: Chest Width	
F: Seat to Axilla	
G: Leg Length	



Product Specifications	Size 1	Size 2	Size 3
Average Age	6 Months - 2 Years	1 - 6 Years	4 - 10 Years
Base Size	50cm x 50cm	54cm - 56.5cm	60cm - 65.5cm
Weight Capacity	25 kg	35 kg	55 kg
Seat Width	33cm	39cm	49.5cm
Seat Depth	35cm	37cm	46.5cm
Seat to Top Of Backpad Height	31.5cm - 42.5cm	44cm - 62cm	50cm - 68cm

Corner seat

Bright, comfortable and ideal first seat for any child, which also suits infants and young children.

	Size 1 (Select colour below) <input type="radio"/>	Size 2 (Select colour below) <input type="radio"/>	Size 3 (Select colour below) <input type="radio"/>
Age range	6 mths - 2 yrs	1 - 6 yrs	4 - 10 yrs
Max user weight	25kg	35kg	55kg
Back pad height	31.5 - 42.5cm	44 - 62cm	50 - 68cm
Seat depth	35cm	37cm	46cm
Base size	40 x 50cm	52 x 54cm	60 x 63cm
Code	<input type="radio"/> Lime: 3140-0041-054	<input type="radio"/> Lime: 3140-0042-054	<input type="radio"/> Lime: 3140-0044-054
	<input type="radio"/> Blue: 3140-0041-066	<input type="radio"/> Blue: 3140-0042-066	<input type="radio"/> Blue: 3140-0044-066
			

Triangular Cushions

For when a flat back with wings is preferred. Fit all seat sizes. Can also be used beside legs to prevent external rotation.

	Lime	Blue
Code	3140-2450-054	3140-2450-066
Price		



Nursery Table

Effective when used together with the Corner Seat.

	Size 0	Size 1	Size 2
Code	3140-7101-000	3140-7102-000	3140-7103-000
Base size	60.5 x 38cm	60.5 x 38cm	70 x 42 cm
Table top size	51 x 60.5cm	51 x 60.5cm	51 x 70.5cm
Table height	19-28cm	24.5-34.5cm	29-44 cm
Tray top cut-out	22 w x 16 cm d	22 w x 16 cm d	32 w x 16 cm d
Table top angle range	+40° to -40° from horizontal		
Fits corner seat size	Size 1	Size 2	Size 3
Price			



medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: CornerSeatNurseryTable_20210907