

# Jenx Ly-on

Modern, child-centred styling, combined with a range of accessories make Ly-on the ideal platform to build a structured yet safe and comfortable posture for children with complex postural or medical needs.

## Feeding, Education and Access Station for users up to 70kg

**Activity Frame and Tablet Holder** easily clip into place allowing the user to be entertained during long periods of lying.

**Standard with a wide range of height adjustment** – supporting children to engage in low to the floor activities whilst also being able to be correctly positioned for even the tallest of carers

**15° tilt recline and forward** allows for weight distribution to be altered to aid in pressure management.




**Side Lying Straps** offering controlled positioning at the chest or pelvis.

**Fold down sides**, allowing for easy engagement with the child and also making transfers in and out of Ly-on safe and dignified.

**Removable and washable** at 60°, this accessory sheet allows for precise positioning of the various support options.

### STANDARD PACKAGE INCLUDES

- 4 braking castors - one of which is directional lockable to aid in manoeuvrability
- Hydraulic lift action
- Tilt-in-space mechanism
- Fold-down side rails
- 4 Side Walls in green vinyl - 250mm high
- End bumpers in green vinyl
- Memory foam topped lying surface (covered with grey wipe clean PU material)
- 1 Hook Receptive Sheet
- 10 Hook Tabs to hold Hook Receptive Sheet in place
- Push handle



**Activity Frame/ Tablet Holder**



**Height Adjustable Platform**



**Support Blocks and Straps**



**Hooks and Loop Receptive Sheet**




**Fold Down Sides**



Code	Max. User Weight	Max. User Height	Tilt	Support Platform Dimensions	Side Wall Height (3 sizes available)	Foot Support Dimensions	Height From Floor	Product Weight
7520-0015-054	70kg	160cm	15°- 15°	168 x 56cm	20/25/30cm	196x 67cm	52-101cm	76kg

# Accessories


**Extra Hook Receptive Sheet**  
7520-7000-000




**Extra Hook Tab**  
7520-7050-000



**Side Wall:**  
Size 1: 7520-7211-054  
Size 2: 7520-7212-054  
Size 3: 7520-7213-054



**Side Lying Strap:**  
Size 1: 7520-1871-000  
Size 2: 7520-1872-000  
Size 3: 7520-1873-000



**Rotation Control Strap:**  
Size 1: 7520-1881-000  
Size 2: 7520-1882-000  
Size 3: 7520-1883-000




**Rotation Control Frame**  
7520-1870-000



**Accessory Mount**  
7520-7925-000



**Activity Frame:**  
Size 1: 7520-7911-000  
Size 2: 7520-7912-000



**Tie-on Beam**  
7520-7920-000



**Tablet Holder**  
7520-7930-000




**Triangular Shelf Block**  
7520-7750-000




**Triangular Positioning Block**  
7520-1850-000



**Bridge Block:**  
Size 1: 7520-7441-000  
Size 2: 7520-7442-000  
Size 3: 7520-7443-000



**Half Roll Block:**  
Size 1: 7520-7821-000  
Size 2: 7520-7822-000



**Roll Pillow**  
7520-7800-000



**Flat Oval Pillow**  
7520-7500-000



**Accessory Cart**  
7520-7950-000



**Accessory Box**  
7520-7940-000

