

Accessories (to suit Junior & Multiseat)

Pommel

Size 0 #3615-1650-000
Size 1 #3615-1651-000
Size 2 #3615-1652-000




Leg Control Pads


Size 1 #3610-1861-000
Size 4 #3614-1864-000
Sz 4 Brackets #3610-1854
Size 5 #3614-1865-000



Back Post Extension
75mm #3610-2480-000




Foam Back Pad
#3610-2390-000




Adjustable Ischial bar
For Seat Well #3615-0501



Hip Pad Extensions
#3610-1201-000



Box Style Footplate
#3430-0281-000



Footplate Pad
#3615-5301-000



Sandals

Size 1: #5356-6081-000
Size 2: #5356-6082-000
Size 3: #5256-6083-000



Sandal Raisers

Size 1: #5356-6091-000
Size 2: #5356-6092-000
Size 2: #5356-6093-000



Jenx Multiseat

The larger version of the Jenx Junior, Multiseat is a practical, comfortable, contoured and flexible modular seating system, designed to suit the needs of older children with mild to complex needs.

Versatile modular seating

All support pads and accessories are designed to **QUICKLY & EASILY FLIP AWAY** without being removed

TWO PIECE BACK with independent lumbar support gives support where it's required

Flip-up **ANGLE ADJUSTABLE FOOTPLATE** to accommodate dorsi or plantar flexion

PELVIC CONTROL SYSTEM with independently adjustable contoured seating

Z BASE has 33.5-57 cm Seat-To-Floor height and a footprint of 63 x 83 cm



MULTISEAT



Standard Configuration:

- Base and Seat Unit with Hip Belt and Hip Supports
- Foam Seat Pad and Cover, Lumbar Support, Back Pad and Cover
- Armrest Bracket, Armrests, Tray, Footpost, Footplate, Sandals, Size 2 Laterals with strap and 2-link Multigrip headrest



Highly adjustable and configurable to suit.

Ease of transfer with tilt-in-space and swing-away supports.

Open design for comfort and support.

SIZE	CODE	AGE RANGE	MAX USER WEIGHT	TILT-IN-SPACE	SEAT DEPTH	HIP WIDTH	BACKREST HEIGHT	SHOULDER HEIGHT	LOWER LEG LENGTH
4	3614-0014-011	8-14	100 kg	10-30°	28.5-45 cm	20-44 cm	46.5-59 cm	46.5-59 cm	25-47 cm
5	3614-0025-011	10-18+	100 kg		35.5-52 cm		52.5-66 cm	52.5-66 cm	25-47 cm