

Kool Stride Sports Buggy

This popular stroller for larger children with disabilities now has a unique advanced feature: an extendable rear wheel system. As the user grows, the rear wheels can extend backwards to maintain stability and control.

Lightweight all-terrain buggy for 10 year old to young adults



Dual parking brakes

Lightweight aluminium frame

5-point safety anterior chest support

Weight limit of 60 kg

Carrier basket under seat, and pocket on the back of the seat which can be used to hold drink bottles etc.

5°, cambered large wheels make walking easy on pathways, gravel or grass



Travelling made easy

Quick release wheels for easy folding down, simple storage and transport.



Roll the sunhood side flaps up to let the cool air in



Stow the sunhood back



Folds flat and compact for transportation

EXPLORE...

Serious sports buggy for cruising outdoors!

- Jogging
- Bushwalks
- Beaches
- Farms
- and much more...



ACCESSORIES

Rain Cover
2650-7600-010

UV Cover
2650-7700-260



Code	Seat Depth	Seat Width	Shoulder Height	Seat Back Height	Foot Support Height	Overall Width	Max Load	Rear Wheel Diameter	Approx Stroller Weight
2650-0001-000	30 cm	43 cm	50 cm	60 cm	45 cm	68 cm	60 kg	50 cm	15 kg