


# Octopus Positioning System

The Octopus offers comfortable head and anterior upper trunk stability, which mounts to wheelchairs and seat shells. Unique retraction retainers fold backwards and click into place with hinges.

## EFFECTIVELY POSITION THE HEAD AND UPPER BODY



Available with a range of mounting hardware - Your Medifab sales consultant can assist as to the best style to suit your client's needs.



Both arms retract for unrestricted transfers into the chair.

Lower each of the positioning arms and adjust for optimal support.

Fasten clasp to hold both arms together in position.

Size	19 mm (3/4") Ball	Neck Size	A	B	C	D
0	1290-4110-000	Up to 85 mm	200 mm	95 mm	200 mm	35 mm
1	1290-4111-000	Up to 100 mm	280 mm	120 mm	230 mm	50 mm
2	1290-4112-000	Up to 125 mm	340 mm	140 mm	260 mm	50 mm
3	1290-4113-000	Up to 150 mm	400 mm	165 mm	300 mm	50 mm

## MANTA HEAD SUPPORT



- Attachment for the Octopus to support the head position
- Universal size to suit various head circumferences
- Easy mounting to the Octopus
- Width adjustable wings



### Full front support Manta attachment

1290-4810-000



### Side support Manta attachment

1290-4820-000

