

Safe Surround Plus Bed

The Safe Surround Plus beds are for children and adults that are able to pull themselves to kneeling and/or standing positions and require higher bed doors to provide a safe environment.

Safe environment for active people with challenging behaviour

Designed for those who have **restless motor functions** or **seriously disruptive behaviours**

High quality craftsmanship: durable Beech timber with attractive aesthetics



Extremely stable and strong construction to withstand forces

Lower the mattress level to enable **self-transfers**

Standard Features



Roomy internal measurements and door height options to protect almost any active child from harming themselves.




- Height adjustable lying surface
- 5 section profiling sleeping platform



- Standard foam mattress
- External door locking bar

DOOR SYSTEM

Two handed exterior door release system with special locking system that prevents the doors being opened from the inside.




PROFILING SLEEP PLATFORM

Assists physical needs of the user, and height adjustable to achieve ergonomic height for the parent/carer



OPTIONS


Mattress Incontinence Cover
7956-0027-000



Pressure Relief Foam Mattress
7956-0011-051



Internal Padding as Required. Velcro attachment.



Code	Door Height	Weight Limit	Inside Dimensions	Overall Height	Lying Height	Outside Dimensions
7953-0008-000	135 cm	150 kg	200 x 100 cm	184 cm	55 cm - 95 cm platform height (top of mattress)	215 x 117cm (with locking bars in place)
7953-0011-000	170 cm			220 cm		