Vigour Back Support

Lightweight yet tough to meet the daily lifestyle demands of active wheelchair users. Back support angle and seat depth offers exceptional adjustability and can be facilitated with the user in the wheelchair. Ventilated spacer-fabric and shell system increases air-flow through the back offering daily comfort.


VIGOUR LO

New invertible shell design means the Vigour Lo can be used either way up on the wheelchair. Sacral and lumbar height support allowing unrestricted upper body range of movement.

Designed for the individual who requires:

- A solid back support with the ability to self-propel
- Minimal postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- For clients that have mild torso asymmetries

Approx. Vigour Lo weight is from 0.8 to 1.5 kg




VIGOUR MID

Lightweight and designed for efficient propelling, with scapular subtractions for unimpeded arm

Designed for the individual who requires:

- A solid back support with the ability to self-propel
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability
- For clients that have mild torso asymmetries

Approx. Vigour Mid weight is from 0.9 to 2.5 kg



VIGOUR HI

Solid full height back support available in three heights for maximum trunk and shoulder support.

Designed for the individual who requires:

- A solid back support
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability

Approx. Vigour Hi weight is from 1 to 3.3 kg




MICRO-CUBE FOAM UPGRADE

Air-flow breathability for maximum control of temperature and moisture. Micro-cube cells provide zoned weight distribution enhancing comfort.



HARDWARE OPTIONS




QUAD MOUNT (DEFAULT)

Default hardware with 4 points of attachment for ultimate strength. #1259-2400-100



DUAL MOUNT (OPTION)

Increased adjustability, angle and growth. #1259-2200-100



SINGLE MOUNT (OPTION)

For lower height backs, and as a lightweight option. #1259-2100-100



Dual-layer cover reduces SHEAR FORCES with air-mesh side panels to counter moisture and temperature

Add laterals as required from our LARGE MODULAR RANGE

Can be upgraded with DIFFERENT FOAM PADDING OPTIONS for excellent results

SLOTS & HOLE POSITIONS

on shell allow easy installation of head supports, lateral supports and harnesses

W/C Width	Vigour Lo Back Support Height	Vigour Lo Back Support Code	Vigour Lo Spare Cover Code	
10"	6"	1252-1006-00L	1252-1006-00L-C	
12"	7"	1252-1207-00L	1252-1207-00L-C	
14"	8"	1252-1408-00L	1252-1408-00L-C	
16"	10"	1252-1610-00L	1252-1610-00L-C	
17"	10"	1252-1710-00L	1252-1710-00L-C	
18"	11"	1252-1811-00L	1252-1811-00L-C	
20"	11"	1252-2011-00L	1252-2011-00L-C	

Vigour Mid Back Support Height	Vigour Mid Back Support Code	Vigour Mid Spare Cover Code
10"	1252-1010-00M	1252-1010-00M-C
12"	1252-1212-00M	1252-1212-00M-C
14"	1252-1414-00M	1252-1414-00M-C
16"	1252-1616-00M	1252-1616-00M-C
17"	1252-1717-00M	1252-1717-00M-C
17"	1252-1817-00M	1252-1817-00M-C
17"	1252-2017-00M	1252-2017-00M-C

W/C Width	Vigour Hi Back Support Height	Vigour Hi Back Support Code	Vigour Hi Spare Cover Code	Vigour Hi Back Support Height	Vigour Hi Back Support Code	Vigour Hi Spare Cover Code
10"	12"	1252-1012-00H	1252-1012-00H-C	14"	1252-1014-00H	1252-1014-00H-C
"	16"	1252-1016-00H	1252-1016-00H-C			
12"	14"	1252-1214-00H	1252-1214-00H-C	16"	1252-1214-00H	1252-1214-00H-C
	18"	1252-1218-00H	1252-1218-00H-C			
14"	16"	1252-1416-00H	1252-1416-00H-C	18"	1252-1418-00H	1252-1418-00H-C
"	20"	1252-1420-00H	1252-1420-00H-C			
16"	18"	1252-1618-00H	1252-1618-00H-C	20"	1252-1620-00H	1252-1620-00H-C
"	23"	1252-1623-00H	1252-1623-00H-C			
17"	18"	1252-1718-00H	1252-1718-00H-C	20"	1252-1718-00H	1252-1718-00H-C
	23"	1252-1723-00H	1252-1723-00H-C			
18"	18"	1252-1818-00H	1252-1818-00H-C	20"	1252-1820-00H	1252-1820-00H-C
п	23"	1252-1823-00H	1252-1823-00H-C			
20"	18"	1252-2018-00H	1252-2018-00H-C	20"	1252-2018-00H	1252-2018-00H-C
II .	23"	1252-2023-00H	1252-2023-00H-C			