

# Vigour Cushion



## CLINICAL


Designed for the user who:

- Requires pressure reduction and redistribution
- Requires pelvic stability
- Needs moderate control of posture and upper extremities
- Needs moderate control of abduction or adduction propensities

## FUNCTIONAL

- Gel-foam ischial well evenly distributes pressure under the pelvis
- Special foam layers' combination, strategically combined to offer a supportive foundation for a correct seating posture
- Slight thigh elevation built-in for increased balance
- Specific contouring offers a high degree of pelvic stability
- Coccyx channel for enhanced pressure relief (for sizes 16" +)
- Dual-layer cover reduces shear forces, and breathable air-mesh on cushion sides counters moisture
- Inner incontinence cover included
- Customisable for many requirements such as leg length discrepancies

## CUSHION SIZE CHART



- Red squares are standard sizes.
- Non-standard or custom sizes are available on request.

## CUSHION DETAILS

| Width   | Height | Capacity       | Trough Height (mm) |      | Product Weight |              |
|---------|--------|----------------|--------------------|------|----------------|--------------|
|         |        |                | Std                | High | Std            | High         |
| 10"     | 2"     | 200lbs (90kg)  |                    |      |                |              |
| 11"-15" | 2.5"   |                | 30                 | 50   | (16x16)        | (16x16)      |
| 16"-17" | 2.75"  | 350lbs (160kg) |                    |      | <b>1.4kg</b>   | <b>1.5kg</b> |
| 18"-20" | 3"     |                |                    |      |                |              |

## PRODUCT PERFORMANCE


|                   |   |  |   |  |  |
|-------------------|---|--|---|--|--|
| Pressure Care     |   |  | 3 |  |  |
| Positioning       |   |  | 3 |  |  |
| Lateral Stability |   |  | 3 |  |  |
| Offloading        |   |  | 3 |  |  |
| Changing Needs    | 1 |  |   |  |  |

- Note: This is a guide only and ratings are dependent on correct use of product.

# Codes for Standard Sized Cushions

## Our cushion codes are meaningful!

Here's how to generate a product code to place your order.




### EXAMPLE

A standard Contour Spex cushion of width 18" and length 20", would have the code 1105-1820-OSC.

## SPEX & VIGOUR CUSHION MEASUREMENTS GUIDE

This FRONT-VIEW of a Spex/Vigour cushion, shows the different measurements referred to on the subsequent product pages.



## Stand out from the crowd!

Upgrade your cushion with a splash of colour, for no extra cost!

We offer 5 colour options for our Spex cushion outer covers alternative to the standard black fabric. Simply add the upgrade codes to your cushion order.



**Granite Grey**  
1111-0491-090

**Ocean Blue**  
1111-0491-060

**Royal Purple**  
1111-0491-077

**Calypso Pink**  
1111-0491-028



**Chilli Red**  
1111-0491-022

**Mango Orange**  
1111-0491-030

**Bumblebee Yellow**  
1111-0491-044

**Leaf Green**  
1111-0491-054

## Spare cushion covers

All our cushions come standard with an inner incontinence cover and an outer cover, however if you would like a spare incontinence or outer cover, these are available.

See our scripting and ordering guide for ordering codes or get in touch with your local Spex distributor.

