

Vigour High Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	
8" 20 cm	8" 20 cm	1103-0808-OHC	<input type="checkbox"/>
	10" 25 cm	1103-0810-OHC	<input type="checkbox"/>
10" 25 cm	10" 25 cm	1103-1010-OHC	<input type="checkbox"/>
	12" 30 cm	1103-1012-OHC	<input type="checkbox"/>
	XX XX	1103-10XX-OHC	<input type="checkbox"/>
11" 28 cm	10" 25 cm	1103-1110-OHC	<input type="checkbox"/>
	12" 30 cm	1103-1112-OHC	<input type="checkbox"/>
	XX XX	1103-11XX-OHC	<input type="checkbox"/>
12" 30 cm	12" 30 cm	1103-1212-OHC	<input type="checkbox"/>
	14" 35 cm	1103-1214-OHC	<input type="checkbox"/>
	XX XX	1103-12XX-OHC	<input type="checkbox"/>
13" 33 cm	12" 30 cm	1103-1312-OHC	<input type="checkbox"/>
	14" 35 cm	1103-1314-OHC	<input type="checkbox"/>
	XX XX	1103-13XX-OHC	<input type="checkbox"/>
14" 35 cm	14" 35 cm	1103-1414-OHC	<input type="checkbox"/>
	16" 40 cm	1103-1416-OHC	<input type="checkbox"/>
	18" 45 cm	1103-1418-OHC	<input type="checkbox"/>
	XX XX	1103-14XX-OHC	<input type="checkbox"/>
15" 38 cm	14" 35 cm	1103-1514-OHC	<input type="checkbox"/>
	16" 40 cm	1103-1516-OHC	<input type="checkbox"/>
	XX XX	1103-15XX-OHC	<input type="checkbox"/>
16" 40 cm	16" 40 cm	1103-1616-OHC	<input type="checkbox"/>
	18" 45 cm	1103-1618-OHC	<input type="checkbox"/>
	20" 50 cm	1103-1620-OHC	<input type="checkbox"/>
17" 43 cm	XX XX	1103-16XX-OHC	<input type="checkbox"/>
	17" 43 cm	1103-1717-OHC	<input type="checkbox"/>
	18" 45 cm	1103-1718-OHC	<input type="checkbox"/>
18" 45 cm	XX XX	1103-17XX-OHC	<input type="checkbox"/>
	18" 45 cm	1103-1818-OHC	<input type="checkbox"/>
	20" 50 cm	1103-1820-OHC	<input type="checkbox"/>
	22" 55 cm	1103-1822-OHC	<input type="checkbox"/>
19" 48 cm	XX XX	1103-18XX-OHC	<input type="checkbox"/>
	18" 45 cm	1103-1918-OHC	<input type="checkbox"/>
	20" 50 cm	1103-1920-OHC	<input type="checkbox"/>
20" 50 cm	XX XX	1103-19XX-OHC	<input type="checkbox"/>
	18" 45 cm	1103-2018-OHC	<input type="checkbox"/>
	20" 50 cm	1103-2020-OHC	<input type="checkbox"/>
	22" 55 cm	1103-2022-OHC	<input type="checkbox"/>
22" 55 cm	XX XX	1103-20XX-OHC	<input type="checkbox"/>
	18" 45 cm	1103-2218-OHC	<input type="checkbox"/>
	20" 50 cm	1103-2220-OHC	<input type="checkbox"/>
24" 60 cm	22" 55 cm	1103-2222-OHC	<input type="checkbox"/>
	18" 45 cm	1103-2418-OHC	<input type="checkbox"/>
	20" 50 cm	1103-2420-OHC	<input type="checkbox"/>
Up to 20" 50 cm	22" 55 cm	1103-2422-OHC	<input type="checkbox"/>
Up to 25" 63 cm	XX = _ _ _ _	1103-2929-OHC	<input type="checkbox"/>
	XX = _ _ _ _	1103-2939-OHC	<input type="checkbox"/>


PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit



2. SELECT COLOUR

Your colour selection applies to any extra covers ordered



3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
<input type="checkbox"/>	Extra Outer Cover (A)	
<input type="checkbox"/>	Extra Incontinence Cover (B)	
<input type="checkbox"/>	Rigidiser Base Panel (can be inserted inside cover for a rigid surface)	
<input type="checkbox"/>	5° Wedge 1209-2724-300	

4. CUSTOMISE

Refer to Cushion Customisations Page

Cushion Customisations


This page details common customisations, for full custom options please refer to the Spex Custom Cushion Script Form




LEG LENGTH DISCREPANCY

1111-0292-000

MEASUREMENTS




REFERENCE DRAWING




FRONT EDGE CUT-BACK FOR LEG CONTRACTURES

1111-0317-000

MEASUREMENTS




REFERENCE DRAWING



REMOVE MEDIAL THIGH SUPPORT

1111-0335-000


REFERENCE DRAWING



SPECIFY MEDIAL SUPPORT DIMENSIONS

1111-0294-000

MEASUREMENTS




REFERENCE DRAWING




SPECIFY LATERAL THIGH SUPPORT DIMENSIONS

1111-0295-000

MEASUREMENTS




REFERENCE DRAWING



INCREASE REAR HALF HEIGHT TO FULL LATERAL THIGH SUPPORT (SUPERHIGH CONTOUR ONLY)

1111-0304-000


REFERENCE DRAWING



ADDITIONAL PRESSURE RELIEVING OVERLAY

1111-0391-000

REFERENCE DRAWING



SUPRACOR® HONEYCOMB

1111-0431-000