


Walking Rein

Curious individuals with challenging behaviour can suddenly slip out of your sight despite your best intentions. A walking rein will help you keep your loved ones close.

Help control children and adults when out walking



Parent/carer's rein is attached to the side straps which can act as a directional aid, helping to guide the user

May prevent personal harm by offering a certain **degree of restraint**

Chest pad is made from fire-retardant PVC and is easy to **wipe clean**

Particularly good for those with a **limited sense of safety awareness**



Ideal to use as a back-up safety method to keep someone within reach, especially during walks to the park where the dangerous roads are never far or in shopping centres where lots of people are distracting.

It is important for safety reasons that the rein strap is always held behind the child and not pulled from the front. This ensures the strain is on the very strong webbing belts and not on the PVC chest pad.

*Note: Never leave your child unattended. This product is not a substitute for proper adult supervision

VERSIONS AVAILABLE

STANDARD REIN

Has waist belt plus chest strap and detachable rein strap



ESCAPOLOGIST REIN

For the very active person. Fitted with a chest belt, waist belt and shoulder straps. Has adjustable and detachable rein.



CUSTOM REINS

Medifab modifications available to meet your custom requirements. Should you require a specific walking rein to suit behavioural needs of young children to adolescents/adults, contact Medifab to discuss available options for your situation.

| Model | Code | Age/Size |
|--------------------|---------------|---|
| Standard | 4605-0062-017 | Fits most ages, up to 97 cm waist circumference |
| Small Escapologist | 4605-0065-017 | Up to 8 yrs |
| Large Escapologist | 4605-0066-017 | 9 yrs to adult |