

Hero

SAFETY FEATURES

Swivelling special
needs carseat
for forward or
rear-facing



HEIGHT ADJUSTABLE HEADREST ensure this is adjusted as your child grows taller

PADDING ACCESSORIES can be set up to ensure effective support to suit your child's needs

Backrest bar has **TWO ANGLE SETTINGS** - ensure this is adjusted to suit the vehicle

RED INDICATOR STRIPS show when ISOFIX, swivel or floor leg isn't engaged and ready for travel

TILT HANDLE Ensure adequate tilt setting is applied to help protect your child's head and neck

Hero

Hernik Hero car seat is adapted from the Britax range and is specifically modified to accommodate for special needs children.

- The Hero is designed to be used in the back seat of a vehicle and **does not require use of the vehicle seat belt; the Isofix anchorage system providing a secure attachment to the vehicle.**
- **Rear setting must be used for a child of 0-15 months,** although in NZ recommendation is until 2 years old.
- Install the child seat in the car in such a way that it is **not trapped by the front seat or the vehicle doors.**
- The **support leg must be extended to the appropriate length** to suit the vehicle.
- The **belts must not be twisted or trapped** and must be completely **tightened.**
- **Before each journey,** check the Isofix is secured, support leg in place, and swivel locked.
- Ensure child has **not outgrown the weight capacity which is 18kg.** If child is approaching this weight limit, consider what car seat would be next option.
- Periodically **check for loose or missing parts.**



Ensure child has not outgrown the weight capacity limit

Contact the car seat specialists for more information:
0800 543 343
solutions@medifab.com

Distributed by **medifab**[®]