


Kidsflex



Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Basic Packages

Size 2	Size 3
<p>Includes:</p> <ul style="list-style-type: none"> • Blue • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20 degs) • Isofix attachment • Standard head support • Small fixed lateral supports • Slim left turn swivel base • Seat belt guide on the left 	<p>Includes:</p> <ul style="list-style-type: none"> • Grey • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20 degs) • Isofix attachment • Standard head support • Large fixed lateral supports • Slim left turn swivel base • Seat belt guide on the left
4170-0022-060	4170-0023-090

Specifications

Size	Total Height	Seat to Head	Back Height	Seat Width	Seat Depth	Lateral Sizes	Distance of Lateral Supports	Weight	Max User Weight
2	80 cm	73 cm	41-57 cm	35 cm	34-39 cm	26-36 cm	19-32 cm	8.9 kg	36 kg
3	90 cm	82 cm	47 - 65cm	40 cm	40-46 cm	38-48 cm	24-37 cm	12.3 kg	75 kg

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Kidsflex_20230118_AU